

*Committed to supporting independent living!*



# HOME HEALTH CATALOGUE

## 2014 - 2015



**BOWERS**  
HOME HEALTH  
(A DIVISION OF RADION LABORATORIES)

# BOWERS HOME HEALTH LOCATIONS



## BRITISH COLUMBIA

Unit 7 - 3691 Viking Way  
Richmond, BC V6V 2J6  
Tel: 604.278.7566  
Toll Free: 1.800.663.0047  
Fax: 604.278.7525



## ALBERTA

11783 - 186th STREET  
Edmonton, AB T5S 2Y2  
Tel: 780.454.1666  
Toll Free: 1.800.561.3687  
Fax: 780.451.6762


### Hours of Operation:

#### BC

Monday - Thursday 8 am - 5:30 pm  
Friday 8 am - 4:30 pm  
*Closed weekends & statutory holidays*

#### Alberta

Monday - Friday 8 am - 4:30 pm  
*Closed weekends & statutory holidays*



Email us at:

[HomeHealth@bowersmedical.com](mailto:HomeHealth@bowersmedical.com)

## How to Order:

Placing an order with Bowers Home Health is easy and convenient. We have customer service agents across Canada waiting to serve you by phone, fax, email, online, B2B and EDI. Bowers Medical Supply has over 40 years of experience in the medical distribution industry and we are now pleased to provide the everyday community with their home care needs.

No account required or minimums to meet. Some orders might be eligible for prepaid freight.

**Questions? Call or email us!**  
**[HomeHealth@bowersmedical.com](mailto:HomeHealth@bowersmedical.com)**

## Terms & Conditions:

### Returns:

Return requests must be authorized through your local Bowers Home Health location.

All products must be returned within 30 days of receipt of authorization.

Products may be subject to restocking fees.

Custom and special order products are not returnable.

### Defects & Damages:

Claims for defects or shortages must be made within 24 hours of receipt of goods

### Payment:

Equipment and special orders may be subject to prepayment or deposit. All delivered orders must be paid prior to shipping. Applicable freight and taxes are extra. We accept all major credit cards. Debit & cash accepted for pick ups.



Follow us on FaceBook for the latest deals, events and tradeshow.

### We also submit claims for:

- DVA - Veterans Affairs Canada
- NIHB - Non-Insured Health Benefits
- SAIL - Saskatchewan Aids for Independent Living
- AADL - Alberta Aids for Daily Living
- BC Ministry of Housing & Social Disability
- Ontario Disability Program & Ministry of Community & Social Services

# EQUIPMENT RENTALS

**Bariatric Hospital Bed** ..... \$386.00 / month ..... \$90.00 / week

*Includes: Mattress & 1/2 headend side rails*

**Standard Hospital Bed** ..... \$150.00 / month ..... \$75.00 / week

*Includes: Mattress & 1/2 headend side rails*

*Optional Upgrade: 3/4 headend side rails*

Overbed Table ..... \$40.00 / month ..... \$25.00 / week

Bed Assist Rail ..... \$45.00 / month ..... \$25.00 / week

Roho Cushion w/ cover ..... \$90.00 / month ..... \$50.00 / week

Roho Cushion w/out cover ..... \$80.00 / month ..... \$40.00 / week

Walker Folding w/ wheels ..... \$35.00 / month ..... \$17.50 / week

Walker Folding w/out wheels ..... \$25.00 / month ..... \$12.50 / week

**Four Wheel Walker w/ slow down brakes** ..... \$60.00 / month ..... \$30.00 / week

**Wheelchair 18"** ..... \$60.00 / month ..... \$30.00 / week

Transport Chair - ultra light ..... \$55.00 / month ..... \$25.00 / week

Wheelchair Ramp - 3' x 30" ..... \$50.00 / month ..... \$25.00 / week

Crutches - All sizes ..... \$18.50 / month ..... \$7.50 / week

Forearm Crutches ..... \$35.00 / month ..... \$17.50 / week

IV Pole ..... \$20.00 / month ..... \$10.00 / week

**Patient Lift** ..... \$120.00 / month ..... \$70.00 / week

ZONE 1: \$50.00 one way, \$80 both ways - Vancouver, Richmond, Coquitlam, East Vancouver

ZONE 2: \$60.00 one way, \$100 both ways - Burnaby, North Vancouver, West Vancouver Horseshoe Bay and Deep Cove

ZONE 3: \$70.00 one way, \$120 both ways - Surrey/Delta, White Rock, Langley, Cloverdale, Horseshoe Bay to Squamish and Maple Ridge

ZONE 4: \$90 one way, \$160 both ways - Outside of Zones 1, 2 and 3.

A valid credit card is required on all rentals.

**Not sure if you want to purchase an item? Then take advantage of our rent-to-own program! Under this program, 100% of the first six month's rent goes towards buying the equipment.**



# CONTENTS

## INCONTINENCE

Liners, Pads & Briefs.....	1 - 8
Wipes.....	9
Underpads & Mattress Protectors.....	10 - 13

## BABY CARE

Diapers.....	14
Pacifiers.....	15

## ODOR CONTROL

Urine Free.....	16
Odor Eliminators.....	17

## SKIN CARE

Cleansers, Moisturizers & Protectants.....	18 - 26
Wound Cleansers & Skin Preps/Antiseptics.....	27 - 28

## PERSONAL HYGIENE

Mouth Swabs & Cleansing Products.....	29
Nail Products & Patient Care.....	30

## BEDSIDE PRODUCTS

Bed Pans & Urinals.....	31
Wash Basins.....	32
Enema Products.....	33

## CATHETERS

Foley Catheters.....	34 - 35
Intermittent Catheters.....	36 - 37
External Catheters.....	38
Lubricants & Clamps.....	39
Catheter Trays.....	40

## URINARY DRAINAGE

Drainage Bags & Accessories.....	41 - 43
----------------------------------	---------

## STABILIZATION DEVICES

Catheter & Multi Purpose Tube Holders.....	44 - 46
--	---------

## BATHROOM SAFETY

Toileting Aids.....	47
Showers Chairs & Bath Boards.....	48
Tub rails & Grab Bars.....	49 - 50

## AIDS FOR DAILY LIVING

Reachers & Stocking Aids.....	51
TableWare.....	52 - 55
Clothing Protectors & Hair Wash Trays.....	56

## MOBILITY AIDS

Transport Chairs & Wheelchairs.....	57
Walkers & Accessories.....	58 - 64
Crutches, Canes & Accessories.....	65 - 66
Uplift Seats.....	67

## FALL PREVENTION

Patient Alarms.....	68
Fall Mats.....	69
No-Slip Socks & Transfer Belts.....	70 - 71
Hip Protectors.....	72 - 74

## PRESSURE RELIEF

Seat Cushions & Mattress Overlays.....	75 - 77
Heel Relief.....	78 - 81
Roho Cushions.....	82
Bed Wedges.....	83



# CONTENTS

## BEDS

Electric Beds & Mattresses.....	84 - 88
Alternating Air Mattresses.....	89 - 90

## EQUIPMENT

Patient Lifts & Slings.....	91
Scales & Privacy Screens.....	92
Recliners, IV Hooks & Overbed Tables.....	93 - 94
Bed Assist Rails.....	95
Pill Crushers.....	96
Autoclave.....	97
Feeding Pumps & Bags.....	98

## SYRINGES & NEEDLES

Syringes without Needles.....	99 - 100
Insulin & Safety Syringes.....	101 - 102
Needles & Winged Infusion Sets.....	103 - 104
Tourniquets.....	105
IV Solutions & Administration Sets.....	106
Sharps Disposal.....	107 - 108

## GLOVES

Exam Gloves - Sterile & Non-Sterile.....	109 - 113
Surgical Gloves.....	114
Glove Box Dispensers.....	115 - 117

## PERSONAL PROTECTION

Shoe Covers, Bouffants & Gowns.....	118
Masks.....	119

## HAND HYGIENE

Manorapid & Dermotan.....	120
Hand Sanitizers & Accessories.....	121 - 126

## DISINFECTANTS

Surface Disinfectants.....	127 - 128
High Level Disinfectants & Detergents.....	129 - 130

## FIRST AID

First Aid Kits, CPR Masks & Guedels.....	131
Hot/Cold Compresses, Bandages, Eye Wash.....	132
Splints.....	133
Adhesive Bandages.....	134

## TAPES

Paper, Transparent & Hypoallergenic.....	135
Athletic & Specialty Tapes.....	136

## WOUND CARE

Gauze (ST & N/S) & Conforming.....	137
Cohesive Bandages.....	138 - 139
Mesalt, Alldress & Mepore.....	140
Mepilex.....	141 - 142
Melgisorb, Mesorb & Mepitel.....	143
Allevyn.....	144
Intrasite, Ichthopaste & Acticoat.....	145
Jelonet, Bactigras, Iodosorb,	
Inadine & Dressing Trays.....	146
Tegaderm, Duoderm, Adaptic & Aquacel.....	147
Promogran.....	148
Silvercel, Actisorb, Packing Strip & Cica-Care.....	149
Tubular Gauze & Opsite.....	150

## CASTING

Gypsona.....	151
--------------	-----

## COMPRESSION THERAPY

Profore, Surepress & Comprilan.....	152
Compression Garments & Accessories.....	153 - 157

## BRACES & SUPPORTS

Elbow Straps.....	158
Wrist & Abdominal Supports & Cervical Collars.....	159
Knee Supports, Hernia Aids & Arm Slings.....	160
Back Supports.....	161



# CONTENTS

## REHAB & FITNESS THERAPY

Exercise Bands, Tubing & Accessories.....	162
Exercise Balls.....	163
Shoulder Pulley & Hand Exercisers.....	164
Flexbar & Pedal Exercisers.....	165
Moist Heat Packs & Covers.....	166
Adhesive Felt.....	167
Conductive Pastes & Gels.....	168
Ultrasound Gels & Warmers, Tens.....	169
Tens Electro-Stimulators.....	170 - 171
Electrodes.....	172 - 173

## RESPIRATORY

Lung Flute.....	174 - 177
Oxygen Tubing & Masks.....	178
Nebulizers.....	179

## DIAGNOSTICS

Pulse Oximeters.....	180
Stethoscopes.....	181 - 183
Automatic Blood Pressure Units.....	184 - 185
Sphygmomanometers & Manual BP Units.....	186 - 190
Vital Signs Devices.....	191 - 192
Diagnostic Sets.....	193
Otoscopes & Ophthalmoscopes.....	194 - 195
Thermometry & Probe Covers.....	196 - 199
Speculums.....	200 - 201

## INSTRUMENTS

Forceps & Scissors.....	202
Scalpel blades & Handles.....	203

## DIABETIC SUPPLIES

Glucose Tablets & Gel.....	204
Glucose Meters & Lancets.....	205

## GENERAL SUPPLIES

Cotton Tipped Applicators.....	206
Tongue Depressors.....	207
Straws & Cotton Balls.....	208
Table Paper & Medication Cups/Spoons.....	209
Paper Products & Dispensers.....	210 - 216



## DAY LIGHT PADS



Pads designed for light to moderate bladder and/or bowel incontinence protection. As a 2-Piece system, these pads provide a more discreet, body-close fit when worn with TENA knit or comfort Pants. TENA Day Light Pads for mild to moderate bladder control protection are designed to be used with washable TENA Reuseable Pants. They are anatomically-designed for comfort and security. TENA Day Light Pads also feature an adhesive strip.

PROM62314

AADL: M242

## TENA® ACTIVE LINERS - REGULAR



Many things can trigger the Unexpected Leak™—pregnancy, childbirth, prescription drugs. Even laughing, coughing or yes, the natural part of aging. TENA® ACTIVE™ Liners are designed for those very light leaks, but they'll protect you better than a regular feminine protection pantiliner so you stay dry no matter what.

PROM56300

## TENA® HEAVY PADS - LONG



Heavy leaks can make it hard to sit through a movie or walk through the park. Get Fearless Protection™ with longer pads that provide extra coverage when you need it most.

PROM47300

AADL: M242

## TENA® ACTIVE™ ULTRA THIN PADS



Finally, an Ultra Thin pad that moves with you! Whether you're out on a run or heading to work, do it with the Fearless Protection™ of TENA®. Don't miss a beat of your active, on-the-go life.

REGULAR - PROM46500

Don't let life's Unexpected Leak™ slow you down. TENA® ACTIVE™ Ultra Thin Long Pads provide a little extra coverage so you have Fearless Protection™ whether you're dancing, driving or doing it all.

LONG - PROM48200

## TENA® ANYWHERE™ ULTRA THIN PADS REGULAR



Enjoy all-day comfort no matter what the day brings or where it takes you. Whether you're in the car or the check-out line, TENA® ANYWHERE™ Ultra Thin Pads lock in the Unexpected Leak™ giving you the Peace of Mind to live your on-the-go life.

PROM52060

## TENA® ANYWHERE™ ULTRA THIN PADS LONG



The Unexpected Leak™ happens to most women at one time or another. And it can happen anywhere. Stuck in traffic? Trapped in a meeting? Have Peace of Mind with the Fearless Protection™ of TENA® and all-day comfort designed to fit your lifestyle.

PROM52070

## TENA® PADS MODERATE REGULAR



Enjoy Fearless Protection™ from life's moderate leaks. Because TENA® Pads fit in your pocket or purse, you never have to pass on plans—like a glass of wine or a cup of coffee.

PROM39748

## TENA® PADS MODERATE WITH ALOE VERA LONG



Finally, you can feel better about Unexpected Leak™. Get the protection you need and the comfort you want, thanks to a soothing hint of aloe vera.

PROM80600



**TENA® MODERATE LONG**

A little extra coverage can make a big difference. Take a walk, take a nap and take another drink of water knowing you have the Fearless Protection™ of TENA® Pads Moderate Long.

PROM39759

**TENA® ULTIMATE PADS**

Worrying about the Unexpected Leak™ is a thing of the past! TENA® has you covered with Ultimate Pads to protect you from moderate to heavy leakage.

PROM50000

**TENA® OVERNIGHT PADS**

Does the Unexpected Leak™ sometimes keep you from getting a good night's sleep? Rest assured. TENA® has you covered with our top selling Overnight Pads to protect you from moderate to heavy leakage.

PROM57400

**CLASSIC PLUS BRIEF**

Competitively positioned briefs provide comfort and dignity, with secure protection and a trim fit for moderate to heavy bladder and/or bowel protection.

Moderate to heavy absorbent protection and security against leakage, and suitable for fecal incontinence. Hook Tab fastening allows unlimited refastening to ensure a secure, comfortable fit around the body to minimize leakage. TENA Classic Plus also features the InstaDri Skin-Caring System™ technology for outstanding dryness and dual core technology that rapidly absorbs urine and retains it away from the body. Soft nonwoven backsheet for wearer comfort and dignity.

PROM67713 MEDIUM  
PROM67714 REGULAR

**ULTRA BRIEF**

Moderate to heavy absorbent protection and security against leakage, and suitable for fecal incontinence. Now features hook tabs for unlimited refastening to ensure a secure, comfortable fit around the body. Can be continually readjusted without tearing the product. TENA Ultra Briefs also feature the InstaDri Skin-Caring System™ technology for outstanding dryness and dual core technology that rapidly absorbs urine and retains it away from the body. Soft non-woven backsheet for wearer comfort and dignity.

PROM67200 MEDIUM  
PROM67201 REGULAR AADL: M203  
PROM67300 LARGE AADL: M207  
PROM68010 EXTRA LRG

## TENA ULTRA & SUPER BRIEFS



- Fully breathable
- Superabsorbent microbeads help reduce odour and lock in liquid
- Blue target absorption zone locks in liquid and odour where more protection is needed
- Soft, nonwoven backsheet is gentle and less noisy
- Moisture-proof lining protects clothing and linen and is color coded for product selection
- Curved leg elastics for better fit, more protection and comfort
- Wetness indicator to alert caregivers
- Unique TENA® "Sure Stay" wide hook system delivers unlimited, fast, secure fastenability—first time, every time

### ULTRA BRIEFS

PROM67200	MEDIUM	40/PK 2PK/CS
PROM67201	REGULAR	40/PK 2PK/CS
PROM67300	LARGE	40/PK 2PK/CS
PROM68010	XLARGE	15/PK 4PK/CS

### SUPER BRIEFS

PROM67401	MEDIUM	28/PK 2PK/CS
PROM67405	REGULAR	28/PK 2PK/CS
PROM67501	LARGE	28/PK 2PK/CS
PROM68011	XLARGE	15/PK 4PK/CS

## TENA STRETCH ULTRA & SUPER BRIEFS



- Fully breathable
- Stretch sides allow for comfort, accurate sizing and 'stay in place' fit
- Superabsorbent microbeads help reduce odour and lock in liquid
- Blue target absorption zone locks in liquid and odour where more protection is needed
- Soft, nonwoven backsheet is gentle and less noisy
- Moisture-proof lining protects clothing and linen and is color coded for product selection
- Curved leg elastics for better fit, more protection and comfort
- Wetness indicator to alert caregivers
- Unique TENA® "Sure Stay" wide hook system delivers unlimited, fast, secure fastenability—first time, every time

### STRETCH ULTRA BRIEFS

PROM67802	M/R	36/PK 2PK/CS
PROM67803	L/XL	36/PK 2PK/CS
PROM61390	2XL	32/PK 2PK/CS

### STRETCH SUPER BRIEFS

PROM67902	M/R	28/PK 2PK/CS
PROM67903	L/XL	28/PK 2PK/CS

## TENA PROTECTIVE UNDERWEAR EXTRA



Designed for both men and women, for full coverage and extra security, with heavier bladder leakage protection. New Body-Close Fit™ design moves with the body for more comfort and greater wearer dignity.

Designed for heavy urine loss, TENA Protective Underwear Extra provides dignified protection that is secure and worn like normal underwear. The exceptional security is due to the InstaDri Skin-Caring System™ technology, which rapidly absorbs urine into the core and away from the skin. So, even when seated or lying down for long periods of time, the wearer will stay dry and comfortable, and the risk of skin irritation is reduced.

PROM72131  
PROM72231

SMALL  
MEDIUM

PROM72331  
PROM72412

LARGE  
XLARGE

Designed for both men and women, for full coverage and extra security.

New Body-Close Fit™ design moves with the body for more comfort and greater wearer dignity.

Designed for moderate to heavy urine loss, this product is extra secure and worn like normal underwear.

PROM72239  
PROM72339

MEDIUM  
LARGE

AADL: M203  
M207

## TENA PROTECTIVE UNDERWEAR PLUS



## TENA PROTECTIVE UNDERWEAR SUPER PLUS MEN



PROM81780  
PROM81920

M/L  
XL

TENA® MEN™ Heavy Protection Underwear are anatomically designed for men with moderate to heavy bladder weakness.

Specifically designed for the needs of Men, for a body-close, comfortable fit

The efficient, highly absorbent target zone with micro beads ensures protection where it's needed most, while Advanced Odour Protection technology fights odours.

The comfortable waist band is soft, stretchable and narrow like your regular briefs, while the cottony-soft material offers maximum comfort and discretion.

TENA® MEN™ Heavy Protection Underwear comes in Super Plus absorbency and is available in M/L and XL sizes.

## TENA FOR MEN

Specifically designed for the light bladder leakage needs of men. All-day comfort, dryness and security that is designed specifically for men, this cup-shaped product is designed to follow the contours of the male body for a snug, secure fit. The adhesive strip secures the product in its optimum position within close-fitting brief style underwear. TENA for MEN™ is discreet and comfortable protection for men with light urinary incontinence. The unique design locks in liquid to keep it away from the body, plus there's the added reassurance of odour control for feel-good freshness.



PROM50600  
AADL: M242

## BELTED BRIEFS



Multi-Layer Core; Soft, Stay-Dry Topsheet; Elastic, Reusable Straps; Multi-Strand Elastic Leg Gathers; Soft, Cloth-Like Barrier

PROM62900  
AADL: M203

## YOUTH BRIEFS



The TENA Youth Brief offers moderate to heavy bladder and bowel control protection.

PROM61166  
AADL: M200

## COMFORT EXTRA PADS

Used for restorative care  
Pads designed for light to moderate to maximum bladder and/or bowel incontinence protection. As a 2-Piece system, these pads provide a more discreet, body-close fit when worn with TENA knit or comfort Pants.



PROM62321

## COMFORT PADS



InstaDri™ technology locks moisture away & Sure Stay adhesive tabs secure the pad in place for use with regular underwear or TENA® Reusable Knit Pants. Contoured for a more comfortable, body close fit. Gathered sides provide enhanced leakage protection for added security Wetness indicator to alert caregivers

### COMFORT DAY PLUS

PROM62620 46/PK 2PK/CS

### COMFORT NIGHT SUPER

PROM62630 30/PK 2PK/CS

Pads designed for moderate to heavy bladder and/or bowel incontinence protection. As a 2-Piece system, these pads provide a more discreet, body-close fit when worn with TENA Knit or Comfort Pants.

This protection comes in four absorbency levels that are designed to be used with washable TENA Reuseable Pants. They are anatomically-designed for comfort and security, and the Dry-Fast Dual Core™ ensures long-lasting dryness. A wetness indicator on the outside makes it easy to determine whether changing is required and reduces the need for intrusive checking.

## DAY PLUS PADS



PROM62618

AADL: M240

## NIGHT SUPER PADS



Pads designed for moderate to heavy bladder and/or bowel incontinence protection. As a 2-Piece system, these pads provide a more discreet, body-close fit when worn with TENA Knit or Comfort Pants.

TENA Pads for heavy protection come in four absorbency levels designed to be used with washable TENA Reuseable Pants. They are anatomically-designed for comfort and security, and the Dry-Fast Dual Core™ ensures long-lasting dryness. A wetness indicator on the outside makes it easy to determine whether changing is required and reduces the need for intrusive checking.

PROM62718

AADL: M240

## COMFORT PANTS

Designed to support and hold TENA® Pads in place. These washable reusable pants ensure that the TENA Pads stay securely in place. The soft, elastic edges and latex-free material guarantee comfort and can be washed several times without losing their shape. Suitable for both men and women.



### COMFORT PANTS

PROM36044	S/M
PROM36055	L/XL
PROM36066	XXL/XXXL

### SEAMLESS PANTS

PROM63311	S/M
PROM63322	L/XL
PROM63344	XXL/XXXL



## SMALL BRIEFS



Moderate to heavy absorbent protection and security against leakage, and suitable for fecal incontinence. Now featuring hook tabs for unlimited refastening to ensure a secure, comfortable fit around the body. Can be continually readjusted without tearing the product. TENA Small Briefs also feature the InstaDri Skin-Caring System™ technology for outstanding dryness and dual core technology that rapidly absorbs urine and retains it away from the body. Soft nonwoven backsheet for wearer comfort and dignity.

PROM66100

AADL: M200

## BARIATRIC BRIEFS

TENA Bariatric Briefs are designed to meet the specialized needs of larger/ taller individuals, helping to support maximum comfort and dignity with a customizable fit range.

Moderate to heavy absorbent protection and security against leakage, and suitable for fecal incontinence. Now featuring hook tabs for unlimited refastening to ensure a secure, comfortable fit around the body. Can be continually readjusted without tearing the product. TENA Bariatric Briefs also feature the InstaDri

Skin-Caring System™ technology for outstanding dryness and dual core technology that rapidly absorbs urine and retains it away from the body. Soft nonwoven backsheet for wearer comfort and dignity.



PROM61375

XL

AADL: M207

PROM61390

XXL

PROM61385

XXXL

## TENA FLEX



Fully breathable one piece belted brief style promotes less heat build-up against the skin, to help promote overall skin health for those with moderate to heavy incontinence.

A design breakthrough for wearers and caregivers. Easier and quicker to change than traditional products, TENA® Flex™ is highly rated by caregivers of individuals who are bedridden, require assistance or have sensitive skin. This unique belted brief is designed to minimize the need for lifting and has been proven to reduce the risk of back strain on caregivers to make changing less intrusive for the wearer.

PROM67804

SIZE 8

AADL: M200

PROM67806

SIZE 16 AADL: M207

PROM67805

SIZE 12

AADL: M203

PROM67807

SIZE 20

## WET WIPES

These large, soft, alcohol-free wipes come in a handy softpack dispenser.

The special fibres are impregnated with a fluid designed to gently cleanse and moisturize the skin, while removing surface impurities. The washcloth ensures quick, efficient and gentle perineal care at each incontinence product change, and can also be used for gentle washing of the whole body.

Contains dimethicone, aloe and chamomile.



PROM64608 48/PK

## FLUSHABLE WET WIPES



Our NEW TENA® Flushable Washcloths are the newest member of the TENA® disposable washcloth line-up. These washcloths are the ONLY adult-sized disposable washcloths certified flushable by the NSF. This washcloth supports toileting programs, independence and dignity. The TENA® Flushable Washcloth is specially designed to be over 99% biodegradable and will not clog pipes when used as directed, it is sewer and septic safe.

PROM64388 48/PK

## DRY WIPES

Very soft and very strong wipes designed to quickly and gently cleanse the skin. These multi-purpose dry disposable washcloths are ideal for perineal care, but are also suitable for the whole body, and are more hygienic than textile cloth alternatives. They can be used with TENA Skin-Caring Wash Cream for perineal care, or used with TENA Bodywash and Shampoo for all-over cleaning.



PROM74500 50/PK

## DISPOSABLE UNDERPADS



TENA Underpads for bed or chair protection are designed to cover mattresses in order to absorb leakages, reduce odours and maintain dryness. They consist of a polyethylene layer combined with sealed edges to ensure leakage protection. The highly absorbent material is distributed evenly throughout the core.

TENA Underpads provide far greater levels of absorption and retention than draw-sheets or wadding-filled products. Available in a comprehensive range of sizes and absorbencies, TENA Underpads can be bought in bulk for large institutions or in smaller packs for home users and small institutions.

### REGULAR ABSORBENCY

CASC350	12 pkgs of 25 Pcs/Case	43 x 61 cm (17 x 24")
CASC351	8 Pkgs of 25 Pcs/Case	58 x 61 cm (23 x 24")
CASC352	6 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")

### EXTRA ABSORBENCY

CASC353	12 pkgs of 25 Pcs/Case	43 x 61 cm (17 x 24")
CASC354	8 Pkgs of 25 Pcs/Case	58 x 61 cm (23 x 24")
CASC355	6 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")

### ULTRA ABSORBENCY

CASC357	4 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")
---------	-----------------------	-----------------------

### BARIATRIC

CASC361	10 Pkgs of 10 Pcs/Case	91 x 91 cm (36 x 36")
---------	------------------------	-----------------------

**"SURE STAY" BACKING** keeps absorbent pad in place without damaging surfaces with tape, pins or adhesives.

CASC381	6 pkgs of 25 Pcs/Case	57 x 60cm (22.5 x 23.5")
---------	-----------------------	--------------------------



## TENA INSTADRI AIR UNDERPAD



TENA InstaDri Air™ Underpads for bed or chair protection are designed to cover mattresses in order to absorb leakages, reduce odours and maintain dryness.

TENA InstaDri Air™ Underpads provide far greater levels of absorption and retention than draw-sheets or wadding-filled products through the use of a special pin-core structure to allow for even distribution of fluid to eliminate clumps and blotches.

CASC359	10 Pkgs of 7 Pcs/Case	58 X 91 cm (23 x 36")
CASC360	8 pkgs of 5 Pcs/Case	76 X 91 cm (30 x 36")

## TENA AIR FLOW UNDERPAD

- SPECIALLY DESIGNED FOR USE ON LOW AIR LOSS THERAPY BEDS works in tandem with specialty beds to maximize air circulation.
- NONWOVEN TOPSHEET is gentle against the skin.
- BREATHABLE NONWOVEN BACKSHEET allows air to circulate for healthier skin while protecting the mattress from excessive soilage.
- SUPERABSORBENT MICROBEADS help to improve absorbency for greater security and skin dryness.



CASC370	6 Pkgs of 10 Pcs/Case	58 X 91 cm (23 x 36")
---------	-----------------------	-----------------------

## TENA FRESH UNDERPAD



TENA Fresh Underpads for bed or chair protection are designed to cover mattresses in order to absorb leakages, reduce odours and maintain dryness. They consist of a polyethylene layer combined with sealed edges to ensure leakage protection. The highly absorbent material is evenly distributed throughout the core.

CASC373	6 Pkgs of 25 Pcs/Case	55 X 90 cm (21.5 x 35.4")
CASC376	6 pkgs of 25 Pcs/Case	72 X 76 cm (28.4 x 30")

## VINYL MATTRESS PROTECTOR



This mattress protector is a great economical waterproof protection option. The 9 ½ in depth will accommodate most standard mattresses, and it can be removed quickly and easily when changing or cleaning is required. The superior grade fire retardant white vinyl can be easily wiped clean.

- Contour fit – elastic cord sewn into the edge
- Material: 100% Vinyl
- Provides ideal waterproof protection
- 100% superior grade fire-retardant vinyl barrier
- Wipes clean easily
- Latex-free

AMG745430	Hospital, 36 in x 80 in x 9.5 in (91 cm x 203 cm x 24 cm)
AMG745433	Twin, 39 in x 76 in x 9.5 in (99 cm x 193 cm x 24 cm)
AMG745437	Double, 54 in x 76 in x 9.5 in (137 cm x 193 cm x 24 cm)
AMG745443	Queen, 60 in x 80 in x 9.5 in (152 cm x 203 cm x 24 cm)
AMG745447	King, 78 in x 80 in x 9.5 in (198 cm x 203 cm x 24 cm)

## VINYL MATTRESS PROTECTOR ZIPPERED

This mattress protector covers the entire mattress, providing excellent waterproof protection. Made from hospital grade white vinyl, it features a surface that wipes clean easily.  
Zipper fit, hospital.



AMG745405	Hospital, 36 in x 80 in x 6.5 in (91 cm x 203 cm x 16.5 cm)
AMG745450	Twin, 39 in x 76 in x 6.5 in (99 cm x 193 cm x 16.5 cm)
AMG745455	Double, 54 in x 76 in x 6.5 in (137 cm x 193 cm x 16.5 cm)
AMG745460	Queen, 60 in x 80 in x 6.5 in (152 cm x 203 cm x 16.5 cm)
AMG745465	King, 78 in x 80 in x 6.5 in (198 cm x 203 cm x 16.5 cm)



## QUILTED UNDERPADS

Washable and reusable quilted underpads for bed or chair protection. Available in various sizes.

ABS1817	17 x 24"
ABS1882	34 x 36"



## PATIENT POSITIONER



Made of vinyl/100% Nylon  
Slider sheets

ABS-TL90713	98 X 210 CM
-------------	-------------

## TLC PAD

Three-layer pad absorbs and holds urine

- Protects bed linens
- Wicks urine away from sensitive skin
- Soft Cozy Cloth surface reduces friction and shear
- Stands up to repeated launderings
- Strong web handles facilitate positioning and transfer
- 40" x 48"

New CozyCloth cover for reduced friction and shear.



SKIL555012

## ABENA BAMBOO BABY DIAPERS

*Bambo Nature puts both the interests of your child and the environment as our top priorities.*

Bambo Nature nappies are just like ordinary disposables, but they have a much lower eco footprint. Bambo Nature nappies are both ultra-thin and ultra-absorbent with a unique layer that draws wetness away from baby's skin, keeping baby dry and comfortable even with repeat wetting. The nappies have been designed in conjunction with parents to assure the product has a practical design and are user friendly.


Bambo Nature has produced a range of baby diapers offering minimum environmental impact combined with maximum freedom of movement and comfort for the child. The breathable back sheet allows air to circulate to the skin, thereby reducing moisture. For the first time, state-of-the-art raw material technology makes it possible for ABENA to produce eco-labeled, environmentally friendly baby diapers that are just as slim and comfortable as ordinary baby diapers.

### Quality Control

Our diaper production is subject to strict 24-hour quality control in our state of the art laboratory. This ensures that all Abena products meet and exceed international standards.

### Panel Test on Children's Diapers

Abena test samples from every production run to ensure that the product meets its strict specifications.



## GUMDROP PACIFIERS



***The Only GumDrop  
Used in Hospitals!***

Trusted by hospitals, GumDrop® offers the newest design in one-piece pacifiers with the same nipple shape as Soothie®. The contoured pacifier stays clear of baby's nose, safely connects to pacifier attachers and the lightweight construction features the nipple shape babies love.

### Features

1. One-piece design made of durable, latex-free silicone
2. Offers an easy teething solution for front, middle and back teeth
3. Offers your baby a variety of textures to choose from
4. Can be chilled for added relief on tender gums
5. Easily attaches to a pacifier attacher in two different places
6. Phthalates and BPA free



CHS1040013	PREEMIE GREEN	100/CS
CHS1040011-104	NEWBORN MANGO	100/CS
CHS10400124	NEWBORN PINK	100/CS
CHS10400125	NEWBORN BLUE	100/CS

## CHRISAL ALL PURPOSE CLEANER

A New Generation of Cleaning Products for Healthcare

Despite intensive disinfection procedures, recurring outbreaks of pathogens are a major concern. Sanitizers and disinfectants are becoming less effective due to bacteria developing chemical resistance. Once chemical cleaners dry, they leave the surface open to invasion by harmful micro-organisms

CHRISAL probiotic cleaning products provide sustainable solutions to these problems

1. Deep cleaning to the microscopic level and the breakdown of biofilm
2. Continues to work on surfaces for up to 5 days
3. Elimination and prevention of unpleasant odours by excluding odour-causing bacteria
4. Provides a safe and healthy microflora on cleaned surfaces to block re-growth of harmful bacteria

### FEATURES

- Highly concentrated
- Biodegradable
- Enriched with probiotics
- Contains surfactants & enzymes
- No perfumes or scents

### BENEFITS

- Cleans all types of hard surfaces
- Deodorizes while cleaning
- Removes biofilm
- Works for up to 5 days.
- Safe & easy to use
- Creates a healthy environment

### APPLICATIONS

Hard Surfaces

- Flooring
- Plastic & Vinyl
- Counters & Tables,
- Sanitary Fixtures
- Stainless Steel,
- Glass & Mirrors



CHRISL60140-1G

1 GALLON

CHRISL60140-32Z

32 OZ

## METRIMIST



- A freshly scented deoderizer spray for eliminating odors.
- Natural Aromatic Deodorizer
- Environmentally friendly alternative to aerosol air fresheners
- Eliminates offensive odors by cleaning the air, not masking the smell

### Features:

- Biodegradable and nontoxic
- Nonflammable
- Cost effective

METRMX1158      8OZ    12/CS

## ODORGON



Odorgon has a long history of success with no allergy complaints, or reactions. Still, there are situations where any inherent smell of a product is not acceptable, so Quest fills this need. Additional uses include: hunting, fishing & security needs. It's an interesting phenomenon when people use odor eliminators, they want to smell something! So using products with no scent may require training to prevent excessive use.

CHEM83006      114ML      12/CS



## DERMASOFT WITH ALOE HAND AND BODY LOTION



### APPLICATIONS:

- After shower or bath
- Ideal for use in hospitals and nursing homes

### FEATURES:

- Provides moisture to dry skin
- Contains Aloe Vera Gel
- Non-greasy
- Ideal for daily use

ADERM-119	4L	4/CS
ADERM-115	237ml	48/CS

## SUPERSOFT DRY SKINCARE LOTION



### APPLICATIONS:

- For face, hands and body
- After bath or shower
- After shaving

### FEATURES:

- Moisturizes
- Soothes dry, sensitive skin
- 100% fragrance free

ADERM-120P	480ML	24/CS
------------	-------	-------

## GENTLE WASH BODY WASH/SHAMPOO



### APPLICATIONS:

- In the shower or bath
- As a body wash/shampoo
- For handwashing

### FEATURES:

- Multipurpose lotionized skin cleanser
- Enriched with Vitamin E
- Promotes soft, moist and smooth skin
- Gentle enough for everyday use
- Leaves no soap scum or oily residue

ADERM-224	800ML	24/CS
-----------	-------	-------

## MANUAL DISPENSER

- Wall dispenser for 800ml cartridge

ADERM-DISP



## PERISHIELD BARRIER OINTMENT & PROTECTANT



### APPLICATIONS:

- As a barrier ointment
- As a skin protectant to seal out moisture
- As a rash treatment
- For use in hospitals and nursing homes

### FEATURES:

- Antimicrobial skin protectant
- Antiseptic — Protects against infection
- Enriched with Vitamins A, D & E and Aloe
- Zinc Oxide and Chloroxylenol additives help prevent and treat diaper rash due to incontinence
- Baby-fresh scent
- Easy to apply and remove

ADERM-500

120G

24/CS

## PROTECTIVE OINTMENT WITH VITAMINS A&D



### APPLICATIONS:

- For general skin care
- For use in pharmacies, hospitals and nursing homes
- On geriatric patients and children
- For home healthcare

### FEATURES:

- Enriched with Lanolin
- Enriched with Vitamins A & D

ADERM-165

120G

24/CS

## BOWERS SINGLE STEP WASH CREAM

- A gentle all in one cleanser and moisturizer
- Ideal as a massage cream, perineal cleanser or for regular skincare
- No rinsing or drying necessary
- Apply to damp cloth, wash area and wipe away residue



BOWSC000	1 GALLON
BOWSC001	940ML
BOWSC002	250ML

## BOWERS SINGLE STEP ALL IN ONE PERINEAL CLEANSER & MOISTURIZER

- Reduces perineal odour and enhances proper skin care
- No rinsing or drying necessary
- Spray affected area than gently clean skin and wipe away residue



BOWSC003	500ML
BOWSC004	1 GALLON

## BOWERS SINGLE STEP SHAMPOO & BODY WASH

- A gentle blend of quality ingredients that cleans your hair as well as being mild to the skin when used as a body wash



BOWSC006	16OZ
BOWSC007	1 GALLON

## TENA WASH CREAM



Our very popular no-rinse cleanser creates a simple, effective routine for caregivers. This gentle alternative to soap and water cleanses, restores and protects even the most delicate skin. Ideal for incontinence care when frequent cleansing is required, TENA Wash Cream is a mild emulsion designed to cleanse, restore and protect the skin. Documented in epicutaneous tests under dermatological control, the cream has a 20-year history of safe and frequent use around the world.

PROM64331	250ML
PROM64340	500ML
PROM64351	1 L
PROM64355	1 L Fragrance free

## BODYWASH & SHAMPOO



TENA Bodywash and Shampoo is designed to cleanse and condition fragile skin and hair. The extra-gentle formula consists of mild yet effective ingredients that make the hair soft and shiny, and the skin clean and fragrant. The anti-static and conditioning agents are particularly suited for combing the hair and sensitive scalps of elderly people.

PROM64365	500ML
PROM64375	1L

## SECURA MOISTURIZING CLEANSER



Antimicrobial skin cleanser, for perineum or body, aids in the removal of urine and feces or other foreign material. No-rinse, one-step cleansing for perineum and body saves time. Contains aloe and glycerin to condition and moisturize skin, pH balanced. Helps control odor and pediatrician tested.

SMTH59430879 118ML  
SMTH59430979 236ML

## SECURA DIMETHICONE PROTECTANT



Skin protectant. Helps treat and prevent diaper dermatitis. Provides protection from exposure to urine and feces. Protects minor skin irritation associated with diaper dermatitis and helps protect from wetness. Helps prevent and temporarily protects chafed, chapped, cracked or windburned skin or lips. Contains 5% dimethicone.

SMTH59432279 114G

## SECURA PROTECTIVE CREAM



Skin protectant. Helps treat and prevent rash associated with diaper use or continued exposure to feces, urine or both. Protects skin from minor irritation associated with diaper dermatitis and exposure to feces and urine, and helps seal out wetness. Contains 10% Zinc Oxide, Clove oil helps mask odor, Aloe, Vitamin E and Allantoin moisturize and condition the skin. Pediatrician tested.

SMTH59431379 40G  
SMTH59431279 78G

## SECURA PROTECTIVE OINTMENT



Skin protectant. Helps treat and prevent rash associated with diaper use or continued exposure to feces, urine or both. Protects skin from minor irritation associated with diaper dermatitis and exposure to feces and urine, and helps seal out wetness. Contains 98% Petrolatum, Clove oil helps mask odor, Vitamins A, D and E to soothe and condition sensitive skin. Pediatrician tested.

SMTH59431579 70G

## SECURA EXTRA PROTECTIVE CREAM



Skin protectant that helps treat and prevent rash associated with diaper use or continued exposure to feces, urine or both. Protects skin from minor irritation associated with diaper dermatitis and exposure to feces and urine and helps seal out wetness. Contains 30% Zinc Oxide.

SMTH59432479 90G

## SECURA PROFESSIONAL CARE LOTION



Skin moisturizer. Overall body moisturizer that helps maintain skin's natural moisture. Relieves itchiness and discomfort caused by dry skin. Non-greasy, easily absorbed, conditions and soothes dry skin.

SMTH80236 360ML



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687



## PROSHIELD PLUS



### SKIN PROTECTANT

- The only barrier cream indicated for use on intact and injured skin associated with incontinence
- Dimethicone based skin protectant with copolymer bio-adhesives
- Smooth, greaseless, fragrance free barrier
- Adheres to moist and dry skin
- Provides protection from urine, faeces and incontinence associated dermatitis
- Provides protection from friction and shearing forces
- Indicated for clinically dry skin

### FOAM & SPRAY CLEANSER

UNI421750

115G

OINTMENT

UNI421743

235ML

SPRAY

- No-rinse formula
- Gentle and pH balanced
- Minimises irritation of intact or injured skin associated with incontinence
- Formulated to help eliminate odour and break down dried stool
- Moisturisers leave skin supple and hydrated

## BAZA PROTECT



Zinc Oxide and Dimethicone Skin  
Protectant Cream

Special formulation of ingredients that aids  
in the protection, relief and/or treatment of  
diaper dermatitis (rash) from incontinence.  
Safe to use with incontinence briefs.

BAZPRO1880 142G

## ATRAC-TAIN



With Urea and Alpha Hydroxy Acid (AHA)  
Advanced moisturizing relief for severely afflicted  
dry, scaly, cracked skin.  
Ideal for heels, elbows, knees, lower  
leg areas. Cream with 10% Urea  
and 4% AHA cream.

Cream with 10% Urea  
and 4% AHA cream.

Lotion with 5% Urea and 2.5%  
AHA.

SWEEN1802

60G CREAM

SWEEN1814

142G CREAM

SWEEN1806

240ML LOTION

## CRITIC-AID



Clear moisture barrier ointment adheres  
to intact or denuded skin and manages  
moisture. Easy to apply and remove.  
Provides protection, relief and/or  
treatment of diaper dermatitis (rash)  
caused by light to severe or chronic  
diarrhea and enzymatic drainage.  
Safe to use with incontinence briefs.

SWEEN7573 71G

SOLabs has developed the only antimicrobial products specifically designed to prepare, protect, and preserve body art.

SOLabs products contain NO parabens or fragrances. PurKlenz, PurTect & PurSan work effectively and contain the same proven ingredient: Chloroxylenol (PCMX)

Recommended Usage Protocols developed in conjunction with Health Educators Inc. are available at [www.solabs.ca](http://www.solabs.ca)

## PURKLENZ

PurKlenz

PurKlenz topical antimicrobial skin cleanser is formulated to cleanse and sanitize prior to a tattoo or piercing.

SOLAB-SOL-PPKUM 900ML 4/CS



## PURTECT

PurTect

PurTect topical antimicrobial ointment is designed for use during the tattoo to protect Clients and Artists against dangerous contact borne and topical pathogens. Excellent for short-term aftercare.

SOLAB-SOL-PT118ML 118ML 24/CS  
SOLAB-SOL-PT20ML 20ML 40/CS



## PURSAN

PurSan

PurSan topical antimicrobial skin cleanser is recommended for optimal aftercare to clean and preserve body art.

SOLAB-SOL-PSAN118 118ML 34/CS



## REMEDY SKIN REPAIR CREAM



- Remedy® Phytoplex Nourishing Skin Cream helps restore skin's natural moisture balance and barrier properties
- Formulated without mineral oil or other petrochemicals to remain breathable
- Absorbs readily into the skin upon application with no greasy residue
- Moisturizes with a blend of natural emollients including safflower oleosomes and carrageenan
- Soothes the skin with oat kernel and licorice root extract
- Ideal for use on hands, face and body
- Enhanced with Phytoplex, additional botanicals and essential fatty acids

MEDLMSC092404

4 oz

## GLAXAL BASE

New Glaxal Base® Moisturizers combine the natural foundation of earth's ancient essence with scientific innovation to create therapeutic protection for your skin. This specially formulated hypoallergenic premium emollient helps seal in moisture to protect the skin from harsh weather, pollution, and everyday stress which upset the skin's delicate balance and result in the appearance of premature aging. Glaxal Base® soothes and rehydrates very dry and sensitive skin.



TARO-006369100770

450G

## CALMOSEPTINE OINTMENT



The ingredients in Calmoseptine Ointment:

Provide a physical moisture barrier, keeping feces, urine and wound drainage from intact and injured skin.

Act as a counter-irritant, and temporarily relieve discomfort and itching from rashes, eczema, insect bites, and contact dermatitis such as poison ivy, oak and sumac.

Menthol has been shown to cause vasodilation. This may increase circulation to an area and aid in healing.

Mild antiseptic ingredients may help prevent secondary bacterial and fungal infections of the skin.

UNI130906

113G



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687

**SKIN PREP**

SMTH59420425 WIPES 50/BX  
SMTH420279 SPRAY 118ML

**NO STING (ALCOHOL FREE)**

SMTH59420600 WIPES 50/BX  
SMTH66800709 SPRAY 1OZ

**REMOVE ADHESIVE  
REMOVER**

SMTH403120 WIPES 50/BX  
SMTH403379 SPRAY 220ML

**AVEENO LOTION**

Aveeno Active  
Naturals Daily  
Moisturizing  
Lotion.

JOJOV4001 354ML

**CAVILON NO- STING  
BARRIER FILM**

3M3344E WIPES 30/BX  
3M3346E SPRAY 28ML

**ADHESIVE REMOVER**

PDIB16400 WIPES 100/BX

**MELITA'S OINTMENT**

Formulated for incontinence: it repels moisture by acting as a buffer to reduce skin breakdown and infection while healing any existing redness or irritation. It protects and helps relieve chapped or cracked skin. Helps protect skin folds where perspiration may occur.

JAM1460 118ML  
JAM1468 237ML

**SEA-CLENS WOUND CLEANSER**

Cost-effective, saline-based wound cleanser designed to loosen debris from the wound bed without interfering with the healing process

- Saline-based cleanser provides safe and effective irrigation
- Soften eschar to facilitate debridement
- Non-toxic, non-ionic, pH-balanced formula

SWEEN1063 178ML  
SWEEN1061 355ML

**SALINE MODUDOSE**

NACL 0.9% for irrigation.

BALLARD116  
15ML 24/BX

**SALINE SOLUTION**

NACL 0.9% (Normal saline)  
for irrigation.

BAXTJF7633P 500ML  
BAXTJF7634 1000ML

**STERILE WATER**

Sterile water for irrigation.

BAXTJF7623P 500ML  
BAXTJF7624 1000ML

**ALCOHOL PREP PADS**

BOW3001 MED 100/BX  
BOW3002 LRG 100/BX

**BENZALKONIUM WIPES**

Benzalkonium  
antiseptic towelettes.

PDID35164 100/BX



## POVIDONE IODINE SWABSTICKS



PDIS41125 3'S 25/BX

## IODINE PREP PADS



Iodine Prep pads 10%

PDIB40600 100/BX

## PVP 10% ANTISEPTIC SWAB SOLUTION



Each pad is saturated with 10% W/V Povidine Iodine solution, USP which is equivalent to 1% available Iodine  
For external use only  
Air dries quickly requiring no wiping  
Latex-Free  
For cleansing minor scrapes, cuts and burns  
Prepares skin by cleansing dead skin, dirt, or oil  
Helps reduce bacteria

LORIS-MEDPVP 200/BX

## POVIDONE IODINE PREP 10% SOLUTION



### BETADINE SOLUTION

UNI00452243 100 ML

UNI00184002 500 ML

LERNA109-09 115 ML

LERNA109-08 500 ML

**DENTASWABS**

Dentaswab Disposable Oral Swabs with polyester foam tip, are individually wrapped and treated with dentrifice. A refreshing alternative when brushing is not a good option. Provides for proper oral hygiene when on the go or confined to a bed or chair. Choice of Mint or unflavored swabs.

DENT12245 250/BX  
DENT12240 250/BX

WITH DENTIFRICE  
UNTREATED

**LEMON-GLYCERIN  
SWABSTICKS**

PDIS16225 3'S 25/BX

**READYBATH  
SHAMPOO CAPS**

MEDLMSC095230 1/PK 30/CS  
MEDLMSC095231 1/PK (F/F) 30/CS

Each disposable cap is pre-moistened with a rinse-free shampoo & conditioner to cleanse and condition the hair

Works on most lengths and types of hair

This one-step system requires no clean-up

Formula is pH balanced, hypoallergenic and alcohol free

Available with green apple scent or fragrance free

Product can be used at room temperature or warmed for patient comfort

Single patient use

Latex-free

**READYBATH  
BASIC F/F**

MEDLMSC095109 5/PK

Super-soft, medium weight cloths are gentle on the skin

Cloths are pre-moistened with rinse-free formula that cleans, moisturizes, conditions and soothes the skin

Formula is pH balanced, hypoallergenic and alcohol free

Available with a light, gender-neutral scent or fragrance free

Single use cloths help reduce chance of cross-contamination

Eliminates the need for basins, soaps, linens and lotions

Product can be used at room temperature or warmed for patient comfort

Latex-free

**NAIL CLIPPERS****FINGER**

MEDLMEC37224

NO FILE

ALMA4-242

WITH FILE 2"

**TOE**

ALMA4-238

WITH FILE 3"

**DIAMOND DEB FILE**

ALMA4-324

8"

**NAIL CLEANER BRUSH**

- Open back hand grip
- 16 x 7 rows of 1/2" (1.3 cm) bristles

Scrub brushes are made from polypropylene and can withstand temperatures of up to 80°C or 176°F.

AMG118140

12/BX

**DISPOSABLE RAZORS**

BIC720

10/PK 72PK/CS

SCHICK0188

100/BX 10BX/CS

**COMBS**

DON-DC5

5" BLACK 12/BG

**EPSOM SALTS**

UNI626556

450G

UNI00626200

900G

**MOUTH RINSE**

No Alcohol

4 oz 60/cs

CHES000413

## ADULT BEDPAN CONVENTIONAL



**VOLL00080**

2 1/4 QT  
AUTOCLAVABLE  
TURQUOISE

## PONTOON BEDPAN



**PREMH120-07**

2 QT  
TURQUOISE

## FRACTURE BEDPAN



**PREMH100-07**

34OZ  
TURQUOISE

## SITZ BATH



**PREMH990-07**

2 LT BAG WITH 60"  
TUBING, TURQUOISE

**PREMH995-10**

ELONGATED  
2LT BAG WITH 60"  
TUBING, DUSTY ROSE

## MALE URINAL WITH COVER



**VOLL00095**

1 QT (1000ML)  
BLUE, AUTOCLAVABLE

## DELUXE MALE URINAL



**PREM140D-01**

1 QT (1000ML)  
TRANSLUCENT

## MALE SPILL-PROOF URINAL



**EKJO10**

## FEMALE URINAL



**AMG760655**

35OZ (1 LT)

## WASH BASINS - SPONGE SOLUTION STAINLESS STEEL



VOLL87420

1 7/8QT (1.8LT)

VOLL87340

7QT (6.6LT)

VOLL87360

8 1/2QT (8LT)

## WASH BASINS/SPONGE ROUND



AUTOCLAVABLE

VOLL00100

7QT (6.6LT)

## WASH BASINS ROUND



AUTOCLAVABLE

VOLL00041

5QT (4.7LT)

## WASH BASINS FOOT TUB - RECTANGULAR



PREMH362-07

7QT (6.6LT)

TURQUOISE



PREMH360-10

6QT (5.7LT)

DUSTY ROSE

## PITCHER WITH 9OZ CUP COVER



VOLL00110

1QT (0.95LT)

BLUE

AUTOCLAVABLE

## EMESIS BASIN

VOLL00060

16OZ (500CC)

TURQUOISE, AUTOCLAVABLE



PREMH300-07

16OZ (500CC)

TURQUOISE



## ENEMA BAGS



Dover™ Enema Bag and Tube, Pre-Lubricated Tip,  
24 Fr/Ch (8.0 mm) x 60" (1.5 m)

KEND145540	48/CS
MDI50-800	50/CS


## GENT-L-TIP ENEMA



The Gent-L-Tip Enema carries the National Brand Equivalent (NBE) emblem. It features a soft, prelubricated tip for patient comfort. A specially designed diaphragm valve prevents bottle reflux. Ready to use.

CHES101503BG 48/CS

## BARD LUBRICATH FOLEY CATHETERS



The BARDEX® LUBRICATH® latex Foley catheter provides for smoother insertion and enhanced patient comfort through the hydrogel coating's unique ability to absorb fluid, thus creating a hydrophilic "cushion" between the catheter surface and the urethra.

BARD365714	14FR	5CC	12/CS
BARD365716	16FR	5CC	12/CS
BARD365718	18FR	5CC	12/CS
BARD365720	20FR	5CC	12/CS
BARD365722	22FR	5CC	12/CS

## BARD SILASTIC LATEX FOLEY CATHETERS

Smooth, non-stick silicone exterior coating reduces calcification build-up. Firm, but flexible, latex-based construction with large drainage lumen.

AADL: M341



BARD33612CE	12FR	5CC	10/CS
BARD33614	14FR	5CC	10/CS
BARD33616	16FR	5CC	10/CS
BARD33618	18FR	5CC	10/CS
BARD33620	20FR	5CC	10/CS
BARD33622	22FR	5CC	10/CS
BARD33624	24FR	5CC	10/CS

## UNCOATED SILICONE FOLEY CATHETERS

The BARD® 100% latex-free all-silicone Foley catheter is available in a variety of French sizes and balloon sizes, including pediatric.

AADL: M345



BARD165814	14FR	5CC	12/CS
BARD165816	16FR	5CC	12/CS
BARD165818	18FR	5CC	12/CS
BARD165820	20FR	5CC	12/CS

## AMSure® Silicone Coated Latex Foley Catheters



Two-Way, Silicone Coated Latex Foley Catheters, Individually Packed, Sterile

AADL: M341

AMSAS41012	12FR	5CC	12/BX
AMSAS41014	14FR	5CC	12/BX
AMSAS41016	16FR	5CC	12/BX
AMSAS41018	18FR	5CC	12/BX
AMSAS41020	20FR	5CC	12/BX

AMSAS42012	12FR	30CC	12/BX
AMSAS42014	14FR	30CC	12/BX
AMSAS42016	16FR	30CC	12/BX
AMSAS42018	18FR	30CC	12/BX
AMSAS42020	20FR	30CC	12/BX

## AMSure® 100% Silicone Foley Catheters



Two-Way, 100% Silicone Foley Catheters, Individually Packed, Latex-Free, Sterile

AADL: M345

AMSAS41012S	12FR	5CC	12/BX
AMSAS41014S	14FR	5CC	12/BX
AMSAS41016S	16FR	5CC	12/BX
AMSAS41018S	18FR	5CC	12/BX
AMSAS41020S	20FR	5CC	12/BX

AMSAS42012S	12FR	30CC	12/BX
AMSAS42014S	14FR	30CC	12/BX
AMSAS42016S	16FR	30CC	12/BX
AMSAS42018S	18FR	30CC	12/BX
AMSAS42020S	20FR	30CC	12/BX

## ALL PURPOSE ROBINSON/NELATON PVC



This is an all purpose PVC catheter with 2 staggered eyes. It is color coded and is 40 cm (16") long. It is for single use. In an effort to reduce costs, this economical catheter is packaged in proper but simple packaging.

RUSCH238500-12	12FR	100/BX	WHITE
RUSCH238500-14	14FR	100/BX	GREEN
RUSCH238500-16	16FR	100/BX	ORANGE
RUSCH238500-18	18FR	100/BX	RED
RUSCH238500-20	20FR	100/BX	YELLOW
RUSCH238500-22	22FR	100/BX	LILLAC

AADL: M329

## ROBINSON PLUS CLEAR CATHETERS



Robinson Plus Clear low friction plastic, w/connector, 16", 2 polished eyes, Sterile.

MEDRX675110	10FR	100/BX
MEDRX675112	12FR	100/BX
MEDRX675114	14FR	100/BX
MEDRX675116	16FR	100/BX
MEDRX675118	18FR	100/BX

## CLEAR PLASTIC CATHETERS



Clear Plastic Cath. w/ Connector, 16", 2 eyes, Sterile.

MEDRX605012	12FR	100/BX
MEDRX605014	14FR	100/BX
MEDRX605016	16FR	100/BX
MEDRX605018	18FR	100/BX

AADL: M328

## VINYL URETHRAL CATHETERS



Latex-Free, Individually Packed in Paper Poly Pouch, Sterile

### FEMALE

AMSAS860612	12FR	6"	50/CS
AMSAS860614	14FR	6"	50/CS
AMSAS860616	16FR	6"	50/CS

### MALE

AMSAS861612	12FR	16"	50/CS
AMSAS861614	14FR	16"	50/CS
AMSAS861616	16FR	16"	50/CS

## SELF-CATH INTERMITTENT CATHETERS

Siliconized catheters with fire-polished eyelets for intermittent catheterization.

100% Latex-free. Sterile.



MENT212	12FR	6" FEMALE STRAIGHT TIPPED	30/BX
MENT240	12FR	6" FEMALE STRAIGHT TIPPED, LEUR END	30/BX

## SELF-CATH PLUS INTERMITTENT CATHETERS

Siliconized catheters with fire-polished eyelets for intermittent catheterization. The hydrophilic coating activates with exposure to water for low friction catheterization. The Self-Cath Plus has an uncoated grip zone for better control and handling.

100% Latex-free. Sterile





## Rochester Original All-Silicone Male External Catheters

The first self-adhering male external catheter silicone sheath. Skin-friendly, non-sensitizing adhesive integrated into the silicone sheath lining.



Offers built-in adhesive in a length and placement found with traditional male external catheters. The soft, breathable sheath conforms to the skin.

### ULTRAFLEX

ROCH33101	SMALL (25MM)	100/BX
ROCH33102	MEDIUM (29MM)	100/BX
ROCH33103	INTERMEDIATE (32MM)	100/BX
ROCH33104	LARGE (36MM)	100/BX
ROCH33105	XLARGE (41MM)	100/BX

AADL: M301



Combines forward adhesive technology with a 1.5" sheath. The maximum adhesive area on a short sheath provides men with a catheter for their specific needs.

### POP-ON

ROCH32101	SMALL (25MM)	100/BX
ROCH32102	MEDIUM (29MM)	100/BX
ROCH32103	INTERMEDIATE (32MM)	100/BX
ROCH32104	LARGE (36MM)	100/BX
ROCH32105	XLARGE (41MM)	100/BX



Provides 70% more adhesive area than traditional self adhering designs. Combined with patented forward placement of the adhesive, this design helps prevent accidental catheter detachment.

### WIDEBAND

ROCH36101	SMALL (25MM)	100/BX
ROCH36102	MEDIUM (29MM)	100/BX
ROCH36103	INTERMEDIATE (32MM)	100/BX
ROCH36104	LARGE (36MM)	100/BX
ROCH36105	XLARGE (41MM)	100/BX



Ideal for patients who change external catheters several times a day for intermittent catheterization. Each box comes with three reusable foam straps.

### NATURAL

ROCH38301	SMALL (25MM)	30/BX
ROCH38302	MEDIUM (29MM)	30/BX
ROCH38303	INTERMEDIATE (32MM)	30/BX
ROCH38304	LARGE (36MM)	30/BX
ROCH38305	XLARGE (41MM)	30/BX

## LUBRICATING GEL



Ideal for both medical procedures and personal use, MedPro Lubricating Gel is a quality lubricant that is:  
Sterile

- Product is sterile until seal is opened, so as to ensure product safety for medical use.

Water Soluble

- Washes off of skin quickly and easily.

Non Greasy

- Will not stain furniture, clothes or equipment.

Unscented

- For patient or personal sensory comfort.

AMG118305

12/BX



AMG118311

3.5G SACHET

AMG118313

5G SACHET

BOWSUJ58

5G SACHET

## CORD/CATH CLAMP



Individually Packaged  
Latex-Free, Sterile

AMSUCC100


50/CS

## CATHETER TRAYS



BOWTRCA102108A 20/CS

- 1 Ea - Fenestrated Drape
- 1 Ea - Under Pad
- 1 Ea - Pair Vinyl Gloves
- 1 Ea - Basin
- 3 Ea - Pvp Swabs
- 2 Ea - 5gm Lubricant Gel
- 1 Ea - 10cc Pre Filled Syringe
- 2 Ea - Metal Forceps



BOWTRCA102109A 20/CS

- 1 Ea - Fenestrated Drape
- 1 Ea - Under Pad
- 1 Ea - Prep Cup
- 1 Ea - Basin
- 5 Ea - Rayon Balls
- 2 Ea - 5gm Lubricant Gel
- 1 Ea - 10cc Pre Filled Syringe
- 2 Ea - Metal Forceps

***\*\*Custom trays are available upon request!***

## CONVEEN SECURITY LEG BAG WITH NON-LATEX FABRIC LEG BAG STRAPS



Lightweight leg bag includes pre-attached, adjustable, anti-kink tubing, soft cloth backing, anti-reflux valve, and centre weld for even distribution of urine.

Latex and PVC free



SWEEN5167

26oz (750ml)

Large w/20" Tubing

10/bx

## CONVEEN SECURITY LEG BAG WITH NON-LATEX FABRIC LEG BAG STRAPS



Deluxe leg bag with pre-attached adjustable, standard tubing, soft cloth backing and anti-reflux valve.

Latex free



SWEEN5174

28oz (800ml)

Extra Large w/18" Tubing

10/bx

## AMSINO LEG BAGS



### PUSH-PULL DRAIN VALVE

AMSAS306N 600ML 48/CS

AMSAS607A 900ML 48/CS

AADL: M354

### TWIST-TURN DRAIN VALVE

AMSAS308N 600ML 48/CS

AMSAS309N 900ML 48/CS

AADL: M354

## LEG BAG EXTENSION TUBE

PVC Extension Tube for Leg Bag  
18" Flexible, Kink Resistant  
With Adapter and Protective Cap  
Latex-Free, Sterile



AMSAS311

50/CS

## DOVER URINE LEG BAG

KEND145514	17OZ (503ML)	50/CS
KEND145516	25OZ (730ML)	50/CS

### URINE BAG WITH 12" EXTENSION TUBING

KEND3433	17OZ (500ML)	50/CS
----------	--------------	-------



## BARD DISPOZ-a-BAG

BARD® DISPOZ-a-BAG® Leg Bags are made from heavy duty vinyl for extended wear, odor containment and leak resistance. Each leg bag contains an anti-reflux valve. Available with or without the FLIP-FLO™ drainage valve which opens and closes with a simple flip. Reusable, sterile

BARD150180	19OZ 18" TUBING	50/CS
BARD150181	32OZ 18" TUBING	50/CS



## BARD DRAINAGE BAGS

BARD's Urine Drainage bags feature a rounded, teardrop shape that ensures more efficient drainage.

Latex-free and bacteriostatic drainage bags are also available. The bacteriostatic collection bag is impregnated with a novel antimicrobial agent that resists pathogens commonly associated with CAUTIs.

Our drainage bags are available in a wide range of sizes including larger capacities.

BARD153504C	2000ML	20/CS
BARD153561-S	2000ML	20/CS

AADL: M352



## KENGUARD DRAINAGE BAGS

Dover™ Urine Drainage Bag, Needle Sampling,  
Anti-Reflux Chamber, Drain Tube, Poly Bag, 2000 mL

KEND3512                      2000ML                      20/CS

Dover™ Urine Drainage Bag, Needle Sampling,  
Anti-Reflux Chamber, Drain Tube Hook and Loop  
Hanger, Poly Bag, 2000 mL

KEND3512V                      2000ML                      20/CS



KEND6206                      2000ML                      20/CS

## CURITY ANTI-REFLUX DRAINAGE BAGS



Curity drainag bags with Splashguard drain spout,  
which helps prevent backsplash of urine.  
SAFEGUARD® needleless sampling port, which is  
accessed with a standard Luer slip syringe or blunt  
cannula.


KEND8206                      2000ML                      20/CS



## MULTI-PURPOSE TUBE HOLDER

CATH-SECURE® is the original multi-purpose tube holder. It provides an efficient, standardized system for securing catheters and medical tubing. The device is easy to adjust and may remain in place for approximately one week without altering skin integrity. CATH-SECURE® is made with alcohol-soluble adhesive.

Unit of Measure: 2 1/2" Long Tab & 3" x 1 3/4" Wide Base




MCJO54452 50/BX  
AADL: M312

## MULTI-PURPOSE TUBE HOLDER

The difference between the CATH-SECURE® 5445-3 & 5445-2 is the length in the tab that wraps around the lumen (see spec sheet). They both provide an efficient, standardized system for securing catheters and medical tubing. The device is easy to adjust and may remain in place for approximately one week without altering skin integrity. CATH-SECURE® is made with alcohol-soluble adhesive.

Unit of Measure: 2 3/4" Long Tab & 3" x 1 1/4" Wide Base



MCJO54453 50/BX  
AADL: M312

## DUAL TAB MULTI-PURPOSE TUBE ANCHORING DEVICE

The CATH-SECURE DUAL TAB® is designed to extend from the center, wrap around the tubing, and anchor it tightly to the tape using a press-on hook-and-loop fastening system. The CATH-SECURE DUAL TAB® improves patient comfort and mobility, is lightweight and flexible, and most of all CDC compliant.

Unit of Measure: (2) 3" Wrap Tabs x 3" Wide Base




MCJO54454 50/BX  
AADL: M312



## WATER-RESISTANT / BREATHABLE TUBE HOLDER

The CATH-SECURE Plus® is made using a breathable, water-resistant, hypoallergenic material that features a butterfly design. The device above is DEHP-Free and Latex Free. This product line improves overall clinical outcomes, lowers healthcare costs, increases nursing productivity all while providing a better quality of care for patients and meeting the healthcare industry standard of care.

Unit of Measure: 2 1/2" Long Tab & Butterfly Base



MCJO54456

50/BX

AADL: M312

## CATH-SECURE FOR KIDS™

The CATH-SECURE For Kids™ is a multi-purpose tube holder. It provides an efficient, standardized system for securing pediatric catheters and medical tubing. The device is hypoallergenic, DEHP-FREE and Latex Free. The CATH-SECURE For Kids™ lowers healthcare costs and increases nursing productivity while providing quality care for patients and meeting the healthcare industry standard of care. The CATH-SECURE For Kids™ can be purchased through your preferred medical distributor.

Unit of Measure: 1 3/4" Long Tab & 1 1/2" x 3/4" Wide Base



MCJO73342

50/BX

AADL: M312

## STATLOCK® FOLEY

STATLOCK® Stabilization Devices are a more effective alternative to tape in helping improve clinical outcomes, quality of care and economic efficiencies.

The STATLOCK® Foley stabilization device accommodates latex 8-22 Fr. and silicone 8-26 Fr. catheters for the ultimate in versatility. Available in adult and pediatric sizes



BARDFOL0102 25/CS

AADL: M312 UPGRD

## STATLOCK® UNIVERSAL PLUS

The StatLock® Universal Plus Stabilization Device sets the new standard of excellence for suture-free stabilization of biliary and nephrostomy percutaneous drainage, and peritoneal lavage catheters.

Advanced technologies now combine with the StatLock® Stabilization Device design sophistication to deliver skin-safe, suture-free, releasable stabilization of even the most lubricious catheters.

The StatLock® Universal Plus Stabilization Device comes in four sizes, to fit virtually all percutaneous drainage catheters:

### Available Sizes:

small (6-8.5 Fr)

medium (10-12 Fr)

large (12-14 Fr)

extra large (14-16 Fr)



BARDVUPD68

MEDIUM 25/CS

## STATLOCK® PICC PLUS

The StatLock® PICC Plus Stabilization Device offers great versatility and ease of use.

Combining the latest in design and engineering technology, StatLock® PICC Plus Stabilization Device releasably secures most winged-catheters on the market.

Just lift the transparent "gull-wings," place catheter wings over the posts, and close.

Anchor pad available in breathable, nonabsorbent tricot polyester or closed cell foam.

Sterile and latex-free.



BARDPIC0222 50/CS

## RAISED TOILET SEAT WITH OR WITH OUT LID



The AquaSense® Raised toilet seat with lid makes it easier to get on and off the toilet.

- With lid for more discretion
- Raises the toilet 2" higher than a regular toilet seat
- Easy-to-clean surface
- Attaches securely to most toilets
- 300 lbs (136 kg) weight capacity
- Easy tool-free installation

AMG770625 2" WITH LID

AMG770626 4" WITH LID

AADL: B220

AMG770620 2" WITHOUT LID

AMG770621 4" WITHOUT LID

## RAISED TOILET SEAT WITH CLAMP & REMOVABLE ARMS



Use without handles for maximum space and easy access.

Use one handle for easy lateral access and a secure grip.

Use both handles for more security and support getting up and down

AMG770618

AADL: B230

## TOILET SAFETY FRAME



Added support for getting safely on and off the toilet

The AquaSense® Toilet Safety Rails make it safer and easier to get on and off the toilet, is aesthetically nice and extremely easy to assemble.

AMG770665

AADL: B255

## COMMUNE



A convenient and safer toilet alternative for persons with reduced mobility. Also use for support over a toilet or as a raised toilet seat.

AMG770315

## ADJUSTABLE BATH SEAT WITH BACK



- Seat size: 20" x 12" (51 cm x 30,5 cm)
- Back size: 20" x 12" (51 cm x 30,5 cm)
- Slip resistant textured surface
- Corrosion resistant aluminum frame
- Sturdy rubber tips - Latex free
- Easy, tool-free assembly: 9 blue knobs
- Weight capacity: 300 lbs (136 kg).

AMG770510  
AADL: B043

## ADJUSTABLE BATH SEAT WITHOUT BACK



AMG770500  
AADL: B024

## TRANSFER BENCH



AMG770406  
AADL: B070

Knocked-Down Blow Molded Bathtub Transfer Bench is built with a light weight anodized rust-resistant frame. Nine (9) adjustments in 0.47 inch increments accommodates different tubs. A removable push-button armrest supports transferring in and out of the tub while a textured molded seat and back prevents slipping and provides better water drainage. A reversible seat / back rest accommodates left or right transfers. Slip resistant tips provide great traction and secure the bench to the floor. Easy, tool-free, push button assembly. Latex-free.

## BATH BOARD

This ergonomically designed bath board fits most bathtubs quickly and easily without the use of tools. What's more, it provides strong support for up to 150 kg (330 lbs) and includes a handle for added security.

- Solid handle for added security
- Fits bathtub inner widths between 16" and 26<sup>3</sup>/<sub>4</sub>"
- Quick and easy tool-free installation
- Supports up to 150 kg (330 lb)



AMG770420  
AADL: B055



## CHROME GRAB BARS



Durable steel grab bars with non-slip knurled chrome surface. Rust resistant. Complete with screws and instructions. Clearance from wall: 1 3/4" (4,4 cm)

- Knurled grip for added security
- Rust-resistant chrome finish
- Stainless steel hardware included

AMG785100	12"
AMG785102	16"
AMG785104	18"
AMG785106	24"
AMG785108	32"

AADL: B116

*\*Also available without rotating flange*



## CHROME BATHTUB HANDRAIL



Features chrome-plated steel construction and a screw-tight clamp mechanism with rubber padded contact points. Knurled tips to reduce hand slippage.

AMG785300  
AADL: B180

## TRI-GRIP GRAB BAR



AMG785304  
AADL: B180

## MULTI-ADJUST TUB RAIL



The new AquaSense® multi-adjust bath safety rail allows for both parallel and perpendicular orientation to the bathtub wall, and features four (4) different height settings. It provides a feeling of security, stability and comfort when entering or exiting the bathtub.

AMG785340  
AADL: B180 Upgrade

## BATH SAFETY RAIL



This AquaSense® Safety Rail provides safety and comfort while getting in and out of the bathtub. Constructed with heavy gauge polyethylene construction, it wipes clean easily and installs securely without the use of tools.

AMG785370  
AADL: B180



## SUCTION GRAB BARS



Provides sturdy support and confidence in and around the tub and shower.  
Helps prevent dangerous slips and falls.  
Grip area has a ribbed surface for extra help with grip.  
Installs instantly.

DRIVE-INC1268 12"  
UNI02181196 20"

## BATH SAFETY RAIL



This AquaSense® Safety Rail provides safety and comfort while getting in and out of the bathtub. Designed to fit a wide variety of tubs, it offers two gripping surfaces. It installs easily and securely without the use of tools and is manufactured from rust-resistant steel. It features an easy to clean, white vinyl coated surface, and non-marring rubber pads on all

AMG785350  
AADL: B180

## BATH MAT WITH MASSAGE ZONES

Full-length bath mat with invigorating massage zones

- Safety

A full-length bath mat with over 400 suction cups to provide an active, secure grip that conforms to your tub, making it safer for everyone

- Hygiene

Antimicrobial formulation helps control mildew and bacteria. Sturdy construction allows for convenient machine washing.

- Versatility

Innovative design accommodates specialized bath seating without compromising stability.



AMG785500

## SHOWER SPRAY



AMG770980

- 3 spray settings can be activated by an easy turn spray dial
- Choices of spray: narrow, wide or standard
- Sturdy on/off pause switch provides easy water flow control at your finger tip
- Ultra-long 80", stainless steel, tangle-free hose allows you to shower comfortably while sitting
- Tool-free setup
- Corrosion resistant
- Easy-to-turn 3-spray selector
- Shower hook included

## REACHERS



AMG704175  
AMG704170  
AMG704185  
AMG704178

20" REGULAR  
26.5" STANDARD  
30" W/SUCTION CUPS - LONG  
32.5" EXTRA LONG

## LONG SHOE HORN 24"

PARSON-16D017

A one piece white epoxy coated steel shoe horn 24" (61 cm) long with a handy hang up hole.



## STOCKING AID

PARSON-16D005

A flexible sock aid with nylon webbing straps. Curved shape holds stocking on the aid while pulling it up your leg.



## RUBBER INVALID RING

AMG745174 16"  
AMG745176 18"

Convenient for home, car or workplace, this invalid ring allows you to sit comfortably for long periods of time. Inflation level can be adjusted to maximize comfort and the open center helps distribute weight evenly. Made from heavy duty rubber with a metal air release valve.



## CUTLERY WITH BUILT UP HANDLE

Cutlery with a grip designed to give independence at meal time. These utensils allow maximum user control with minimum effort. Useful in both institutions and at home. The comfort grip has a soft vinyl coating with finger indentations for a secure grip. Of special note is the Parsons knife blade which has been designed to be useful, both as a regular dinner knife and as a rocker knife. It has a curved cutting edge making it a suitable knife for the one handed person. Suited for persons with poor grasp with a 7/8" (22mm) diameter grip. Light weight, 85 grams (3oz) each. Dishwasher Safe. Latex Free. Blue Grips: All comfort grips are available with blue grips for the visually impaired, when ordering add "BL" to the end of the part number.

PARSON-16T001	KNIFE
PARSON-16T002	FORK
PARSON-16T003	TEASPOON
PARSON-16T004	SOUP SPOON



## LEFT & RIGHT HANDED CUTLERY



Cutlery with a grip designed to give independence at meal time. These utensils allow maximum user control with minimum effort. Useful in both institutions and at home. The comfort grip has a soft vinyl coating with finger indentations for a secure grip. Of special note is the Parsons knife blade which has been designed to be useful, both as a regular dinner knife and as a rocker knife. It has a curved cutting edge making it a suitable knife for the one handed person. Suited for persons with poor grasp with a 7/8" (22mm) diameter grip. Light weight, 85 grams (3oz) each. Dishwasher Safe. Latex Free. Blue Grips: All comfort grips are available with blue grips for the visually impaired, when ordering add "BL" to the end of the part number.

PARSON-16T005	LEFT HAND FORK	PARSON-16T006	RIGHT HAND FORK
PARSON-16T007	LEFT HAND SPOON	PARSON-16T008	RIGHT HAND SPOON
PARSON-16T025	LEFT HAND TEASPOON	PARSON-16T026	RIGHT HAND TEASPOON

## PLASTIC FOOD BUMPER



Food bumpers adapt a standard plate into a plate for self feeding similar to an inner lip plate or a scoop plate. The vertical wall provides an area to push food onto your utensil for scooping and aids in keeping the food on the plate. Fits 8" - 11" (20-28cm) plates. This sanitary snap-on food bumper is molded in one piece. Wash by hand. Microwaveable.

PARSON-16T070

WHITE

PARSON-16T070BL

BLUE

## PARSONS SCOOP PLATE

This 9" (23cm) scoop plate promotes independence at meal time. It allows the user to push food against the plate's high back wall and onto their cutlery. The entry side is low for ease of use. The wide outer rim is useful to grasp the plate while scooping food onto cutlery. The plate has a skid resistant base which is helpful for the individual with one hand or for anyone with poor coordination. Breakproof plastic. Microwaveable and top rack dishwasher safe.

ifferent colour plates are ideal for differentiating food groups, such as meat, vegetarian, etc. Darker plates are useful for the visually impaired, especially if used on a light coloured tablecloth.



PARSON-16T115



PARSON-16T115BL

## INNER LIP PLATE



The original Inner Lip Plate. The inner lip holds food on the plate while the user brings the fork or spoon to the edge of the plate. Ideally suited for one handed persons. Break proof and microwavable. 9" (23cm) diameter, 3/4" (19mm) deep.

PARSON-16T118  
PARSON-16T118BL

SANDSTONE  
BLUE

## SECTION PLATE WITH LID

This sectioned plate has three compartments to keep food separated. Each compartment is about 1.5" (3.8cm) deep. The two smaller sections hold up to 6.75 oz (200ml), while the larger one has a 12.6 oz (375ml) capacity. The high walls of each compartment can be used to help push food onto forks and spoons. This stackable dish comes with a clear plastic cover for easier transporting and storing of food. Diameter measures 8.25" (20cm). Plate and lid are dishwasher and microwave safe.



PARSON-16T167

## CONVALESCENT NO - SPILL CUP



The convalescent cup has a lid with a mouth piece, which has a hole large enough to accommodate a straw. Opposite the mouth piece there is a small air inlet hole in the lid. By putting a finger over the air hole you can regulate the flow of the liquid. 8 oz. (237ml) capacity. Ideal for drinking while lying in bed.

PARSON-16T143



## PARSONS MUG WITH LID

This crystal clear plastic mug has double handles to make drinking safer and easier for people with a weak grasp or tremor. The handles will accommodate four adult fingers and the lid has a drinking spout to prevent spills. The transparent design instantly shows the level of contents. Wide base to resist tipping. Dishwasher safe to 194 F (90 C). Microwaveable. 10 oz (295ml) capacity.



PARSON-16T123

PARSON-16T125 NO-SLIP

## KENNEDY CUP (SPILLPROOF CUP)



Sturdy screw on lid with internal nozzle eliminates accidental spills. The Kennedy Cup will not spill, so it can be used in a reclining or sitting position. Easy to grip handle is perfect for patients of all ages, left or right handed. Lightweight and durable, can be picked up by even the weakest hands. The Kennedy Cup holds 7oz of hot or cold liquid. Accepts any standard straw. Designed especially for rehabilitation and as an accessory for daily living. Helps restore independence and dignity. Institutional dishwasher safe. Blue colour. Straw not included. Not recommended for the microwave.

PARSON-16T149

## HANDYCUP

Handycup was developed in Denmark for patients with coordination difficulties. Thanks to its oblique shape and solid handles drinking from the Handycup is an easy task. Handycup is produced in FDA-approved transparent San to make it easy to see how much liquid is in the cup. This cup also has ml indication so a caregiver can track the amount of liquid ingested. Handycup is household dishwasher and microwave oven safe.



PARSON-16T133A



## TERRY CLOTHING PROTECTOR

This terry cloth bib is 17" X 30" (43 X 76 cm). Attaches with tie strings.

PARSON-16T082



## CLOTHING PROTECTOR WITH SNAPS



A colourful bib to brighten dinner time. 18" wide x 35" long (46cm x 89cm). Poly/cotton front with waterproof backing. Neck opening closes with snap buttons. Comes in a variety of colours and patterns.

PARSON-16T083-2  
PARSON-16T083

WITH TIES  
WITH SNAPS

## COMFORT HAIR WASH TRAY

This plastic tray mounts on the shoulders and is held in place via a 35" (89 cm) strap which hangs in front of the user for them to hold. It can be used in any home sink, bathroom or kitchen, or with any shampoo basin. The 5" (13cm) wide curved neck rest provides a comfortable support. Water poured on the hair drains into the sink or basin. The tray is 17 3/4" (45cm) long x 13" (33cm) wide with 1 1/2" (4cm) high sides.



PARSON-16B096A

## COMFORT HAIR WASH BASIN



PARSON-16B098A

Inflatable shampoo basin makes it possible to wash a person's hair without removing them from bed. The patient's neck and shoulders are cushioned in the inflatable rings. Handy for home or institutional use. Inflates by mouth or by air pump. Made of heavy duty vinyl. Complete with 4" drain hose, plug and instructions. Inflates to 28" (71cm) wide x 24" (61cm) long x 6" (15cm) deep.

## LIGHTWEIGHT TRANSPORT CHAIRS



- Anatomically correct seat-to-armrest height for a naturally erect posture
- Detachable full flip-back armrests for effortless transfer
- Fold-down back for easy transportation and storage
- Comfort-Grip 8" wheels, for in- and outdoor use and easy maintenance
- Removable push-in swing-away footrests with tool-free adjustment of length for greater comfort.
- Easy-adjust footrests with heel-loops for added support and comfort. Heel-loops and seatbelt included as standard features for increased safety.
- Dual locking rear wheels for safe, secure transfers.
- Chip-resistant finishing and no sharp edges on hardware .
- Lightweight and compact weighing only 9.3 kg.
- Weight capacity: 300 lb (136 kg).
- Each chair is assigned a unique serial number located on the bottom crossbar for effective after-sale product support.
- Comes complete with bilingual instruction manual and a tool package.

AMG700846

BLACK

AMG700848

PLAID

## WHEELCHAIR WITH FULL DETACHABLE ARMS AND SWING-AWAY FOOT RESTS



This affordable wheelchair features dual axle positions that allow for varying seat-to-floor heights, a heavy-duty chrome-plated steel frame, solid snap-on tires, adjustable brakes, sealed bearing assemblies, padded armrests, and flame-retardant black naugahyde upholstery. Parts as well as accessories are stocked in Canada for fast and easy availability.

- Institutional tested durability & proven quality
- Proprietary contact points are easy to clean & disinfect
- Dual axle positions allow varying seat heights
- Parts stocked in Canada for fast & easy availability
- Maximum weight capacity: 250 lbs (113,4 Kg)
- Latex free

AMG700522

16"

AMG700622

18"

## CHALLENGER



- The Challenger series is our newest and most versatile walker. With a seat height of 18", the Challenger Mini is a small and compact yet strong walker designed especially for the Japanese market.
- The Challenger Grande comes with a larger seat and a 21" seat height and has a storage bag under the seat.
- All 4 models have 8" wheels and a flexible back support.
- Both the Wide and Low Wide models have a 400 lbs capacity and only weight 16lbs!

AADL: C061

## EVOLUTION



- The Evolution Series is the therapists choice for individuals of higher weight
- The Evolution has a weight bearing capability of up to 400lbs
- The rigid steel frame gives this walker unmatched stability for users suffering from neurological conditions such as Multiple Sclerosis, Parkinson's and stroke conditions.
- The rigid frame allows the addition of the specialized Arm Trough System, for those individuals unable to use a standard walker handle.
- Available with standard rigid back support or flexible back straps
- Wide models have a larger seating area with a bigger soft seat cushion for added comfort.
- Soft seat cushion in every model
- Plastic tray or basket installed tote bag included
- Now with stronger and wider wheels 40mm wide for ease of use on gravel roads, wet grass, deep carpet, etc.

AADL: C064

## HORIZON



- Similar design to the best selling Piper, but with with a steel frame that accommodates heavier weight individuals up to 400 lbs.
- New 40mm wide spoke wheels with long lasting PU rubber.
- Ideal for users with some neurological conditions that require a very solid walking platform.
- Wide models have a larger seating area with a bigger soft cushion for added comfort.

AADL: C064

## PIPER



- Strongest most stable light-weight walker on the market with fully welded 1" aircraft aluminum frame for superb stability
- One hand folding system easiest walker to handle by both user & caregiver
- Short wheel base allows for smaller turning radius, ideal for smaller living spaces
- Narrow frame allows user to go through small door openings like bathrooms, etc
- New safety accessories with reflective tapes visible from every angle
- Redesigned front handles with convenient hook for hanging purses, shopping bags, etc
- Easy locking brake system, ideal for arthritis patients
- New 40mm wide wheels with long lasting PU rubber
- New mesh baskets

AADL: C061

## SPRITE



- Economical light weight walker
- 6" wheels models now come with wider wheels, ideal for indoor use
- 8" wheel model for comfortable outdoor use, (Sprite Grande)
- Rigid contoured back support for secure seating
- New rear fork system on regular Sprite & Sprite Mini
- New Slow Down Brake available as an option on Sprite & Mini
- New safety reflective attachments on all Sprites

AADL: C061

## XPRESSO



- 'X' style center folding mechanism for compact storage & transportation
- Collapsible basket, walker can be folded with the basket still on
- Cable-free brakes with easy to adjust handle height
- New slow-down brake can be added, quickly & easily
- New wide polyurethane tires, long wearing
- Built in padded, waterproof seat cushion
- Comfortable form fitting rubber handgrips
- Available in 10 different flavors: Red, light blue, kiwi green, black, pearl white, orange, champagne, cappuccino, yellow, silver

**\*\*XPRESSO LITE also available**

AADL: C061



## EXCURSION



### Available Heights & Colours Include:


AMG700-926	X18 Hemi Height	Pearl Black
AMG700-927	X18 Hemi Height	Cranberry
AMG700-920	X20 Lightweight	Pearl Black
AMG700-921	X20 Lightweight	Pacific Blue
AMG700-922	X20 Lightweight	Cranberry
AMG700-923	X20 Lightweight	Lime
AMG700-924	X20 Lightweight	Red Rocket
AMG700-917	X23 Lightweight	Pacific Blue
AMG700-928	XWD Lightweight	Pacific Blue
AMG700-929	XWD Lightweight	Pearl Black

The lightest side-fold rollator around!

At just 15.4 lb (6.98 kg) the Airgo® eXcursion™ weighs 20% less than the competition while still providing peace of mind with industry leading stability.

- Easy to fold by just a pull of the strap
- Ultra-light frame makes lifting and storage simple
- Soft basket doesn't need to be removed when folding
- Ergo-form handgrips fit comfortably and place your hand at a natural angle for better brake access
- Padded seat and backrest for added comfort
- Locks in open position for added safety
- Reflective trim for added visibility in the dark
- Easy to activate hand brakes
- 8" (20 cm) Comfort-Grip non-marring wheels are ideal for indoor & outdoor use
- Easy set-up
- Net weight (without accessories): 15.4 lbs (6.98 kg)
- Maximum weight capacity: 300 lbs (136 kg)

AADL: C064



## AIRGO FUSION SIDE-FOLDING ROLLATOR & TRANSPORT CHAIR

The Airgo Fusion was designed to allow you to remain independent and live your life to the fullest. Its lightweight design and unparalleled stability provide the confidence to walk safely again. And if you ever find that you've ventured out too far, the Airgo Fusion converts easily to a transport chair so that a loved one can help bring you back home in comfort and style. Find out more about why so many people love their Airgo Fusion!



- New 18" (45.7 cm) seat height
- Lightweight side-folding rollator and transport chair in one design
- Folds sideways easily by just a pull of the strap!
- Patent pending Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest yours arms in either mode
- Smooth bearings and an improved center of gravity for easy pushing
- Unique ergonomic handgrips fit naturally for easier brake access
- Features our new senior friendly tool-free brake adjustment knob
- Lightweight frame (17.1 lb net) designed for superior stability
- Reflective soft basket can be removed and is machine washable
- 8" (20 cm) Comfort-Grip wheels are ideal for indoor & outdoor use
- Net weight (without accessories): 17.1 lb (7.76 kg)
- Weight capacity of 300 lb in both modes

AADL: C061

AMG700-938  
AMG700-935

Hemi Height Side-Folding Rollator & Transport Chair  
Side-Folding Rollator & Transport Chair

Cranberry  
Pacific Blue



## AIRGO ULTRA-LIGHT 6 ROLLATOR

The lightest rollator available: a breeze to lift up and take anywhere!

- Durable, ultra-light frame enables you to lift the rollator easily and take it anywhere.
- Ergo-form handgrips fit more comfortably and place your hand at a natural angle for comfort and better brake access.
- Multi-set height technology: Six height adjustment.
- Padded seat and backrest.
- Soft basket doesn't need to be removed when folding.
- Easy-to-activate brakes.
- Locks in open position for added safety and in folded position for convenience.
- 6" non-marring wheels are ideal for indoor use.
- Reflective piping for added visibility in the dark.
- The handgrips come with built-in reflectors for more safety in the dark.
- Tool-free set-up.
- Net weight (without accessories): only 4.7 kg (10.5 lb)
- Patents: 7828305, 2587989, 200710103964.9
- Maximum weight capacity: 250 lb (113 kg).



AMG700950

## NEXUS WALKER



The neXus 2 walking aids is the perfect walker for the active user who wants to enjoy life. Complete with basket, the neXus 2 features an ultra-compact "X" folding system that allows it to stand independently when folded. Its lightweight frame makes getting the neXus 2 out of the car so much easier and its solid seat structure includes a contoured, cushioned seat that makes a quick rest during a busy day extremely comfortable.

### Features

- Push-to-Lock braking system with arthritic ball for ease of use
- All the same great features found on the neXus 3
- Large mesh basket
- Stands independently when folded

\*\*Nexus 3 also available with cable-free braking system.

## CLEVER-LITE LS WALKER WITH SEAT, PUSH DOWN BRAKES



DRIVE-10233

- Rear Glide Push down brakes.
- Sturdy 1" diameter anodized, extruded, aluminum construction.
- Comes with flip up seat with built in carry handle.
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame.
- Easily folds with dual lever side paddle release.
- Soft, flexible backrest provides comfort and security while seated.
- 5" fixed front wheels.

Seat Height 19.5" - 23.5"

32" Handle Height

Product Weight 11 lbs

## CLEVER-LITE WALKER WITH SEAT AND LOOP LOCKS

- Front wheels can be set in either "swivel" or "fixed" position.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Sturdy 1" diameter anodized, extruded, aluminum construction.
- Comes with flip up seat with built in carry handle.
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame.
- Easily folds with dual lever side paddle release.
- Soft, flexible backrest provides comfort and stability while seated.
- Attractive blue finish.
- Limited Lifetime Warranty.



DRIVE-10230

5" Casters

Seat Height 19.5" - 23.5"

36" Handle Height

Product Weight 13 lbs



DRIVE-10243

8" Casters

Seat Height 23" - 26"

39" Handle Height

Product Weight 16 lbs

## ADJUSTABLE ALUMINUM FOLDING WALKERS



AMG770120  
AADL: C030

This economical and lightweight walker features 1" diameter anodized aluminum tubing, geometric design, and angled front legs. The walker front bar also acts as a support rail for use over a commode bowl. The one-button folding mechanism makes the walker easy to open and close.

Available in adult or small adult sizes, silver or bronze colors.

Maximum weight capacity: 300 lbs (136 Kg)



AMG770160  
AADL: C030

## WALKER ACCESSORIES

Attach to the bottom of walker legs in order to help the walker glide.



DRIVE-RTL10119  
Tennis Ball Glides



DRIVE100014  
Skis



AMG780010  
Skis



MEDLG07900  
Glide Caps



AMG780000  
Skis



AMG780118  
Glide Caps

## FOREARM CRUTCHES

Constructed of heavy-gauge aluminum tubing, these forearm crutches feature vinyl coated, front opening forearm cuffs, steel handles with plastic grips, and x-large rubber crutch tips. Double height adjustments allow for independent adjustment of both hand grip and cuff height.

AMG720855	28" to 37" (71.1 cm to 94 cm)
AMG720865	36" to 45" (91.4 cm to 114.3 cm)
AMG720875	20" to 29" (50.8 cm to 73.6 cm)

AADL: C186



## FOREARM CRUTCHES EURO STYLE

Constructed from lightweight anodized aluminum tubing, the Euro-Style Forearm Crutches feature a one-piece molded plastic cuff and hand grip assembly, as well as a secure fit adjustment pin assembly.

Adjustable from 28" to 38" (71.1 cm x 96.5 cm), in one inch increments.

AMG720850



## ADJUSTABLE ALUMINUM CRUTCHES



AMG720790	YOUTH 4'6" - 5'2"
AMG720785	MEDIUM 5'2" - 5'10"
AMG720780	TALL 5'10" - 6'6"

## REPLACEMENT PARTS FOR CRUTCHES



AMG725150  
AMG725110  
AMG725160

HAND GRIPS  
CRUTCH TIPS  
UNDERARM PADS

## COMFORT-PLUS CANE

- 12 individual height settings offer an unparalleled range of adjustments that fits most people 5' - 6'5".
- Cane height adjusts from 28" to 39".

AMG730442	BLUE
AMG730444	BLACK
AMG730446	CHARCOAL
AMG730448	OLIVE
AMG730450	BURGUNDY
AMG430452	CORAL



## COMFORT-PLUS ADJUSTABLE FOLDING CANES



AMG730490	BLUE
AMG730491	BLACK
AMG730492	CHARCOAL
AMG730493	OLIVE
AMG730494	BURGUNDY
AMG730495	CORAL

## CANE REPLACEMENT TIPS



AMG735001

## COMFORT-PLUS ADJUSTABLE QUAD



- K-shaped quad base provides greater stability
- Can be easily rotated to allow for left or right hand use without getting in the way.
- 11 individual height settings offer a range of adjustments that fits most people 5' - 6'5".
- Cane height adjusts from 28" to 39".

AMG770842	BLUE	AADL: C086
AMG770844	BLACK	
AMG770846	CHARCOAL	

## HUGO QUADPOD



The new & innovative Hugo QuadPod ultra stable cane tip is 54% lighter & 80% more compact than a large base quad cane. It's unique design recenters itself for maximum stability and allows your cane to stand independently so you can go anywhere you want to go!

AMG781030



## UPLIFT SEAT

The Uplift Seat Assist with Memory Foam and Waterproof cover is a self-powered lifting cushion that converts an armchair or wheelchair into an automatic lifting chair. Uplift safely provides the extra lift needed to help you get in and out of chairs on your own. Allowing you to remain independent and active, without having to rely on others. No switches, electric wires, latches or special installation is necessary. The Lift Actuator lets you choose one of six weight settings 80 to 230+ lb. (37 to 104kg) that's best for you. A heavier weight limit is available on the 16H019A Uplift Plus ranging from 200 to 350 lb. (90kg to 160kg). When you sit down on Uplift's comfortable cushion, your downward motion is "hydraulically dampened" so you settle gently into the chair. No more "falling" into the chair due to muscle weakness when you use the Uplift. Uplift has a very comfortable cushion and adds only one inch to the height of the chair when you're seated. When you want to get up from your chair, Uplift automatically provides the extra lift needed to get you to a standing position. Uplift lifts up to 80% of your weight 10" (25cm) high! The pneumatic Lift Actuator gives a smooth counterforce, following your body's natural motion, while it supplies the needed boost to help you stand up. No more struggling to get up from a chair. uplift provides the lifting strength you need.

PARSON-16H018A



## UPEASY SEAT

UPEASY SET This lifting cushion is very similar to the Uplift seat, except this one uses regular poly foam and a polycotton cover.

WEIGHT CAPACITY: 95/220 LBS

PARSON-16H015A



## ALARM SETS



**1 YEAR WARRANTY  
ON ALARMS!**

- Available for beds, under mattress, wheelchairs, chairs & toilet seat pad.
- Some models compatible with Nurse Calls
- Available as code, silent or voice alarm

## MAT SETS



### MATS:

- Alarm sounds when resident steps on safety mat
- Heavy Duty PVC Safety Mat is slip resistant on both sides
- Latex-free, mildew resistant & flame resistant
- 2 sizes to pick from

### ALARMS:

- Always on for resident safety
  - 80 dB
  - 8' connecting cord
  - Tamper-resistant covers conceal all settings
  - Low battery indicator sound & light continues to alert until batteries are changed
  - Light bar easily identifies activated alarm
  - Holder & batteries included • Screw driver
- Additional SMA-120 Alarm Features:
- Optional AC power supply saves battery life
  - 3 decibel levels 80, 90, 120 dB

## SEAT BELT ALARM



- 31"- 54" adjustable seat belt has quickrelease, non-restraint activation
- 5' connecting cord
- Attaches easily using two existing screws on back of wheelchair
- No tools or mounting materials needed to attach alarm to wheelchair frame
- 85dB alarm with auto reset is always on for resident safety
- The **ONLY** Seat Belt Alarm with metal connecting pin...no breakable plastic parts!

## TRI-FOLD BEDSIDE FLOOR MAT



PSC-SBSM-1

30" x 66" x 2"

- Latex-free
- Tri-fold for easy storage
- Mildew-resistant
- High-density, high-impact foam
- Flame retardant
- Slip-resistant on both sides
- Easy to clean

## BEVELED EDGE SAFETY MAT

- Low one-inch profile makes it easy for residents to step on and off mat
- Beveled edges minimizes resident's potential to trip on mat
- Beveled edge allows wheelchairs and medical equipment to easily roll over mat
- Slip-resistant on both sides
- Latex-free
- Mildew-resistant
- Flame retardant
- High-impact absorption EVA foam



PSC-SBEV-1

38" x 70" x 1"

## ROLL-ON SAFETY MAT



SKIL911548

68" X 26" X 1 1/2"

- Safer for residents who might stumble on standard floor mats
- High-density foam mat absorbs the impact of a fall from bed
- Durable vinyl cover wipes clean for easy care
- Slip-resistant bottom keeps mat safely in place
- Folds in half for easy storage
- Includes strap that secures mat to footboard or allows it to be hung from wall hook.

## SAFESIDE FALL MAT

- Shock-absorbing bonded foam cushions fall impact
- Step-on stability prevents tumbles when resident walks on mat
- Provides safer footing –won't "bottom-out" when walked on
- Slip-resistant bottom keeps mat securely in place
- Low-profile/tapered edge design reduces resident or staff stumbling over mat
- Heat-sealed vinyl cover protects against moisture and bacteria—important for infection control
- Available 24-inch- and 36-inch-wide bi-fold and non-folding models



SKIL911577

68" X 36" X 1"

## SECURESLEEVES™ FOR ARMS

- Protects sensitive skin from tears and abrasions
- Unique blend of Chitosan rayon, nylon, cotton & Spandex
- Anti-microbial properties protect against odor causing bacteria often accompanying deep wound and skin problems
- Also protects against the growth of bacteria, molds and fungi
- Can also be used to cover I.V. and wound sites
- Machine washable and Latex-free

Available in sizes Small-Extra Large



## SECURESLEEVES™ FOR LEGS

- Protects sensitive skin from tears and abrasions
- Unique blend of Chitosan rayon, nylon, cotton & Spandex
- Anti-microbial properties protect against odor causing bacteria often accompanying deep wound and skin problems
- Also protects against the growth of bacteria, molds and fungi
- Can also be used to cover wound sites
- Machine washable and Latex-free

**Available in sizes:** Large/Extra Large  
Small/Medium XX Large



## SECURESOCKS™ FOR FALL MANAGEMENT



- All around non-skid pattern
- One Size Fits Most
- Length - 17"
- All around non-skid pattern to assist in preventing patient/resident falls
- Ultra soft 100% cotton for maximum comfort and warmth
- Longer to help prevent them from coming off a patient/resident's foot

## ACTI-TRED SOCKS




- Double printed socks to help prevent falls
- Various sizes available to fit almost everyone.

MEDAC99942	AQUA	TODDLER
MEDAC99943	CRANBERRY	SMALL/YOUTH
MEDAC99944R	BLUE	MEDIUM
MEDAC99945	BEIGE	LARGE
MEDAC99946	NAVY BLUE	XLARGE
MEDAC99947	GRAY	XXLARGE
MEDAC99948	BLUE	BARIATRIC
MEDAC99949Y	YELLOW	UNIVERSAL

48 PRS/CS

## GAIT BELT WITH METAL BUCKLE

- Yellow Gait Belt for fall risk
- Designed to aid in lifting, walking or transferring residents
- Choice of EZ release YKK® plastic buckle or metal buckle
- 60" or 80" EZ Release YKK® Plastic Buckle – weight capacity: 280 lbs.
- 60" or 80" Metal Buckle – weight capacity: 340 lbs.
- 100% Cotton
- Can be laundered and dried
- Includes loop to hold excess length
- One year warranty



PSC-SGBM-80Y 80"

## PADDED GATE BELT

These padded transfer belts will assist patient transfers with comfort, dignity and safety for both care giver and patient. The padded areas are lined with Bucktex for soft, breathable comfort. The pads are 5.5" wide (13 cm). Four handle loops are spaced around the belts for convenient, safe grip points. Colour coded for quick identification in the rehab department. Available with plastic side release (SR) buckle or the extra safe auto buckle.



PARSON-17S140S	SMALL
PARSON-17S140M	MEDIUM
PARSON-17S140L	LARGE



## REMOVABLE TAILBONE PAD, SEWN IN HIP PADS



PSC-SHP-XS Extra Small: Waist (24"-27") Hips (31"-34")  
 PSC-SHP-S Small: Waist (28"-31") Hips (35"-38")  
 PSC-SHP-M Medium: Waist (32"-35") Hips (39"-42")  
 PSC-SHP-L Large: Waist (36"-38") Hips (43"-45")  
 PSC-SHP-XL Extra-Large: Waist (39"-42") Hips (46"-49")  
 PSC-SHP-XXL Extra-Extra-Large: Waist (42"-46") Hips

- Waterproof high impact foam pads
- Memory foam pads contour to body
- Worn as underwear or outside clothing
- Breathable 94% Cotton / 6% Spandex
- Latex free

## REMOVABLE HIP & TAILBONE



PSC-SHP-RP-XS Extra Small: Waist (24"-27") Hips (31"-34")  
 PSC-SHP-RP-SW Small: Waist (28"-31") Hips (35"-38")  
 PSC-SHP-RP-M Medium: Waist (32"-35") Hips (39"-42")  
 PSC-SHP-RP-L Large: Waist (36"-38") Hips (43"-45")  
 PSC-SHP-RP-XL Extra-Large: Waist (39"-42") Hips (46"-49")  
 PSC-SHP-RP-XXL Extra-Extra-Large: Waist (42"-46") Hips (49"-53")

## OPEN FRONT



PSC-SOHP-S/M Small/Medium 28"-34"  
 PSC-SOHP-L/XL Large/X-Large 34"-42"

## REPLACEMENT PADS



Available  
individually

PSC-SHP-PAD

## HIPSAVER QUICKCHANGE



HipSaver QuickChange is a functional garment designed for people living within a long term or acute environment who are dependent on staff to change their incontinence briefs. The drop down crotch panel assists staff with quick changes of incontinence products.

## HIPSAVER WRAP & SNAP



HipSaver Wrap&Snap is the another innovation from HipSaver developed in response to requests from nursing staff for a hip protector that is both functional and easy to manage when a resident needs a complete change of clothing or because balance issues make changing a prone resident easier. To change a resident, simply unsnap the drop down crotch panel followed by the front panel. The snaps are reinforced because caregivers told us they wanted a quick, one-handed operation. HipSaver Wrap&Snap is made from a comfortable polycotton fabric in an attractive light coffee colour. The hip pads are 1/2" thick and measure approximately 7 1/2" in circumference.

## HIPSAVER OPEN-BOTTOM



HipSaver Open-Bottom is the latest innovation from HipSaver. Requested by nursing staff who saw a need for residents to have a garment that permits easy toileting, the open crotch design is perfect for residents with urge incontinence or challenges in getting clothing on and off. The HipSaver Open-Bottom design retains the discretion of the SlimFit models making it perfect for community dwelling residents as well as those residing in assisted living facilities.

The pant is made from comfortable polycotton fabric and contains no latex or skin irritants. The hip pads are oval in shape measuring roughly 7.5" long by 6.5" wide and 1/2" thick tapering to 1/16" at the bevelled edge. They are heat sealed into a high temperature tolerant waterproof pouch. The hip pads are odourless and contain no skin irritants or latex. HipSavers are comfortable enough to wear 24 hours a day unlike hard shell hip protectors.



## HIPSAVER SLIMFIT



HipSaver SlimFit is a fashionable garment designed for people living independently in the community. So discreet, nobody will know that it's being worn under most clothing. The pant is made from comfortable polycotton fabric in an attractive light coffee colour and contains no latex or skin irritants. The SlimFit model will accommodate incontinence liners comfortably. The hip pads are oval in shape measuring roughly 7.5" long by 6.5" wide and 1/2" thick tapering to 1/16" at the bevelled edge. They are heat sealed into a high temperature tolerant waterproof pouch. The hip pads are odourless and contain no skin irritants or latex. HipSavers are comfortable enough to wear 24 hours a day unlike hard shell hip protectors.

## HIPSAVER NURSING HOME



HipSaver Nursing Home is a utilitarian garment designed for people living in a long term or home care environment. Wide leg openings ensure comfort when worn over top of incontinence briefs. The pant is made from comfortable polycotton fabric in an attractive light coffee colour and contains no latex or skin irritants. The soft hip pads are 1/2" thick and measure approximately 7" x 7". They are heat sealed into a high temperature tolerant waterproof pouch. The hip pads are odourless and contain no skin irritants or latex. HipSavers are comfortable enough to wear 24 hours a day unlike hard shell hip protectors.

## HIPSAVER SOFT SWEATS



Perfect for both women and men, this all soft hip protector has the same proven protection as our other models but in a comfortable, casual style. Features a loose ankle hem - no binding elastic - and six hip sizes and two lengths for a perfect fit. Machine wash and dry. HipSaver SoftSweats are pre-shrunk to fit properly.

Now available in Navy, Grey, Khaki and Camouflage

## SIZING CHART FOR HIPSAVER BRAND

Hips:	Inches	Size	Hips:	Centimeters
28-32	XS		70-81	
32-36	S		81-92	
36-40	M		92-102	
40-45	L		102-114	
45-51	XL		114-129	
51-57	XXL		129-145	

## WAFFLE® SEAT CUSHION

The Original and Standard WAFFLE Seat Cushions are recommended for patients up to 300 pounds. The Pediatric Seat Cushion is recommended for individuals under 100 pounds. All cushions may be used for comfort and treating and preventing pressure ulcers. Each cushion features a unique low profile design to provide safety and convenience for getting in and out of chairs, and the air venting holes allow heat and moisture to easily dissipate for a higher level of comfort. Most importantly, our seating cushions utilize the same static air technology of our top-selling mattress overlays to aid in the prevention and treatment of pressure ulcers up through Stage IV.

Most WAFFLE Cushions are shipped inflated. However, you may need to adjust the inflation from time-to-time. In order to do so, we recommend the WAFFLE L.A.D.(Lower Air Delivery) Pump to inflate your WAFFLE Brand Cushion.

A fluid resistant cover is available with the Original and Standard Adult WAFFLE Seat Cushion (200WCIC and 218WCIC). SINGLE PATIENT USE - 3 MONTH WARRANTY

EHOB214WCI	14 X 14 X 1"	PRE-INFLATED (PEDIATRIC)
EHOB200WCI	17 X 17 X 1.5"	PRE-INFLATED
EHOB218WCI	19 X 19 X 1.5"	PRE-INFLATED



## WAFFLE® BARIATRIC SEAT CUSHION



EHOB240WCI      22 X 28 X 2"  
PRE-INFLATED

The WAFFLE Bariatric Cushion is designed with a diamond-shaped pattern to evenly distribute an individual's weight and minimize material stress and is capable of supporting patients up to 700 pounds. The Bariatric Cushion's low air profile reduces pressure and helps prevent and treat pressure ulcers in most seating applications. It is also recommended for pain management and comfort. It is easy to clean, lightweight and shipped inflated ready for use. SINGLE PATIENT USE

## WAFFLE® M.A.D. PUMP

The WAFFLE M.A.D. Pump provides non-electrical, exact inflation for all WAFFLE Brand Mattresses and the WAFFLE Chair Pad(207GDGP) only. It ensures therapeutic effectiveness and maximizes the medically-proven benefits of static air support.

EHOB5400-MADX-000



## WAFFLE® CHAIR PAD WITH M.A.D. PUMP

Perhaps the most dangerous of all seating applications are the wooden day chairs. Our WAFFLE® Chair Pad is a perfect fit to protect the patient's entire body. The long form of the pad contours and cradles the body and provides a layer of static air between all of the body's major pressure points and the chair. Its low profile design allows for easy entry and exit, and adjustable straps secure pad and provide proper fitting to the chair.

\*3 Month Warranty



EHOB207GDCP 72 X 21 X 3"

## WAFFLE® MULTI-CARE PAD

The WAFFLE Multi-Care Pad is a double cushion offering relief to the back and the seat area in a wheelchair. It can also be used laterally to cradle and stabilize the patient.

\*3 Month Warranty



EHOB201WP

36 X 15 X 2"

## WAFFLE® L.A.D. PUMP

The WAFFLE L.A.D. Pump provides non-electrical, exact inflation for all WAFFLE Brand Cushions, Foot WAFFLE Custom and Heel Elevator Custom only. It ensures therapeutic effectiveness and maximizes the medically-proven benefits of static air support.



## WAFFLE® ECONO EXTENDED CARE PLUS MATTRESS OVERLAY

The WAFFLE® Mattress Overlay is made to go right on top of your existing mattress. The air overlay cradles the body, providing protection and treatment for pressure ulcers. And it is also recommended for pain management. The WAFFLE® Mattress Overlay is easy to transport, cleans with soap and water and is extremely comfortable. This product features a patented design to eliminate 'bottoming out', and air vents to dissipate heat and moisture, and a low profile design for greater stability and patient safety. A WAFFLE® M.A.D. hand pump accompanies the overlay to ensure precise inflation.

\*6 Month Warranty



EHOB1056ECP 76 X 34 X 3"

## WAFFLE® ECONO-CARE PLUS OVERLAY WITH PUMP



The WAFFLE EconoCare is our most affordable solution for comfort and pressure ulcer prevention. Designed primarily for hospital use, the EconoCare features a 30-day warranty and is intended for short term use.

Flex-O-Seal Design.

\*30 Day Warranty

EHOB1025ECP 76 X 34 X 3"

## WAFFLE® MATTRESS EXPANSION CONTROL OVERLAY

Recommended for Hospitals -Three Month Warranty. Pressure Ulcer Prevention and Treatment Through Stage IV.

The WAFFLE® Mattress Overlay is made to go right on top of your existing mattress. The air overlay cradles the body, providing protection and treatment for pressure ulcers. And it is also recommended for pain management. The WAFFLE® Mattress Overlay is easy to transport, cleans with soap and water and is extremely comfortable. This product features a patented design to eliminate 'bottoming out', air venting holes to dissipate heat and moisture, and a low profile design for greater stability and patient safety. A WAFFLE® M.A.D. hand pump accompanies the overlay to ensure precise inflation.



EHOB1004ECP 75 X 30 X 3"

## WAFFLE® FOOT CUSHION - LINED

The Foot WAFFLE® Custom surrounds the lower limb in air and is designed to float the heel while sustaining the integrity of the foot, ankle, calf and Achilles area. The air-filled straps and foot pad offer pressure ulcer prevention, treatment and patient comfort. It can be used conveniently at home and on O.R. tables, stretchers, ICU beds, wheelchairs and general ward beds.

\*36 Month Warranty



EHOB503FW	SMALL	CALF 10" TO 11.5"
EHOB504FW	MEDIUM	CALF 12" TO 14"
EHOB505FW	LARGE	CALF 14.5" TO 18"

## WAFFLE® HEEL ELEVATOR

The Waffle Heel Elevator is used for prevention and healing through Stage IV pressure ulcers. It suspends the heel completely off the surface and reduces plantar flexion and foot drop. The cloth lining and air venting holes keep the patient cool and dry.



EHOB345CX      ONE SIZE FITS MOST PATIENTS

## WAFFLE FOOTHOLD®

A new positioning device for the lower extremity that promotes proper alignment of the ankle, knee and hip, to prevent and treat heel ulcers. The patented Hydro-Cushion floats and lifts the calf while unloading the heel when recumbent. Ideal for use in hospital, long term care and home care.


\*30 Day Warranty



EHOBFH10-SMXX	SMALL
EHOBFH10-MDXX	MEDIUM
EHOBFH10-LGXX	LARGE
EHOBFH10-XLXX	EXTRA LARGE



The HeelZup® Cloud Heel Suspension Cushion is a new concept for the prevention and treatment of heel pressure ulcers. Until now, conventional pillows or contour cut foam cushions were the main choices available for elevating the lower legs and suspending the heels to eliminate pressure. The HeelZup Cloud Heel Suspension Cushion was designed to utilize the advantages and eliminate the disadvantages of both pillows and foam cushions.



## HEELZUP® CLOUD SOLO

For single patient and home use:

Manufactured of high-density, anti-microbial, hypoallergenic VE foam and thermobonded polyester fiber. The standard width is 30". The Solo version is covered with a lightweight coated nylon fabric with a Velcro button enclosure. The cover can be wiped down or removed and laundered in cold water on a gentle cycle.

INT-CLOUD-SD	5'4" AND UNDER
INT-CLOUD-MD	5'5" TO 5'10"
INT-CLOUD-LD	5'11" AND UP



## HEELZUP® UP CLOUD REUSABLE

Manufactured of high-density, anti-microbial, hypoallergenic VE foam and thermobonded polyester fiber. The standard width is 30". A moisture proof, reusable coated four-way stretch fabric cover with No-Slip bottom is sewn onto the cushion. A waterfall flapped zipper enclosure prevents moisture ingress through the zipper. The cover can be cleaned with standard cleaning procedures, using properly diluted solutions.

INT-CLOUD-SR	5'4" AND UNDER
INT-CLOUD-MR	5'5" TO 5'10"
INT-CLOUD-LR	5'11" AND UP



## HEELZUP THERAPEUTIC HEEL ELEVATING CUSHION WITH POPLIZONE THIGH SUPPORT

- Elevates the heels while supporting the leg from the thigh to the ankle.
- Eliminates pressure on the Popliteal Fossa.
- Raised side bolsters prevent legs from falling off of the sides of the cushion.
- Stays neatly in place and does not require constant maintenance and adjustment.
- Coated nylon cover reduces friction and shear, while the No-Slip bottom material helps the cushion to stay neatly in place. No need for constant maintenance and adjustment.
- The PopliZone Thigh Support can be retrofitted to previously purchased HeelZup cushions that are not equipped with the PopliZone Thigh Support.
- Hinged seam between the PopliZone Thigh Support and HeelZup Cushion allows the entire system to articulate.



INT-TSH14R  
INT-TSH14RB

STANDARD  
BARIATRIC

## HEELZUP REUSABLE

Manufactured of high-density, anti-microbial, hypoallergenic foam. The calf bed is wide enough to support both legs while the raised bolsters prevent the legs from falling off of the sides of the cushion. A moisture proof, reusable, coated nylon fabric cover is sewn onto the cushion. The cover can be cleaned with standard cleaning procedures, using properly diluted solutions. Also available in extra large bariatric sizes. Carries a three year limited, non-prorated warranty!



INT-HZ14R



The Heelzup cushion comes compressed for easy shipping.

## HEEL FLOAT

- Floats the heel over an air cavity for zero pressure
- Provides air circulation to aid healing of existing ulcers
- Foam cushioned sides protect lateral and medial malleoli
- Water-based gel insert lowers pressure beneath ankle
- Cozy Cloth cover reduces skin shear and wicks away moisture
- Anti-foot-drop construction controls dorsal flexion
- Launderable, reusable, and guaranteed for six months



SKIL503034

SMALL 3" WIDE

SKIL503035

MEDIUM 4" WIDE

SKIL503036

LARGE 5" WIDE

## ULTRA-SOFT HEEL CUSHION



- Filled with fluffy spun fiber
- Covered with low-friction brushed polyester fabric to reduce skin shear
- Cover wicks away moisture to eliminate skin-damaging perspiration
- Comfortable, extra-wide cloth strap won't irritate skin
- Velcro® closure keeps heel cushion in place
- Reusable and launderable
- Meets California 117B flammability standards

SKIL503030

UNIVERSAL SIZE 1 PR

## THIN-LINE ECONO-GEL CUSHION

- Top and bottom layers of foam provide comfort
- Specially designed gel pod relieves pressure and increases stability
- Safety ties keep cushion in place
- Incontinence-proof vinyl cover wipes clean for easy care
- Available polyester outer cover reduces skin shear
- Two-tone covers serve as a reminder to reposition resident



### VINYL WITH POLYESTER COVER

SKIL751170

16" X 16" X 1 1/2"

SKIL751171

18" X 18" X 1 1/2"

## ROHO® SEAT

The ROHO Group is Shape Fitting Technology. Our technology is unique in its ability to conform to the body's shape. Like fluid, it continuously adjusts with every movement. A person can immerse into the soft cells, and feel almost completely buoyant. Gels can't do it. Foam can't do it. Even alternative air products cannot do what ROHO technology can do. No other product can more evenly distribute the forces of body weight or reduce the effects of friction or shear upon tissue. ROHO is the optimal choice for preventing or treating pressure sores (ischemic ulcers) and postoperative wounds.

## ROHO HIGH PROFILE CUSHION

The ROHO HIGH PROFILE Cushion is made of 4" (10 cm) interconnected air cells that provide an adjustable, low deformation, counter-pressure environment that assists in the healing of ischemic ulcers. Appropriate for the high risk client.

\*Also available: Dual Compartment can be adjusted independently to increase positioning and stability for either side-to-side or front-to-back control.

The ROHO HIGH PROFILE Cushion comes in a wide range of sizes, to fit chair widths 12" to 28" / 30cm to 71cm.

Approximate weight of the HIGH PROFILE is 3.8 lbs. / 1.7 kg.

24 - Month Limited Warranty

U.S. Medicare Code: E2622/2623



## ROHO LOW PROFILE CUSHION

The ROHO LOW PROFILE Cushion is made of 2" (6cm) interconnected air cells that maximize stability for the client while protecting from a potential risk of tissue breakdown. Lightweight and portable, this cushion is ideal for the active client.

\*Also available: LOW PROFILE Dual Compartment can be adjusted independently to increase positioning and stability for either side-to-side or front-to-back control.

The ROHO LOW PROFILE Cushion comes in a wide range of sizes, to fit chair widths 12" to 22" / 30cm to 56cm.

Approximate weight of the LOW PROFILE is 2.5 lbs. / 1.1 kg.

24 - Month Limited Warranty

U.S. Medicare Code: E2622/2623



## BED WEDGES



- Designed for patients requiring head, back, or leg elevation
- Conveniently and economically replaces pillows or rolled towels while positioning patients securely
- Provides maximum support and patient comfort
- Comes with removable, washable, navy poly / cotton cover

17 X 19 X 7"

24 X 24 X 7"

17 X 19 X 10"

24 X 24 X 10"

17 X 19 X 12"

24 X 24 X 12"

## NO-SLIP WEDGE

The No-Slip Wedge completes the ITI Positioner Line by adding a highly reliable turning solution. The No-Slip Wedge is manufactured of Pressure Reducing, 1.85 pound density foam that is soft and comfortable, yet supportive enough to hold the torso in a turned position.

- High density foam provides support for patients up to 300 pounds while providing optimal pressure redistribution.
- Consistently and reliably maintains a 30 degree side lying position.
- Contour cuts promote air circulation between the patient and the No-Slip Wedge.
- Available in both covered reusable and bare foam disposable versions. Also available in an optional "All Sides No-Slip" cover design for the ultimate in slip resistance. Call Intensive Therapeutics for details.
- Manufactured of high-density, anti-microbial, hypo-allergenic and latex-free foam.
- Available in a 20" Standard version and a 26" Extra-Long version.
- Reusable version is covered with a skin friendly, moisture proof friction and shear reducing cover with a panel of "No-Slip" fabric on the bottom.
- No-Slip fabric on the bottom of the reusable version prevents the wedge from migrating away from the patient.

Manufactured of high-density, anti-microbial, hypoallergenic foam.

The dimensions are 20" W x 11" Deep x 7" Tall

INT-NSW-R

20" STANDARD





## REXX SERIES



The M.C. ReXX, features a travel range of 7 7/8" to 26" and provides tremendous value with its mobility at all heights, 450 lb safe working load and range of options.

- 7 7/8" low to 26" high
- Mobile or locked at any height
- Travels virtually straight up and down
- Variety of options
- Orthopedic grid deck
- Honeycomb roller bumpers and corner mattress stops

Collapsible 3/4 length  
side rails



Transport  
Storage  
dollies



1/2 length head end  
side rails



3 position  
assist rail



Fixed assist rail

## ADVANTAGE HIGH/LOW BED

Choosing the right bed for your facility can be difficult. Let the Advantage Bed Configurator help you get the right bed for your facility and the features and accessories you need. All while using a top quality simulator to get a detailed and accurate view of the bed you design!

Start with the Advantage High/Low bed frame and see exactly how it works. The virtual control panel allows you to experiment with the different positions and functions that the Advantage frame offers.

Add your choice of four caster options ranging from basic at the head of the bed to foot activated locking double casters at the head and foot of bed.

Continue by selecting which Series and style of head and foot boards you prefer and then select the finish that suites the facility best.

Then choose the right rails for your patients. Customize your rail selection from a range that includes  $\frac{3}{4}$  length collapsible side rails,  $\frac{1}{2}$  length side rails, fixed assist rails and three position assist rails.



Assist Rail



$\frac{1}{2}$  Length Rail




3 Position Assist rail

All beds come standard with a foot board mounted control panel. Finish your customization with a selection of accessories. Choose a pendant and select from various pendant mounts to get the controls right where you need them. Add IV brackets and an IV pole and finish your bed off with a trapeze.



Do all this with the click of a button and view the changes and variations at every step in the viewing window! Always know exactly what to expect when buying a bed with the Advantage Bed Configurator



## BARIATRIC ELECTRIC BED

The BA8000 Bariatric Bed features a sleeping surface that is comprised of 4 aesthetically pleasing smooth sided tubular 'box' sections that adjust in seconds from 35" to 48" widths. A clear span under bed provides unobstructed access for lifts and over bed tables. The bed travels straight up and down avoiding caster movement and damage to walls during operation. A 9" low position allows for use as standard low bed. Four locking, hospital grade, 5" casters come standard in this easily portable bed.

BA8000

## PIVOTING REHAB BED



Transform a standard bed into the Total Solution! The PR1000 Pivoting Rehab features a low height for patient safety, restraint free pivoting assist bars, maintenance free embedded controls, special pressure redistribution sleep surface, four total swivel lock casters, special ergonomic assisting bars, and unimpeded under bed clearance. The perfect bed for safe egress and ingress, cardiac chair position. Great for in room physical therapy, applying lift slings to patient, and for a good position to assist feeding and 360 degree access to patient.

Warranty: 15 year limited on frame, 3 year on drive system, 1 year on wood products and casters, lifetime on welds.

PR1000

## MEDLITE BEDS

MedLite Beds feature premium Linak™ motors that provide smooth reliable operation. Recessed headspring supports and foot assembly reduce the chance of pinching and injuries. Lightweight beds make transporting easier. Strong beds have an extra steel plate welded to the middle of the frame, painted frames with enamel finish and clear coat for a chip and scratch-resistant finish, plastic wood grain end panels and heavy-duty remotes that resist breakage and moisture.

Easy to use and clean

Full-Electric, Semi-Electric and Full-Electric Low Beds are available.

### Specifications Values

Bariatric	No
Bed Style	Full Electric
Deck Style	Split Spring
Includes Staff Control	No
Latex Free	Yes
Length Inches	80 in
Max Height Inches	22.5 in
Max Length Inches	80 in
Max Width Inches	36 in
Min Height Inches	14.5 in
Min Length Inches	80 in
Min Width Inches	36 in
Product Use	Homecare
Weight Capacity Pounds	450 lb



MEDLMDR107003L

## INVACARE BED PACKAGE



INVA-CBED23-1000

Invacare's Home Care Bed offers a number of innovations that translate into long-term savings over the product lifecycle.

5411 IVC Full-Electric Bed

6630 Reduced Gap Half Length Rails

Meets Fire Regulation Code 16 CFR 1633 and 16 CFR 1632

SPS1080 Solace Prevention Foam Mattress

Back Height: 15" min. - 23" max. (Measurement from floor to bed deck.)

Dimensions: Dimensions Overall: 36"W x 88"L

Sleep Surface: 36"W x 80"L

Product Weight Capacity: 350 lbs. patient weight


450 lbs. inclusive of all accessories

Warranty: All Welds: 5 years Mechanical/electrical: 2 years



## MODEL A FOAM REPLACEMENT

The Model A replacement mattress is a 6" multi layered foam mattress. Features a 1.5" Memory foam top layer, a 4.5" fire retardant foam core, 2-3" foam side rails to prevent slipping. The zippered cover has a stretch fabric cover, anti-slip bottom, and colour choice of blue or burgundy.




## ROHO SELECT-AIR MATTRESS SYSTEM

The Roho SelectAir low air loss mattress systems are a complete line that offers the most cost-effective alternative for skin protection and pressure sore prevention. SelectAir is reliable, whisper quiet during operation, portable, lightweight, low on energy consumption, with a user friendly control panel, microprocessor controls and a firm mattress mode. Features a timer to facilitate nursing care during dressing changes and patient repositioning.



## PRISM MATTRESS SYSTEMS

Fatigue tests support VPF's strong cell consistency compared to tests done on conventional foam. Gatching notches allow mattress to move with bed frame for better patient positioning. Various configuration and grades of Reflex Plus™ brand foam available as Prism® inserts subject to your patients needs: All components, including foam insert and cover are replaceable. Product carries a 7 year warranty. Cover is antimicrobial Staph-CHEK® with built in handles for easy transport. VELCRO® fasteners secure Prism® to standard bed frames.





## SUPRA DPS MATTRESS REPLACEMENT SYSTEM



### MATTRESS LOW AIR / ALT PRESS, DIGITAL

Medline's state-of-the-art Supra DPS mattress is our most feature-driven mattress replacement system.

Functions include alternating pressure, low air-loss, static float, auto firm and seat inflation

Cell-on-cell design prevents bottoming out during a power outage

Digital pump with choice of four alternating pressure cycle times

Upgraded low-friction, anti-shear stretch cover

Low- pressure alarm with audio and visual indicators

400-lb. weight capacity

Also available with 4" raised air bolsters to help reduce the risk of patient falls and entrapment

MEDLMDT24SUPRADPS

## ALTERNATING AIR PRESSURE PUMP



AMG707200

The Alternating Air Pressure Pump is a popular solution for minimizing and preventing pressure sores and decubitus ulcers. This compact unit, CSA approved, features a hospital grade plug, an easy-to-use pressure control dial, an on-off switch, and a rippled air mattress made up of over 150 individual air bubbles that are inflated and deflated alternately by the control unit. The rippled air mattress is constructed from bacteriostatic medical grade P.V.C. and inflates to approximately 31" x 74" x 1 1/2" (78,4 cm x 188 cm x 3,8 cm).

- Power: 120V / 60 Hz
- 10-foot long SJT Power Cord
- Includes a 5 minute timer
- Pressure Range: 50-105 mmHg
- 2 standard air outlets (4/16")

## HYDRAULIC STANDARD PATIENT LIFT WITH SIX POINT CRADLE



Silver vein steel construction provides maximum strength.  
6 point cradle design.  
Safely raises or lowers individuals from any stationary position.  
High performance hydraulics raise or lower individuals gradually and safely.  
Comes standard with 5" casters.  
Easy to operate caster brakes provide additional security.  
Adjustable width base.  
Limited lifetime warranty (frame).  
1 year limited warranty (pump)

DRIVE13023-SV

### STANDARD SLING



Sling Points: 2  
13019 - Chain  
Design: Mesh  
Material: Polyester  
Dimensions: 39" x 25"  
Cradle Points: 2 or 6  
Weight Capacity: 330 lbs

### SLING WITH COMMODE OPENING



Sling Points: 2  
13019 - Chain  
Design: Mesh  
Material: Polyester  
Dimensions: 39" x 25"  
Cradle Points: 2 or 6  
Weight Capacity: 330 lbs

**\*\*Call for more sling options**  
*Chains sold separately*

## PAINTED HYDRAULIC LIFT



INVA9805P

Six-point swivel bar attachment easily adapts to all sling styles and positions  
20" - 64" range allows for patient to be picked up from lying position on floor  
Pump handle can rotate from side to side for convenience of attendant  
Offset mast and boom style provide better lift path, maximizing full travel range  
Elongated handle design serves a broad range of operator heights and offers improved handling ability

### STANDARD SLING



### SLING WITH COMMODE OPENING



### SLING CHAINS & STRAPS



**\*\*Call for more sling options**  
*Chains sold separately*

Overall Height: 20" - 64"  
Clearance: 6.5"  
Base Length: 46.75"  
Back Width: Open: 42.25", Closed: 22"  
Caster Options: Front: 5", Rear: 5"  
Product Weight Capacity: 450 lbs.  
Product Weight: 68 lbs.  
Warranty: 3 years on lift, 1 year on hydraulic pump

## HEALTH O METER® PROFESSIONAL UPRIGHT SCALE



The Health-O-Meter® Model 402KL\* is a 350 pound/160Kg capacity balance beam column scale with height measuring rod from Health-O-Meter®, one of the world's leading manufacturers of personal and professional scales. The scale includes counterweight capacity extenders to increase it's capacity to 500 pound/227 Kg.

Tough and reliable, this balance beam scale is built to last for years in even the most demanding weighing applications. No need for batteries or A/C adapters. The scale's easy to read dual poise balance beam weighing mechanism makes weighing fast and easy and with the height measuring rod you can measure height and weight at the same time. If you do not need to measure height consider the Model 400KL, an identical performing scale without the height measuring rod.

HEALTH402KL	FIXED BEAM
HEALTH450KL	ROTATING BEAM

## 3 PANEL MOBILE PRIVACY SCREEN

- Caster wheels allow screen to be moved easily.
- Folding panels facilitates storage.
- Opaque white curtain panels enhance privacy.
- Lightweight anodized aluminum construction.

Dimensions:

53.3" (135 cm) height

75.6" (192 cm) width

Ground Clearance: 19" (48 cm)

Meets NFPA 701 Standards Flammability Testing

AMG770882



## AMG 3-POSITION RECLINER

- May be adjusted to 3 different ergonomic positions.
- High weight capacity of 300 lb (136 kg).
- Includes full armrest, elevating footrest and activity tray.
- Quality high density foam in the chair seat and back provide patients with optimal support.
- Full length armrests provide patients with additional comfort.
- Easy to transport and maneuver with swivel casters, and to set in place due to its rear locking casters.
- Made of easy to clean upholstery to facilitate maintenance.



Total Height: 116 cm  
Seat Width: 51 cm  
Seat Depth: 51 cm  
Armrest Width: 9 cm  
Net Product Weight: 38 kg  
Weight Capacity: 300 lb

May be set to 3 positions:  
Upright Angle: 117 degrees  
Second Angle: 129 degrees  
Third Angle: 142 degrees

Available in:  
Blue, Rosewood and Jade

AMG770050	BLUE
AMG770051	ROSEWOOD
AMG770052	JADE



## ECONOMY 2-HOOK I.V. STAND

This I.V. stand has a low profile, four point base, twin wheeled swivel casters, and a circular "Ram's Horn" design. As well, the Economy 2-Hook I.V. Stand features a positive friction grip, chrome-plated steel pole, stable base design, and friction grip lock adjustment.

AMG775751



## AUTOMATIC OVERBED TABLE

This adjustable overbed table is easily assembled. Features include four 2" swivel casters, and a chrome-plated "H" base.



Item Number	Weight	Table Top (L x W)	Adjustable Height Range	Color
<b>770-870</b>	11.5 kg	31.5" x 15" 80 x 38.1 cm	30" - 44" 76.2 x 111.8 cm	Walnut
<b>770-872</b>	11.1 kg	31.5" x 15" 80 x 38.1 cm	29" - 43" 73.7 x 109.2 cm	Walnut
<b>770-873</b>	11.1 kg	31.5" x 15" 80 x 38.1 cm	29" - 43" 73.7 x 109.2 cm	Light brown

AMG770870

WALNUT

TILT-TOP

AMG770872

WALNUT


REG TOP

AMG770873

LIGHT BROWN

REG TOP

## ASSISTA-RAIL




The AssistaRail™ is an "inverted U" style bed rail for home beds, providing a unique combination of occupant safety and in-bed mobility.

### Simple Installation.

Slide the Assista-Rail™ between the box spring and mattress, secure to bed frame with mandatory safety strap. Traction grips help to prevent rail from shifting during use.

### Assista-Cane™ Add On Kit (sold separately).

A sturdy handle that mounts to the Assista-Rail™ and can be fixed at any angle to provide a safe weight bearing surface. Minimizes the risk of falls by reducing body twisting.



### Assista-Cane™ Leg Extender (sold separately).

Increases the Assista-Cane™ height from 32.5" / 83 cm to 35.5" / 90 cm for taller users.

### Dual Rail System (sold separately).

Conversion kit allows an Assista-Rail™ to be placed on both sides of a bed.

HCAST-S-1

## SUPERPOLE



### Maintain Confidence and Independence.

Place support exactly where it's needed for valuable assistance beside a bed, bath, toilet or chair.

### Quick, Easy, Secure Installation.

Simply turn the screw at the bottom of the pole and the pole will expand between the floor and ceiling, creating unparalleled sturdiness and reliability.

### Versatile.

Need to remove or relocate the pole? Not a problem! Rubber pads on the floor and ceiling plates prevent marking, even on stipple ceilings.

### Add On Kits.

The SuperPole™ also serves as a mounting support for a complete range of options that can be varied to meet changing needs: The SuperBar™, SuperTrapeze™, SuperTray™ & Safety Island™.

HCSP-UF

## SILENT KNIGHT PILL CRUSHER



MEDLNONSK0100

- Extremely quiet and works with very little effort on the operator's part
- Special 8mL pouches help eliminate cross-contamination
- Lifetime warranty
- With 50 pouches

### Pill Crusher Bags

DON-PCS8133C 50/BG 20/BX

## DIE CAST CHROME PILL CRUSHER



MEDLNON131000

- Die-cast pill crusher pulverizes even the hardest tablets with very little effort
- One hand operation
- Heavily chromed components and stainless steel fittings to prevent corrosion
- Rarely needs cleaning because pills are crushed between two standard 3/4-oz souffle cups

## EZ9 - AUTOMATIC AUTOCLAVE



With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer EZ9 fully automatic autoclave will fill, sterilize, exhaust and dry at the touch of a button.

### EZ9 Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying.
- Hepa filter 0.2 µm air filter (replaceable).
- Pre-loaded cycles for: unwrapped instruments, wrapped/packs, glassware and additional drying cycle.
- All program parameters can be customized and stored.
- Automatic shut off at the end of both the sterilization and dry cycles.
- Audible and visual cycle interruption alert.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir.
- Constant monitoring of temperature and pressure for sterilization.
- Power outage recovery system.
- Door design prevents steam from coming in contact with control panel.
- Optional internal printer documents: date, time, temperature and pressure.
- 510k clearance for dental handpieces.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- International certifications.
- 2 year parts and labor warranty.

TUTTEZ9P WITH PRINTER

## KANGAROO™ EPUMP ENTERAL FEEDING PUMP



The Kangaroo™ ePump system is the latest innovation in enteral feeding pump technology. When patients require nutrition through continuous feeding, intermittent feeding or feeding and flushing, Kangaroo™ ePump enteral feeding pump delivers in one compact, easy-to-use device.

KEND4824005


## KANGAROO™ 924 ENTERAL FEEDING PUMP

The Kangaroo™ 924 Enteral Feeding Pump features pre-set volume and overflow infusion safeguards to ensure the correct and consistent delivery of nutrition. The easy-to-use controls make set-up and operation of the unit intuitive. *The Kangaroo™ 924 pump is used exclusively with the Kangaroo™ 924 pump sets.*



KEND392455

## ENTERALITE INFINITY ADMINISTRATION SETS



MAQ-INF0500 500ML

The Enteralite Infinity 1200 ml delivery set (INF1200) comes with pre-attached pump set, patented in-line occluder for automatic free-flow protection, built-in hanger, and a large top-fill opening with leak-proof cap to minimize formula spills and waste.



MAQ-INF1200 1200ML

All Moog enteral delivery sets are DEHP and latex-free.

Non-IV Compatible

All mating adapters meet AAMI guidelines to prevent potentially hazardous situations in connecting to the patient's feeding tube.



## SYRINGE TIP CAPS



Sterile, single-use, multiple packaged syringe tip cap; polyolefin plastic formulation; fits luer lock or regular luer syringe tips, covers entire luer skirt of syringe.

SHER8881682085

## ECCENTRIC LEUR TIP

An eccentric tip has the nozzle at the end of the syringe not in the centre of the syringe but at the side. This causes the blade attached to the syringe to lie almost in line with the walls of the syringe itself and they are used when the blade needs to get very close to parallel with the skin (when injecting into a surface vein or artery for example).



20CC  
60CC

25/BX  
25/BX

TERUSS20ES  
TERUSS60ES

## CURVED TIP SYRINGE



Monoject curved tip syringe are ungraduated but ideal for pressure irrigation.

12CC

50/BX

MONO412

## TUBERCULIN SYRINGE

TB syringes have clear, bold markings to ensure accurate dosing and control.

1CC

100/BX

TERUSS01T



## CATHETER TIP SYRINGES



35CC	30/BX	SHER8881535770
60CC	25/BX	TERUSS60C

The catheter tip, which is essentially a slip tip but longer and tapered, making it good for pushing into things where there the plastic taper can form a tight seal.

## LEUR-LOK SYRINGES



1CC	100/BX	BD309628
3CC	200/BX	BD309657
3CC	100/BX	TERUSS03L
5CC	100/BX	TERUSS05L
10CC	100/BX	BD309604
10CC	100/BX	TERUSS10L
12CC	80/BX	SHER8881512878
20CC	40/BX	BD309661
20CC	50/BX	TERUSS20L2
30CC	25/BX	TERUSS30L
30CC	40/BX	BD309650
60CC	25/BX	TERUSS60L
60CC	40/BX	BD309653
60CC	100/BX	SHER8881560125

Luer-Lok fittings are securely joined by means of a tabbed hub on the female fitting which screws into threads in a sleeve on the male fitting.

## SLIP TIP SYRINGES



3CC	100/BX	TERUSS03S
3CC	100/BX	SHER118030055
5CC	100/BX	TERUSS05S
5CC	100/BX	BD301603
10CC	100/BX	TERUSS10S
20CC	50/BX	TERUSS20S2
20CC	40/BX	BD301625
30CC	25/BX	TERUSS30S

Luer-Slip fittings simply conform to Luer taper dimensions and are pressed together and held by friction (they have no threads).

## MONOJECT™ INSULIN SAFETY SYRINGE



The Monoject™ safety syringe is a proven engineering control and preferred by millions of users each day. Our exclusive design helps keep clinicians safe with maximum clinical flexibility. The Monoject™ safety syringe complies with OSHA Bloodborne Pathogens Standard for Engineering Controls.

100/BX

SHER8881511144	3/10 mL 29 G x 1/2" (0.33 mm x 1.3 cm)
SHER8881511136	1/2 mL 29 G x 1/2" (0.33 mm x 1.3 cm)
SHER8881511110	1 mL 29 G x 1/2" (0.33 mm x 1.3 cm)
SHER8881511344	3/10 mL 30 G x 5/16" (0.305 mm x 0.8 cm)
SHER8881511336	1/2 mL 30 G x 5/16" (0.305 mm x 0.8 cm)
SHER8881511310	1 mL 30 G x 5/16" (0.305 mm x 0.8 cm)

## MAGELLAN™ SAFETY NEEDLE AND SYRINGE COMBINATIONS



Magellan™ safety needle and syringe combinations are needle based safety devices preattached to syringes that protect the healthcare worker from accidental needlestick injuries. The preattached combinations save valuable nursing time and improve clinical efficiencies. The product is easy to use and the robust safety shield ensures clinicians, patients, and downstream workers receive unparalleled protection from a sharps injury.

100/BX

SHER8881811558	1 mL 25 G x 5/8" (0.508 mm x 1.6 cm)
SHER8881811510	1 mL 25 G x 1" (0.508 mm x 2.5 cm)
SHER8881811310	1 mL 23 G x 1" (0.635 mm x 2.5 cm)
SHER8881833210	3 mL 22 G x 1" (0.711 mm x 2.5 cm)
SHER8881833215	3 mL 22 G x 1-1/2" (0.711 mm x 3.8 cm)
SHER8881833310	3 mL 23 G x 1" (0.635 mm x 2.5 cm)
SHER8881833510	3 mL 25 G x 1" (0.508 mm x 2.5 cm)
SHER8881833558	3 mL 25 G x 5/8" (0.508 mm x 1.6 cm)

## TERUMO® INSULIN SYRINGES

Low Dead Space design, for use with U-100 insulin only, with attached needle. Based on extreme needle sharpness and smooth plunger movement, TERUMO® insulin syringes are designed for maximum comfort. Clear, bold scale markings assure dosage control and accuracy.

TERUSS30M3009	3/10 mL	30G x 3/8"	100/BX
TERUSS30M2913	3/10 mL	29G x 1/2"	
TERUSS10M2913	1 mL	29G x 1/2"	
TERUSS05M2813	1/2 mL	28G x 1/2"	
TERUSS10M2813	1 mL	28G x 1/2"	
TERUSS10M2713	1 mL	27G x 1/2"	



## BD SAFETYGLIDE™ SYRINGE

3 mL BD SafetyGlide™ syringe with shielding subcutaneous injection needle, regular bevel, regular wall. Detachable needle.

BD309570 3 mL 25G x 5/8"



3 mL BD SafetyGlide™ syringe with shielding intramuscular injection needle, regular bevel, regular wall. Detachable needle.

50/BX

BD305906 3 mL 22G x 1 1/2"

BD305909 3 mL 21G x 1 1/2"



## BD LUER-LOK SYRINGE WITH NEEDLE

3 BD PrecisionGlide™ needle combination. Regular bevel, regular wall. Detachable needle.

BD309570 3 mL 25G x 5/8"



100/BX

8BX/CS

## TERUMO® DENTAL NEEDLES


TERUMO® Dental Needles are latex free and are available in both short and long sizes for conduction anesthesia and infiltration anesthesia.

TERUDN2541 25G 1 5/8"

TERUDN2722 27G 7/8"

TERUDN2721 27G 13/16" SHORT

TERUDN3021 30G 13/16" SHORT



100/BX

## TERUMO® HYPODERMIC NEEDLES



Terumo's Hypodermic needles use a unique double bevel design and precision manufacturing to provide ultra sharp needles, which translate into greater comfort for the patient.

100/BX

10BX/CS

TERUNN1825R	18G x 1", Thin Wall Needle
TERUNN1838R	18G x 1 1/2", Thin Wall Needle
TERUNN2025R	20G x 1", Ultra-Thin Wall Needle
TERUNN2038R	20G x 1 1/2", Ultra-Thin Wall Needle
TERUNN2125R	21G x 1", Ultra-Thin Wall Needle
TERUNN2138R	21G x 1 1/2", Ultra-Thin Wall Needle
TERUNN2225R	22G x 1", Ultra-Thin Wall Needle
TERUNN2516R	25G x 5/8", Thin Wall Needle
TERUNN2525R	25G x 1", Thin Wall Needle
TERUNN2538R	25G x 1 1/2", Thin Wall Needle
TERUNN2713R	27G x 1/2", Regular Wall Needle
TERUNN2238R	22G x 1 1/2", Ultra-Thin Wall Needle
TERUNN2325R	23G x 1", Ultra-Thin Wall Needle



## MAGELLAN™ SAFETY NEEDLES

Magellan™ safety needles are needle based safety devices that protect the healthcare worker from accidental needlestick injuries. The easy one handed design and robust safety shield ensure the clinician is protected against any exposures from sharps injury.



50/BX

SHER8881850810	18 G x 1" (1.27 mm x 2.5 cm)
SHER8881850815	18 G x 1-1/2" (1.27 mm x 3.8 cm)
SHER8881850910	19 G x 1" (1.067 mm x 2.5 cm)
SHER8881850915	19 G x 1-1/2" (1.067 mm x 3.8 cm)
SHER8881850010	20 G x 1" (0.902 mm x 2.5 cm)
SHER8881850015	20 G x 1-1/2" (0.902 mm x 3.8 cm)
SHER8881850158	21 G x 5/8" (0.813 mm x 1.6 cm)
SHER8881850110	21 G x 1" (0.813 mm x 2.5 cm)
SHER8881850115	21 G x 1-1/2" (0.813 mm x 3.8 cm)
SHER8881850210	22 G x 1" (0.711 mm x 2.5 cm)
SHER8881850215	22 G x 1-1/2" (0.711 mm x 3.8 cm)
SHER8881850358	23 G x 5/8" (0.635 mm x 1.6 cm)
SHER8881850310	23 G x 1" (0.635 mm x 2.5 cm)
SHER8881850558	25 G x 5/8" (0.508 mm x 1.6 cm)
SHER8881850510	25 G x 1" (0.508 mm x 2.5 cm)

## SURFLO® WINGED INFUSION SETS

Take a look at our wings versus the competition and you'll see why SURFLO® is a very well recognized brand. Terumo's soft, flexible interlocking wings readily conform to the body's contour. The wing is angled for ease of positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen. The precisely honed needle with Ultra-Thin wall design assures ease of penetration and continued access with a single insertion. DEHP Free.



50/BX

TERUSV*19BLS	19G x 3/4" Thin Wall Needle, 3 1/2" Tubing
TERUSV*21BLS	21G x 3/4" Ultra Thin Wall Needle, 3 1/2" Tubing
TERUSV*23BLS	23G x 3/4" Ultra Thin Wall Needle, 3 1/2" Tubing
TERUSV*25BLS	25G x 3/4" Thin Wall Needle, 3 1/2" Tubing
TERUSV*25EL	25G x 1/2" Regular Wall Needle, 8" Tubing
TERUSV*27EL	27G x 1/2" Regular Wall Needle, 8" Tubing
TERUSV-18BLK	18G x 3/4" Thin Wall Needle, 12" Tubing
TERUSV-19BLK	19G x 3/4" Thin Wall Needle, 12" Tubing
TERUSV-21BLK	21G x 3/4" Ultra Thin Wall Needle, 12" Tubing
TERUSV-22BLK	22G x 3/4" Ultra Thin Wall Needle, 12" Tubing
TERUSV-23BLK	23G x 3/4" Ultra Thin Wall Needle, 12" Tubing
TERUSV-25BLK	25G x 3/4" Thin Wall Needle, 12" Tubing

## DISPOSABLE LATEX-FREE TOURNIQUETS, DISPENSER BOX (ROLL)

A Skin-Friendly Alternative for Sensitive Skin

Non-Irritating

- Ideal for use on patients with sensitive skin
- Eliminates risk of allergic reaction to latex proteins.
- Its smooth texture makes it feel more comfortable on skin.

Disposable for Improved Hygiene

- Single-use to reduce the risk of cross-contamination between patients
- Protects workers and patients from potentially dangerous allergens

Excellent Elasticity

- Its high degree of stretch makes this tourniquet a suitable alternative to using a latex product

User-Friendly Dispenser Box Available


- Dispenser box with perforated design allows for quick dispensing
- Tourniquets are pre-notched every 18" to facilitate separation



AMG118782 25/RL

## VELCRO TOURNIQUET

This superior quality tourniquet is made of durable blue rubber with a Velcro closure. Latex free.  
Dimensions: 9" x 1" (22.9 cm x 2.5 cm)



AMG118765	CHILD	1 EA
AMG118770	ADULT	1 EA

## IV SALINE SOLUTION




0.9% Sodium Chloride Injection, USP

BAXTJB1322P	250 ML	30/CS
BAXTJB1323	500 ML	24/CS
BAXTJB1324	1000 ML	12/CS

## SECONDARY MEDICATION SETS (10 DROPS / ML)

Secondary Medication Set with  
male Luer lock adapter, hanger  
37" (94 cm)



BAXTJC7453      48/CS

## CONTINU-FLO SOLUTION SETS (10 DROPS / ML)

CONTINU-FLO Solution Set with 2 Luer activated valves, male Luer  
lock adapter  
104" (2.6 m)



BAXTJC8519      48/CS



6.6 L BD™ chemotherapy sharps collector. Red.  
(12/ca)

BD300459



1.4 L BD™ tray sharps collector. Yellow with funnel entry.  
(36/ca)

BD300460



5.1 L BD™ next generation sharps collector. Yellow.  
(20/ca)

BD300974



3 Gallon Bemis Sharps Collector. Rotating cylinder top.  
Uses bracket #423 or #425, or cabinet 1031 for wall mounting.  
(12/cs)

BEMIS333-030



SharpSafety™ Sharps Container, Transportable, Red, 1 Quart  
20/CS

DEVO8303SA



SharpSafety™ Safety In Room Sharps Container Counterbalance Lid, Transparent Red 5 Quart

DEVO8507SA



SharpSafety™ Sharps  
Container Phlebotomy,  
Yellow 1 Quart  
100/CS

DEVO8906



SharpSafety™ Sharps  
Container, Phlebotomy,  
Red, 1 Quart  
10/CS

DEOV8900SA



Sharps-A-Gator™ Tortuous  
Path Sharps Container,  
Transparent Yellow, 3.5  
Gallon  
(10/CS)

GC31300646



Sharps-A-Gator™ Sharps  
Container, Tortuous  
Path, Transparent Red,  
5 Quart  
(30/CS)

GC31144010



SharpSafety™ Sharps  
Container, Horizontal  
Drop, Yellow, 2 Gallon  
(20/CS)

DEVO8967Y



One-hand, one-step,  
no-touch sharps depositing.  
Translucent Yellow  
(30/CS)

GC31300901


2.5 qt

GC31300919

5 qt



## VINYL GLOVES




### DERMASHEER

- Vinyl Medical Examination Gloves
- Powder Free
- Latex Free

BOW1000S	SMALL
BOW1000M	MEDIUM
BOW1000L	LARGE
BOW1000XL	EXTRA LARGE

## LATEX GLOVES




### DERMASENSE

- Latex Medical Examination Gloves
- Textured
- Powder Free

BOW10110	EXTRA SMALL
BOW10111	SMALL
BOW10112	MEDIUM
BOW10113	LARGE

## NITRILE GLOVES



### DERMATUFF

- Nitrile Medical Examination Gloves
- Textured
- Powder Free
- Latex Free

BOW30310	EXTRA SMALL
BOW30311	SMALL
BOW30312	MEDIUM
BOW30313	LARGE
BOW30314	EXTRA LARGE

## MICRO-TOUCH® ELITE®

- Polyvinyl Chloride
- DEHP-Free
- Cream color
- Smooth finish
- Beaded cuff
- Stretch for better fit and feel
- For low-risk applications
- Length: 230 mm (9.1")
- Powder free
- Non-sterile



ANS3091	SMALL
ANS3092	MEDIUM
ANS3093	LARGE
ANS3094	XLARGE

100/BX      10BX/CS

## MICRO-TOUCH® DERMACLEAN®

- Low latex protein\*
- Cream color
- Textured finish
- Beaded cuff
- Length: 241 mm (9.5")
- Powder free
- Non-sterile



ANS3180	XSMALL
ANS3181	SMALL
ANS3182	MEDIUM
ANS3183	LARGE
ANS3184	XLARGE

100/BX      10BX/CS

\*\*XL = 80/BX

## MICRO-TOUCH® AFFINITY®

- Neoprene
- Green color
- Textured fingertips
- Beaded cuff
- Length: 239 mm (9.4")
- Powder free
- Non-sterile
- Easy donning due to polymer coating
- Tested for use with chemotherapy drugs



ANS3770	XSMALL
ANS3771	SMALL
ANS3772	MEDIUM
ANS3773	LARGE
ANS3774	XLARGE

100/BX      10BX/CS

## MICRO-TOUCH® NITRILE TEXTURED

- Nitrile
  - Blue color
  - Textured finish
  - Beaded cuff
  - Length: 240 mm (9.4")
  - Excellent chemical and puncture resistance
  - Tested for use with chemotherapy drugs.
- ASTM D6978



ANS4000	XSMALL
ANS4002	SMALL
ANS4004	MEDIUM
ANS4006	LARGE
ANS4008	XLARGE

100/BX      10BX/CS

## MICRO-TOUCH® NITRILE

- Nitrile
  - Blue color
  - Thinner to enhance sensitivity
  - Textured fingertips
  - Beaded cuff
  - Excellent chemical and puncture resistance
  - Tested for use with chemotherapy drugs.
- ASTM D6978
- Length: 240 mm (9.4")



ANS6034300	XSMALL
ANS6034301	SMALL
ANS6034302	MEDIUM
ANS6034303	LARGE
ANS6034304	XLARGE

200/BX      10BX/CS

## MICRO-TOUCH® PLUS

- Low latex protein\*
- Cream color
- Textured finish
- Beaded cuff
- Length: 241 mm (9.5")
- Easy donning due to polymer coating



ANS6015300	XSMALL
ANS6015301	SMALL
ANS6015302	MEDIUM
ANS6015303	LARGE
ANS6015304	XLARGE

150/BX      10BX/CS

**MICRO-TOUCH® NITRAFREE™**

- Nitrile
- Pink color
- Chemical accelerator-free\*
- Textured fingertips
- Beaded cuff
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs. ASTM D6978
- Length: 240 mm (9.4")



ANS6034510	XSMALL
ANS6034511	SMALL
ANS6034512	MEDIUM
ANS6034513	LARGE
ANS6034514	XLARGE

100/BX

10BX/CS

**SYNSATION®**

- Polyvinyl chloride
- DEHP-Free
- Clear
- Smooth finish
- Beaded cuff
- For low-risk applications
- Length: 230 mm (9.1")



ANS6043300	XSMALL
ANS6043301	SMALL
ANS6043302	MEDIUM
ANS6043303	LARGE
ANS6043304	XLARGE

200/BX

10BX/CS

**WALL BRACKET**

- Metal
- Horizontal or vertical installation
- Fits all sizes of Ansell examination glove dispensers

ANS5701U 24/CS

## DISPOS-A-GLOVE® STERILE



- Copolymer
- Clear
- Smooth finish
- Length: 269 mm (10.6")



JOJOMDG601	SMALL
JOJOMDG701	MEDIUM
JOJOMDG801	LARGE
100/BX	14BX/CS

## VINYL GLOVES



- Vinyl
- Non-sterile
- Pre-Powdered

AMD9993A	SMALL
AMD9993B	MEDIUM
AMD9993C	LARGE
AMD9993D	XLARGE
100/BX	10BX/CS

## ACCUTOUCH NITRILE GLOVES



MEDLMDS192084	SMALL
MEDLMDS192085	MEDIUM
MEDLMDS192086	LARGE
MEDLMDS192087	XLARGE

100/BX      10BX/CS

Our Accutouch Chemo Nitrile Exam Gloves are both latex and powder free.

Nitrile gloves conform comfortably to the shape of the hands and return to original shape after removal.

They feature excellent tactile sensitivity and high tensile strength for durable barrier protection. Textured for a secure grip, with tight cuffs to minimize accidental exposure to blood-borne pathogens and other contaminated fluids.

Approved for use with chemotherapy drugs, formaldehyde and glutaraldehyde.



**CONFORM® PROCEDURE**

- Cream color
- Bisque finish
- Unbeaded cuff
- Length: 312 mm (12.3")

ANS786960	6.0
ANS787065	6.5
ANS787170	7.0
ANS787275	7.5
ANS787380	8.0
ANS787485	8.5

*1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case*

**DERMA PRENE® ISOTOUCH®**

- Polyisoprene
- Green color
- Smooth finish
- Beaded cuff
- Ergonomic for better comfort
- Tested for use with chemotherapy drugs. ASTM D6978
- Length: 306 mm (12")

ANS20685255	5.5
ANS20685260	6.0
ANS20685265	6.5
ANS20685270	7.0
ANS20685275	7.5
ANS20685280	8.0
ANS20685285	8.5
ANS20685290	9.0



*1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case*

**PERRY® STYLE 42®**

- Cream color
- Smooth finish
- Beaded cuff
- Ergonomic for better comfort
- Length: 300 mm (11.8")
- Contains natural rubber latex

PER5711100	5.5
PER5711101	6.0
PER5711102	6.5
PER5711103	7.0
PER5711104	7.5
PER5711105	8.0
PER5711106	8.5
PER5711107	9.0

*1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case*



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

**Richmond: 1-800-663-0047**  
**Edmonton: 1-800-561-3687**

## GLOVE BOX DISPENSER - SINGLE



Glove Box Dispenser - Single with Built-in Flexible Spring  
Holds one box of gloves

Two-way keyholes for vertical or horizontal wall mounting  
White powder coated steel

BOWMAN-GB001

## GLOVE BOX DISPENSER - DOUBLE



Glove Box Dispenser - Double  
Holds two boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting  
White powder coated steel

BOWMAN-GB002

## GLOVE BOX DISPENSER - TRIPLE



Glove Box Dispenser - Triple  
Holds three boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting  
White powder coated steel

BOWMAN-GB003

## GLOVE BOX DISPENSER - SINGLE WITH FLEXIBLE SPRING



Glove Box Dispenser - Single with Flexible Spring  
Holds one box of gloves with built-in flexible spring

Two-way keyholes for vertical or horizontal wall mounting  
Clear PETG plastic

BOWMAN-GP013

## GLOVE BOX DISPENSER - DOUBLE



Glove Box Dispenser - Double  
Holds two boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting  
Clear PETG plastic

BOWMAN-GP014

## GLOVE BOX DISPENSER - TRIPLE



Glove Box Dispenser - Triple  
Holds three boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting  
Clear PETG plastic

BOWMAN-GP015



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687

## GLOVE BOX DISPENSER - QUAD



Glove Box Dispenser - Quad  
Holds four boxes of gloves horizontally

Keyholes for wall mounting  
Clear PETG plastic

BOWMAN-GP340

## GLOVE BOX DISPENSER - SINGLE WITH FLEXIBLE SPRING



Glove Box Dispenser - Single with Flexible Spring  
Holds one box of gloves with built-in flexible spring

Two-way keyholes for vertical or horizontal wall  
mounting  
Stainless steel

BOWMAN-GS004

## GLOVE BOX CADDY - SINGLE



Glove Box Caddy - Single  
Holds one box of gloves

Attaches to the bottom of protection organizer LD-037  
Quartz ABS plastic

BOWMAN-GD011

## SHOE COVERS

Footwear protection with fluid repellant, slip-on covers.

- Sewn seams provide increased tear strength.
- Elastic opening for easy on-and-off convenience.
- Non-conductive, non-skid and plain options.



MEDICOM8006 XL 100/BX 3BX/CS

Anti-skid, water resistant, blue shoe covers.

*Latex free*

BOWW350 UNIVERSAL 100/BX 10 BX/CS



## BOUFFANTS



The classic Surgeons Cap is made of a soft, reinforced material that is absorbent and comfortable to wear. The Nurses Cap is made of a comfortable spunbonded material which is cool and can be worn well for hours. All caps are Latex Free.

DUKAL321W	21" White	100/BX	5BX/CS
DUKAL324W	24" White	100/BX	5BX/CS

## ISOLATION GOWNS



BOWW301	10/PK	10PK/CS	ELBOW/WRIST/YELLOW
DUKAL301	10/PK	10PK/CS	ELBOW/WRIST/YELLOW
BOWW301WH	10/PK	10PK/CS	ELBOW/WRIST/WHITE
BOWW302	10/PK	10PK/CS	IMPERVIOUS/YELLOW



Our 3-Ply Tissue Gowns are constructed of three layers of tissue bonded into one layer for maximum comfort and absorbency.

GRAHAM70221N 50/CS



## EARLOOP MASKS



BOWB2115 50/BX



AMD2115 50/BX

## PROCEDURE EARLOOP W/VISOR



PRIMAGARD®160 Procedure Ear Loop with attached Visor, Anti-fog Foam, Anti-glare, Indigo

PM-PG4-1573  
25/BX 100/CS

## SURGICAL TIE W/VISOR



PRIMAGARD®120 Surgical Tie with attached Visor, Anti-fog Foam, Anti-glare, Denim Blue

PM-PG4-2592  
25/BX 100/CS

## SURGICAL TIE MASK



PRIMAGARD® Surgical Mask with Soft Horizontal Tie, Blue Swoosh

PM-PG4-2560  
50/BX 300/CS

## N95 MASKS



ALPH695  
20/BX 20/CS

## MANORAPID



Manorapid is a highly effective poly alcohol formula that when rubbed on your hands for 20-30 seconds kills bacteria and fungi.

*That is a lot of germs! 99.999.33% is not your average sanitizer.*

Manorapid has a special lanolin-derivative emollient that ensures your skins stays soft and smooth after each application even under the harshest winter conditions found in Canada.

Manorapid is a proven formula used in scores of countries and for over 40 years has helped prevent infections in hospitals and nursing homes.

REMAS-HAND-120  
REMMANO1000P  
REMMANO1000B

120ML WITH FLIP CAP  
1 L WITH PUMP  
BRACKET FOR 1 L

## DERMOTAN



Dermotan is an imported soap from Germany manufactured with caring substances found in high quality cosmetic products making it very gentle on the skin.

The smallest amount of Dermotan will clean the dirtiest of hands, making this a cost effective product.

Dermotan can also be used in bath water cleansing the most fragile skin; from new born babies to the elderly even with difficult skin problems.

Dermotan is used as a pre scrub prior to using Manorapid in the OR for a perfect combination to keep your hands from drying out even in the winter months when the air is dry and cold.

REMDERM1000

1 L



**PURELL® ADVANCED Hand Rub**

The only hand sanitizer permitted for Professional Healthcare, Food Premise and Commercial use based on current Health Canada efficacy standards\*

\* PURELL Advanced Hand Rub gel formulation. Based on current, higher efficacy standards of the Health Canada GUIDANCE DOCUMENT Human-Use Antiseptic Drugs, EN-162959.

## **PURELL® HAND SANITIZER**

### **The Brand People Know And Trust.**

PURELL Advanced Hand Rub is the first and only sanitizer allowed for Commercial Use, Hospital and Healthcare Professional Use and in Food Premises in Canada.\* That's because it is the first sanitizer to be granted an Exemption by Health Canada for Hospital and Healthcare Professional use, Commercial use and in Food Premises based on the 2009 Health Canada GUIDANCE DOCUMENT Human-Use Antiseptic Drugs.\*

It provides proven persistent effect and kills 99.99% of germs\* including:

- Norovirus (MNV)
- Staphylococcus aureus
- Escherichia coli
- Staphylococcus epidermis
- Serratia marcescens
- Candida albicans
- Aspergillus niger
- Rotavirus
- Rhinovirus
- Adenovirus

This green certified hand rub gel formulation contains ingredients made from natural renewable resources. Fragrance free. Meets EcoLogo hand sanitizer standard for environmental leadership and proven performance. Scientifically advanced, patent-pending 70% formulation outperforms other leading hand sanitizers milliliter for milliliter(1). Includes a nourishing blend of four different skin conditioners. Clinically proven to help maintain skin health(2). North America's #1 instant hand sanitizer.

## PURELL® 70 INSTANT HAND SANITIZER



The power of PURELL® in a 70% ethyl alcohol, fragrance free formula. Effective in destroying certain bacteria. Acceptable for frequent use. Contains vitamin E and glycerin.

GOJO9770-12 354ml

## PURELL® 70 INSTANT HAND SANITIZER MOISTURIZING FOAM

PURELL 70 Moisturizing Foam is specially formulated to be mild to the skin, non-drying and increase skin hydration.



GOJO5695-24  
45ML

GOJO5795-04 535ML



PURELL® Personal Gear  
Retractable Clip

GOJO9608-24

## PURELL® VERSAHOLD™ POINT OF CARE BRACKET



Drive hand hygiene compliance by making PURELL Instant Hand Sanitizer available at the point of patient care. The PURELL VERSAHOLD bracket provides flexible options for attaching to a wide variety of fixtures.

PURELL VERSAHOLD Point of Care Bracket

- Addresses a key barrier to hand hygiene compliance: access to hand hygiene products
- Helps meet WHO and PIDAC compliance standards throughout the patient zone – in patient rooms, on mobile transport, in work areas, wherever needed
- Quick and easy on and off mounting capability, no tools required

## PURELL® SANITIZING HAND WIPES

Enhance guest and patient experience and satisfaction with economical towelettes that clean and sanitize. Individually wrapped, “on-the-go” amenity item for a wide variety of markets, including foodservice, healthcare, travel and hospitality. Made with 62% ethyl alcohol – the same active ingredient as original formula PURELL Instant Hand Sanitizer.



GOJO9022-10      100/BX

## PURELL® ALCOHOL HAND SANITIZING WIPES



For personal hand hygiene to help prevent the spread of certain bacteria •Use anytime for quick clean-up when soap and water are not available.

GOJO9031-06

175/TUB



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687



## PURELL® 70 INSTANT HAND SANITIZER MOISTURIZING FOAM

PURELL 70 Moisturizing Foam is specially formulated to be mild to the skin, non-drying and increase skin hydration.

TFX™ 1200 mL Refill

- Fresh dispensing valve with each refill
- Refill snaps in place with audible click

GOJO5395-02



## PURELL® ADVANCED MOISTURIZING FOAM HAND RUB



GOJO1913-02

1200 mL Refill

- GOJO SANITARY SEALED™ refill helps lock out germs
- Durable, recyclable PET bottle uses 30% less material than standard rigid HDPE
- Patent pending GOJO CONTROLLED COLLAPSE technology keeps bottle shape longer while emptying
- Crystal clarity for visual appeal and easier monitoring of fill level
- Fresh dispensing valve with each refill
- Removable pump promotes easier recycling

GOJO1903-02

ADVANCED HAND RUB 1200ML

## PURELL® LTX-12™ DISPENSER - WHITE



GOJO1920-04

- Reliable touch-free dispensing
- Smart, trouble-free electronics eliminate battery changes in most installations
- Large sight window and skylight for at-a-glance product monitoring
- 1200mL high capacity output for fewer changeouts
- Easily converts to locked dispenser
- Lifetime guarantee even includes batteries
- Fully ADA compliant
- Customization available

## GOJO® CLEAR & MILD FOAM HANDWASH

Spa-inspired, green certified foam hand soap offers luxurious foam lather in a fragrance- and dye-free formulation. Enriched with moisturizers and skin conditioners. USDA certified biobased formulation. EcoLogo certified. Delights users while supporting a healthy environment.

GOJO1311-03

700ML



## GOJO® LTX-7™ DISPENSER - WHITE



GOJO1380-04

Compact 700mL dispenser is ideal for tight spaces. Lifetime Performance Guarantee; including batteries. GOJO SMART ELECTRONICS technology minimizes energy use per actuation, eliminating need to change batteries over the life of the system in most installations. Large sight window, skylight and clear refill bottles make it easy to check fill status. CONTROLLED COLLAPSE refills hold their shape as they empty for a more attractive appearance. Easily converts to locking dispenser. Removable pump makes recycling easy. SANITARY SEALED™ refills made of durable, recyclable PET material, using 30% less plastic than rigid HDPE bottles. Fully ADA compliant.

## GOJO® GREEN CERTIFIED FOAM HAND CLEANER

A mild, biodegradable foam soap formula certified by Green Seal™ and EcoLogo® to meet their joint standard (GS-41/CCD-104) for hand cleaners and hand soaps. The certifications mean that the product has been thoroughly evaluated by these unbiased, third-party organizations and meet standards for being environmentally responsible without sacrificing product performance

GOJO1380-04



## ISAGEL HAND SANITIZING GEL

Ethyl Alcohol Antiseptic No-Rinse Hand Gel  
70% Ethyl Alcohol

Effective, fast-acting handrub kills greater than 99.9% of most transient microorganisms (including MRSA and VRE) in 15 seconds or less. Enriched with emollients, Isagel is gentle on the skin. Compatible with gloves and CHGbased products.



SWEEN1644	118 ML
SWEEN1645	621 ML
SWEEN7041	800 ML

## DISPENSER BRACKET



Made of very sturdy plastic.  
Made specifically to hold 621 ml  
sized pump bottles. Theft resistant.

SWEEN0089

## DISPENSER FOR BAG



Safe and compact, easy to  
dispense. Designed for 800 ml  
and 1000 ml sized bags.

SWEEN7251

## DRIP TRAY



Accessory to be use with 7251  
to prevent spillage.

SWEEN7262

## ACCEL CS20



Accel Instrument Products provide the perfect balance between Safety and Efficacy by providing a safer, more responsible solution to disinfection without the use of harmful chemicals such as Aldehydes.

Accel CS 20 is designed for use in Healthcare and Personal Services when sterilization in a rapid contact time is required. Accel CS 20 achieves chemical sterilization and high level disinfection in 20 minutes and has a 14 day re-use claim.

VIROXCS20-1	1 LT
VIROXCS20-4	1 GALLON

## ACCEL CONCENTRATE



Accel Concentrate is a Ready-to-use One-step disinfectant cleaner for use as a daily cleaner and as a disinfectant of environmental surfaces and patient care devices. Accel Concentrate uses the Accelerated Hydrogen Peroxide (AHP) technology and is effective against a broad spectrum of pathogens in only 5 minutes!

VIROXACCDISC4	1 GALLON
---------------	----------

## ACCEL PREVENTION



Accel PREVENTION Hydrotherapy is designed for use as a cleaner and disinfectant for foot baths and hydrotherapy tubs where the requirement for broad-spectrum germicidal performance in a short contact time is required. Accel Prevention is the FIRST and ONLY concentrated 5 minute intermediate level disinfectant available in Canada.

VIROXACCPREXH4	1 GALLON
----------------	----------

## ACCEL TB



Accel TB is a Ready-to-use One-step disinfectant cleaner for use as a daily cleaner and as a disinfectant of environmental surfaces and patient care devices. Accel TB uses the Accelerated Hydrogen Peroxide (AHP) technology and is effective against a broad spectrum of pathogens in only 5 minutes!

### WIPES

### LIQUID

VIROXACCDISR1-TB	1 LT
VIROXACCDISR4-TB	4 LT

VIROXWIP1-TB	6X7"	160/TUB
VIROXWIP1-TB	10X10"	60/TUB
VIROXWIP1TB-SP	7.5X10"	10/PK SOFT PACK

**CAVICIDE**

METRMX1024	24 OZ SPRAY
METRMX1000	1 GALLON
METRMX1025	2.5 GALLON

**CAVIWIPES**

METRMX1100	160/TUB
METRMX1150	66/TUB (XL)
METRMX1224	45/PK FLAT PACK

**As seen on  
BC Aquaria!**

**14 DAY**

METRMX1400	3.78 LT
------------	---------

**28 DAY**

METRMX2800	3.78 LT
------------	---------

**CAVICIDE****Surface Disinfectant/Decontaminant Cleaner**

CaviCide is a convenient, ready-to-use, intermediate-level surface disinfectant which is effective against TB, HCV, bacteria (including MRSA and VRE) and fungi. It is safe for all areas of the facility including NICU, operating rooms, isolation rooms, patient care areas and laboratories. When used as directed, it will also effectively clean and decontaminate critical and semi-critical instrumentation.

CaviCide is safe for use on non-porous surfaces. It is perfect for cleaning and decontamination of environmental and medical device surfaces.

**CAVIWIPES****Low Alcohol (17.20% isopropanol)  
Surface Disinfectants**

Presaturated with CaviCide®, CaviWipes is a cleaner and disinfectant, in one. CaviWipes are durable, nonwoven, nonabrasive towelettes that offer quick, easy-to-use, time-saving convenience. They're recommended for use on non-porous surfaces and fixtures.

Ideal for operating rooms, surgical centers, neonatal units, medical and dental operatories and other critical-care areas. CaviWipes should not be disposed of in the toilet.

**METRICIDE****High-Level Disinfectant and Sterilant**

MetriCide is a 2.6% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant on immersible, lensed and delicate instruments for up to 14 days.

**METRICIDE 28****High-Level Disinfectant and Sterilant**

MetriCide 28 is a long life, 2.5% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant for up to 28 days.





## **METRIZYME**

### **Extra Strength Dual-Enzymatic Detergent**

MetriZyme is a detergent cleaner comprising of protease dual-enzymes. The active ingredients offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid soils and protein soils. The protein-dissolving enzymes in MetriZyme® allow for fast, simple removal of body fluids that often cannot be accomplished manually with soaps and other detergents.

METRMX4000

1 GALLON



## **METRISPONGE**

### **Sponge Saturated with Dual-Enzymatic Detergent**

According to 2008 SGNA guidelines, manual cleaning of endoscopes immediately after patient use is the first and most important step in removing the microbial burden prior to automated or manual disinfection.

METRMX4025

25/BX



## **EMPOWER**

### **Dual-Enzymatic Detergent**

EmPower is a fresh scented, detergent cleaner containing two protease enzymes that offers the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments. Non-corrosive to metals.

METRMX4100

1 GALLON



## **EMPOWER FOAM**

### **Foaming Enzymatic Spray**

EmPower Foam propels thick foam of protease enzymes. Effectively cleans away blood, tissue, mucous and other protein-rich body fluids. Foaming spray is ideal for use in surgery, central processing, operating room and other departments where instrument cleaning is not immediately available.

METRMX4224

24 OZ



## METRIGUARD

### General Purpose Cleaner and Disinfectant

MetriGuard is a ready-to-use, multi-purpose intermediate-level surface disinfectant which is effective against TB, HIV-1, bacteria (including MRSA and VRE) and fungi. It is safe and compatible for use on the surfaces of non-critical inanimate objects. The large 32 oz spray offers an economical advantage.

METRMX7506

32 OZ SPRAY



## METRICLEAN 2

### Industrial Cleaner, Polisher and Cart Washer Agent

MetriClean2 is a multi-purpose instrument cleaner, polisher and cart washer agent. The alkaline pH, low-foaming, concentrated cleaner is ideal for stainless steel instruments and cart washers. It enhances functionality and life of the instruments by helping to remove stains and reduce corrosion. For use in manual cleaning systems, washer sterilizers, washer

METRMX8100

3.8 LT

METRMX8150

5 GALLON

## CIDEX SOLUTIONS



- Cidex activated dialdehyde solution is a fast and effective way to disinfect and sterilize many medical devices and instruments
- Cidex 14 day use destroys 99.8% of mycobacterium buerculosis in 45 minutes at 25 degrees Celsius
- Cidex Plus (high-level disinfectant) is a 3.4% glutaraldehyde solution
- Cidex OPA (ortho-phthalaldehyde) is compatible with many medical devices, requires no activation or mixing, and destroys mycobacterium tuberculosis in 12 minutes at 25 degrees Celsius

JOJOCX250

CIDEX

3.8 LT

JOJOCX785

PLUS SOLUTION

3.8 LT

JOJO20394

OPA SOLUTION

3.8 LT

JOJO20392

OPA TEST STRIPS

60/BTL

## POCKET CPR MASKS



Mask with  
O2 & 1 way  
valve.

WASFA007691

## CPR FACE SHIELD WITH KEYCHAIN



WASFA030502

## AIRWAY GUEDELS



ALM80-8500	SIZE 00	40MM
ALM80-8510	SIZE 0	50MM
ALM80-8520	SIZE 1	60MM
ALM80-8530	SIZE 2	70MM
ALM80-8540	SIZE 3	80MM
ALM80-8550	SIZE 4	90MM
ALM80-8560	SIZE 5	100MM
ALM80-8570	SIZE 6	110MM

\*\*Please note colours may vary between brands.  
Available individually or as a complete set.

## FIRST AID KITS



WASIPF89BN26C	BASIC
WASIPF891N27C	LEVEL 1
WASIPF892N270	LEVEL 2/3

## STERI-STRIPS



3MR1540	1/8" X 3"	5/SL	50/BX
3MR1541	1/4" X 3"	3/SL	50/BX
3MR1542	1/4" X 1.5"	6/SL	50/BX
3MR1547	1/2" X 4"	6/SL	50/BX
3MR1548	1" X 5"	4/SL	25/BX

## BUTTERFLY CLOSURES



ACME1990033

LARGE 100/BX

## POLYSPORIN



UNI417337	30G CREAM
UNI418509	30G OINTMENT

## HOT/COLD COMPRESS



Reusable hot & cold packs. Can be microwaved or frozen for instant relief.

RAP12246-24	4X6"
RAP12259-24	6X9"
RAP12290-12	9X11"

## INSTANT COLD PACKS



Single use instant cold packs. Great for first aid kits, sports teams, medical offices, etc.

RAP31246-50	4X6"
RAP31259-24	6X9"

## ELASTIC BANDAGES



Elastic bandages for minor sprains or strains. Also available in latex free.

**REGULAR**

BOWBDE2	2"
BOWBDE3	3"
BOWBDE4	4"
BOWBDE6	6"

**LATEX FREE**

BOWQ532	2"
BOWQ533	3"
BOWQ534	4"
BOWQ536	6"

## FINGER GUARD



Finger guard plastic protectors help protect injured fingers and toes. Full length immobilizes the entire digit.

AMG118614	SMALL
AMG118616	MEDIUM
AMG118618	LARGE

## TRIANGULAR BANDAGES



MEDICOM5460  
40" X 40" X 56"

### HOT WATER BOTTLES



### HYDROGEN PEROXIDE



## EYE WASH STATION



PICK04  
PICK05

## EYE WASH BOTTLE & SOLUTION



PICK03	32OZ BTL
PICK04A	1 L Solution
PICK04A1	4 L Solution

## ALUMINUM MESH SPLINT

3.75"x24" Aluminum Mesh  
Splint for first aid uses.

STM80-3265-0

1EA



## SAM SPLINT

Sam Splints have the universal capability to immobilize any bone in the body, including the neck. They are durable, waterproof, reusable, lightweight, and are invisible on x-rays. They are an ideal splint for trauma first aid kits.



WASIPF6103301 46CM x 11CM

WASIPS6104301 90CM X 11CM

## PADDED FINGER SPLINT

An aluminum splint used for fractures of the phalanges.

Can be cut to the desired length with scissors and shaped.

Can be used as a baseball splint, finger protector, or finger cot.

DAL100301P

1" X 18"

DAL100303P

1/2" X 18"





## ADHESIVE BANDAGES



- Air vented, non-stick nylon pad
- Sterilized
- Individually wrapped
- 100 per box

AMG018600	3/4" X 3"	100/BX
BOWW018600	3/4" X 3"	100/BX
AMG018605	1" X 3"	100/BX

### ASSORTED SHEER

ACME1260033	100/BX
-------------	--------

## FLEXIBLE FABRIC STRIPS



ACME1595033	1 X 3"	100/BX
-------------	--------	--------

## SHEER PATCHES



ACME1307033	7/8"	100/BX
-------------	------	--------



## FABRIC KNUCKLE



ACME1602033	100/BX
-------------	--------

ACME1308033	1.5" X 1.5"	100/BX
-------------	-------------	--------

## FABRIC PATCHES



ACME1617033	2 X 3"	100/BX
-------------	--------	--------

**MICROPORE TAPE**

3M1530.5	1/2"	24/BX	AADL: M665 M665
3M15301	1"	12/BX	
3M15302	2"	6/BX	
3M15303	3"	4/BX	

**TRANSPORE TAPE**

3M1527	1/2"	24/BX
3M15271	1"	12/BX
3M15272	2"	6/BX

**PAPER TAPE**

DUKALP50	1/2"	24/BX	AADL: M665 M665
DUKALP110	1"	12/BX	
DUKALP210	2"	6/BX	
DUKALP310	3"	4/BX	

**TRANSPARENT TAPE**

DUKALT50	1/2"	24/BX
DUKALT110	1"	12/BX
DUKALT210	2"	6/BX
DUKALT310	3"	4/BX

**MEFIX TAPE**

SANC310250-05	2.5CM X 10M
SANC310500-05	5CM X 10M
SANC311000-04	10CM X 10M
SANC311500-05	15CM X 10M
SANC312000-01	20CM X 10M

**HYPAFIX TAPE**

BSN7144300	2.5CM X 10M
BSN7144301	5CM X 10M
BSN7144302	10CM X 10M
BSN7144303	15CM X 10 M



## LEUKOTAPE

Leukotape® K is an elastic tape used to support the sensory functions. The stretchable cotton material works ideally with the skin's natural elasticity. The support is coated with a skin-friendly polyacrylate adhesive. It is permeable to air reducing the risk of maceration. The wave pattern of the adhesive can support the therapeutic effect. Leukotape® K is very durable, sports and showering present no problem.

### LEUKOTAPE K

BSN7297811 5CM X 5M  
BSN7297812 7.5CM X 5M

### LEUKOTAPE P

BSN7616800 3.8CM X 13.7M

## MEPITAC TAPE

Mepitac is highly conformable and easy to apply and remove. The Safetac layer assures secure fixation and gentle removal from skin and bodily hair. If inspection or repositioning is necessary, Mepitac can be lifted and re-applied.

Mepitac is an ideal fixation solution for patients with fragile and sensitive skin. Mepitac is suitable for patients requiring repeated application and removal of tape over the same area e.g. neonates and dialysis patients.



SANC298300 2CM X 300CM  
SANC298400-00 4CM X 1.5M

## CLOTH TAPE



DUKALC50 1/2" 24/BX  
DUKALC110 1" 12/BX  
DUKALC210 2" 6/BX

## ATHLETIC TAPE



RENF20 1" 12/BX  
RENF30 1.5" 8/BX  
RENF40 2" 6/BX  
RENF50 3" 4/BX  
RENF50 4" 3/BX  
RENF101226 1.5" SECONDS 32/CS  
JAYB300X-1515 1.5" SECONDS 100/CS

## FOAM UNDERWRAP



JAYB50-27530 48/CS

**GAUZE****NON-STERILE**

DUKAL6112	2"X2"	200/BG	
DUKAL6113	3"X3"	200/BG	
DUKAL6114	4"X4"	200/BG	AADL: M715

**STERILE**

DUKAL6122	2"X2"	2/PK	25PK/BX	
DUKAL6123	3"X3"	2/PK	25PK/BX	
DUKAL6124	4"X4"	2/PK	25PK/BX	AADL: M717

**TELFA NON-ADHERENT  
PADS**

KEND1961	2"X3"	100/BX
KEND2132	3"X4"	100/BX
KEND2138	3"X8"	50/BX

**ABD PADS**

KEND9190A	5"X9"	36/TRAY
KEND9194A	8"X10"	18/TRAY

AADL: M730

**COMBINE ROLL**

DUKAL3001 8"X20YDS

**CONFORMING GAUZE**

DUKAL602D	2" CONFORM	12/BG
DUKAL603D	3" CONFORM	12/BG
DUKAL604D	4" CONFORM	12/BG
DUKAL606D	6" CONFORM	6/BG

AADL: M246



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687

## 3M™ COBAN™ SELF-ADHERENT WRAP



5 YDS/RL

A self-adherent elastic wrap that functions like a tape, but sticks only to itself. Available in sterile and nonsterile styles, and in a variety of widths and colors to meet your application needs.

3M1581	1"	30/CS
3M1582	2"	36/CS
3M1583	3"	24/CS
3M1584	4"	18/CS
3M1586	6"	12/CS

## CO-PLUS®/CO-PLUS® LF ELASTIC COHESIVE BANDAGE



Co-Plus® flexible elastic cohesive bandage provides non-slip light support and compression. Because of its unique construction, Co-Plus® adheres to itself and not to skin, aiding in patient comfort. No clips or fasteners required. Available in latex and latex free versions.

### Features and Benefits

- Adheres only to itself
- Non-slip support
- Breathable
- Hand-tearable (Co-Plus® LF only)
- Latex Free (Co-Plus® LF only)

BSN7210014	2.5 cm x 4.5 m* Beige	60 RL/BX
BSN7210016	5 cm x 4.5 m* Beige	36 RL/BX
BSN7210017	7.5 cm x 4.5 m* Beige	24 RL/BX
BSN7210018	10 cm x 4.5 m* Beige	18 RL/BX
BSN7210019	15 cm x 4.5 m* Beige	12 RL/BX
BSN7210033	5 cm x 4.5 m* Mixed Colors	36 RL/BX



## COHESIVE BANDAGES

Our cohesive bandage maintains consistent non-slip compression, yet is soft, flexible, and porous. It sticks only to itself and never leaves residue. Utilizing woven fabric, Dukal Cohesive Bandage tears by hand in both directions. An excellent product when additional support is needed. This durable self adhesive wrap is well suited for use in a variety of applications and comes in a many colors and sizes.



*This product contains latex.*

DUKAL8015AS	1" Non Sterile Cohesive Bandage, Assorted	30 Rolls per CS
DUKAL8015DB	1" Non Sterile Cohesive Bandage, Dark Blue	30 Rolls per CS
DUKAL8015T	1" Non Sterile Cohesive Bandage, Tan	30 Rolls per CS
DUKAL8015W	1" Non Sterile Cohesive Bandage, White	30 Rolls per CS
DUKAL8155AS	1.5" Non Sterile Cohesive Bandage, Assorted	48 Rolls per CS
DUKAL8155DB	1.5" Non Sterile Cohesive Bandage, Dark Blue	48 Rolls per CS
DUKAL8155R	1.5" Non Sterile Cohesive Bandage, Red	48 Rolls per CS
DUKAL8155T	1.5" Non Sterile Cohesive Bandage, Tan	48 Rolls per CS
DUKAL8155W	1.5" Non Sterile Cohesive Bandage, White	48 Rolls per CS
DUKAL8025AS	2" Non Sterile Cohesive Bandage, Assorted	36 Rolls per CS
DUKAL8025DB	2" Non Sterile Cohesive Bandage, Dark Blue	36 Rolls per CS
DUKAL8025T	2" Non Sterile Cohesive Bandage, Tan	36 Rolls per CS
DUKAL8035AS	3" Non Sterile Cohesive Bandage, Assorted	24 Rolls per CS
DUKAL8035DB	3" Non Sterile Cohesive Bandage, Dark Blue	24 Rolls per CS
DUKAL8035T	3" Non Sterile Cohesive Bandage, Tan	24 Rolls per CS
DUKAL8045T	4" x 5yds Non Sterile Cohesive Bandage, Tan	18 Rolls per CS
DUKAL8045AS	4" x 5yds" Non Sterile Cohesive Bandage, Assort	18 Rolls per CS

**MESALT**

Mesalt is intended for the management of heavily discharging and discharging infected wounds in the inflammatory phase and deep cavity wounds such as pressure sores and surgical

SANC285580	5CM X 5CM	30/BX
SANC285780	7.5CM X 7.5CM	30/BX
SANC286080-11	10CM X 10CM	30/BX

**MESALT RIBBON**

SANC285280-11	2CM X 1M	10/BX
---------------	----------	-------

**ALLDRESS**

Alldress is an all-in-one self-adherent absorbent dressing ideal for a wide variety of low to moderately exuding wounds as either primary or secondary dressing. The backing film provides a moist wound environment and protects the wound from water and contamination.

SANC265320-14	10CM X 10CM	10/BX
SANC265340	15CM X 15CM	10/BX
SANC265360-14	15CM X 20CM	10/BX

**MEPORE**

Mepore is the ideal dressing for surgical wounds, and can be used for all wounds that have low to medium exudate levels.

The absorbent pad of Mepore has a wound contact layer, which is designed to give minimum adherence to the wound.

SANC670800	6CM X 7CM	60/BX
SANC670900	9CM X 10CM	50/BX
SANC671000	9CM X 15CM	50/BX
SANC671100	9CM X 20CM	30/BX
SANC671200-17	9CM X 25CM	30/BX
SANC671300	9CM X 30CM	30/BX

**MEPORE PRO**

Showerproof, self-adhesive absorbent dressing

Mepore Pro offers all the advantages of a Mepore dressing, but is also showerproof. Comfortable and easy to work with, Mepore Pro offers gentle yet secure protection while minimizing inconvenience to patients.

SANC670820	6CM X 7CM	60/BX
SANC670920	9CM X 10CM	50/BX
SANC671020	9CM X 15CM	50/BX
SANC671120	9CM X 20CM	30/BX
SANC671220-18	9CM X 25CM	30/BX
SANC671320	9CM X 30CM	30/BX

## MEPILEX



Mepilex is a soft and highly conformable foam dressing that absorbs exudate and maintains a moist wound environment. The Safetac layer seals the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk for maceration. The Safetac layer ensures that the dressing can be changed without damaging the wound or surrounding skin or exposing the patient to additional pain.

SANC294100	10CM X 10CM	5/BX
SANC294200	10CM X 20CM	5/BX
SANC294300	15CM X 15CM	5/BX
SANC294400	20CM X 20CM	5/BX


## MEPILEX LITE



Mepilex Lite is a thin and highly conformable foam dressing that absorbs exudate and maintains a moist wound environment. The Safetac layer seals around the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk for maceration. The Safetac layer ensures atraumatic dressing changes. Mepilex Lite can be cut to suit various wound shapes and locations.

SANC284000	6CM X 8.5CM	5/BX
SANC284100	10CM X 10CM	5/BX
SANC284300	15CM X 15CM	5/BX
SANC284500	20CM X 50CM	2/BX

## MEPILEX AG



Mepilex Ag combines the unique features of Safetac technology with the bacteria reducing power of silver. Mepilex Ag goes to work quickly, inactivating wound pathogens within 30 minutes and for up to 7 days. At dressing removal, Mepilex Ag does not stick to the wound or strip surrounding skin, minimizing patient pain and wound trauma.

SANC287110	10CM X 10CM	5/BX
SANC287210	10CM X 20CM	5/BX
SANC287310	15CM X 15CM	5/BX
SANC287410	20CM X 20CM	5/BX

## MEPILEX HEEL



Mepilex Heel minimizes trauma to the wound and pain to the patient at dressing changes. The Safetac technology prevents Mepilex Heel from sticking to the wound bed yet adheres gently to the surrounding skin allowing easy application of secondary fixation. Mepilex Heel absorbs effectively and ensures a low risk of maceration. Mepilex Heel conforms well to the heel preventing exudate from pooling and causing maceration. Mepilex Heel is ideal for most shapes or sizes of feet.

SANC288100	13CM X 20CM	5/BX
------------	-------------	------

## MEPILEX TRANSFER



Mepilex Transfer is thin and conformable making it easy to keep the dressing in contact with the wound surface and the surrounding skin, even in awkward or uneven areas. The structure of Mepilex Transfer allows the exudate to be transferred vertically into a secondary absorbent pad. The Safetac layer seals around the wound edges preventing the exudate to leak horizontally on to the surrounding skin, thus, minimising maceration and excoriation. Mepilex Transfer maintains a moist wound environment in combination with an appropriate secondary dressing.

SANC294800	15CM X 20CM	5/BX
SANC294502	20CM X 50CM	5/BX
SANC294600	7.5CM X 8.5CM	5/BX
SANC294700	10CM X 12CM	5/BX

## MEPILEX BORDER



Mepilex Border is an all-in-one foam dressing that effectively absorbs and retains exudate and maintains a moist wound environment. The Safetac layer seals the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk for maceration. The Safetac layer ensures that the dressing can be changed without damaging the wound or surrounding skin or exposing the patient to additional pain.

SANC295200	7.5CM X 7.5CM	5/BX
SANC295300	10CM X 10CM	5/BX
SANC295000	12.5CM X 12.5CM	5/BX
SANC295400	15CM X 15CM	5/BX
SANC295600	15CM X 20CM	5/BX
SANC295800	10CM X 20CM	5/BX
SANC295900	10CM X 30CM	5/BX

## MEPILEX BORDER-LITE



Mepilex Border Lite minimizes trauma to the wound and pain to the patient at dressing changes. The Safetac technology prevents Mepilex Border Lite from sticking to the wound bed yet adheres gently to the surrounding skin. Mepilex Border Lite is ideal for those locations where a thin self-adherent dressing is required. It is a very conformable and gentle alternative to other thin adhesive bordered dressings.

SANC281000	4CM X 5CM	10/BX
SANC281100	5CM X 12.5CM	5/BX
SANC281200	7.5CM X 7.5CM	5/BX
SANC281300	10CM X 10CM	5/BX
SANC281500	15CM X 15CM	5/BX

## MELGISORB



Melgisorb® can be used for a wide range of moderately to heavily exuding partial to full thickness wounds such as pressure sores, venous and arterial ulcers, diabetic ulcers, donor sites and post-operative wounds. Infected wounds should be evaluated on a regular basis and treated appropriately.

SANC250600	5CM X 5CM	10/BX
SANC251100	10CM X 10CM	10/BX

## MELGISORB AG



Benefits of Melgisorb Ag

- Excellent absorption for heavily exuding deep wounds
- Provides sustained antimicrobial effect
- High fiber strength allows removal in one piece

SANC256100-00	10CM X 10CM	10/BX
---------------	-------------	-------

## MESORB



Mesorb® is designed for direct use on heavily discharging wounds and open drains - without any need for supplementary dressings in between.

Discharge and exudate is drained and transferred into Mesorb by the wound contact layer and the diffusion layer. The highly-absorbent core absorbs exudate, while simultaneously providing a good protective and ventilating wound cushion. The protective backing on the reverse side of the absorbent core prevents clothes and bedlinen from being contaminated. As the backing covers the sides of the dressing, it also prevents side leakage.

SANC677001	10CM X 10CM	50/BX
SANC677201	10CM X 15CM	50/BX
SANC677401	10CM X 20CM	50/BX
SANC677601	15CM X 20CM	50/BX
SANC677701	20CM X 25CM	30/BX
SANC677801	30CM X 20CM	30/BX

## MEPITEL



Mepitel is a thin, semi-transparent dressing that will conform to the shape of your wound and should be used in combination with an outer absorbent dressing.

Containing the unique and patented Safetac® adhesive on both sides of the dressing, Mepitel prevents the outer absorbent dressing from sticking to your wound (or the skin around your wound). This makes the dressing changes as pain free as possible and also makes sure that the skin around your wound stays dry and healthy.

Mepitel can be left in place for up to 14 days, depending on the condition of your wound, which means that only the outer dressing needs replacing during the 14 day period as and when needed.

SANC290510	5CM X 7.5CM	10/BX
SANC290710	7.5CM X 10CM	10/BX
SANC291010-15	10CM X 18CM	10/BX
SANC292005	20CM X 30CM	5/BX



## ALLEVYN ADHESIVE



ALLEVYN ◊ Adhesive - Hydrocellular Dressing

- Low allergy adhesive
- Tapered construction with micro-thin edges
- Unique trilaminate structure
- Perforated wound contact layer
- Highly absorbent central hydrocellular layer
- Waterproof outer film layer

SMTH66000043	7.5CM X 7.5CM	10/BX
SMTH66000044	12.5CM X 12.5CM	10/BX

## ALLEVYN NON-ADHESIVE



ALLEVYN◊ Hydrocellular Foam Dressing

- Unique trilaminate structure
- Non adherent wound contact layer
- Highly absorbent central hydrocellular layer
- Waterproof outer film layer
- Non adhesive

SMTH66007643	5CM X 5CM	10/BX
SMTH66007637	10CM X 10CM	10/BX
SMTH66000093	15CM X 15CM	10/BX

## ALLEVYN GENTLE



ALLEVYN◊ Gentle Gel Adhesive Hydrocellular Foam Dressing

- Advanced triple-layered construction
- Absorbent hydrocellular pad
- Perforated soft gel adhesive wound contact layer
- Highly permeable waterproof outer film

SMTH66800247	5CM X 5CM	10/BX
SMTH66800248	10CM X 10CM	10/BX

## ALLEVYN GENTLE BORDER



ALLEVYN ◊ Gentle Border Gel Adhesive Hydrocellular Foam Dressing with Border

SMTH66800269	7.5CM X 7.5CM	10/BX
SMTH66800270	10CM X 10CM	10/BX
SMTH66800272	12.5CM X 12.5CM	10/BX

## INTRASITE CONFORMABLE



IntraSite\* Conformable is a conformable hydrogel dressing which combines the advantages of IntraSite Gel with a non-woven dressing to aid in the gentle packing of deep, open or undermined wounds. The non-woven substrate ensures the gel remains in close contact with the entire wound surface. It is suitable for use in necrotic, sloughy and granulating wounds, and helps create a moist wound environment to facilitate healing.

SMTH66000324	10CM X 10CM	10/BX
SMTH66000325	10CM X 20CM	10/BX
SMTH66000326	10CM X 40CM	10/BX

## INTRASITE GEL



SMTH7308	8G
SMTH7311	15G
SMTH7313	25G

## ICHTHOPASTE



SMTH4959 7.5CM X 6M

## ACTICOAT



Acticoat\* (with SILCRYST\*\* Nanocrystals) Dressings are an effective barrier to bacterial penetration. The barrier function of the dressing may help reduce infection in partial and full-thickness wounds including pressure ulcers, venous ulcers, diabetic ulcers, surgical wounds, first and second degree burns and graft and dermal substitutes recipient sites. Acticoat\* dressings may be used over debrided and grafted partial thickness wounds.

### ACTICOAT

SMTH20101	4" X 4"	12/BX
-----------	---------	-------

### ACTICOAT 7

SMTH20141	4" X 5"	5/BX
SMTH20241	6" X 6"	5/BX

## ACTICOAT FLEX



SMTH66800396	5CM X 5CM	5/BX
SMTH66800399	10CM X 10CM	12/BX
SMTH66800409	10CM X 20CM	12/BX
SMTH66800419	20CM X 40CM	5/BX
SMTH66800432	40CM X 40CM	6/BX

## JELONET



### Indications

- Minor burns and scalds
- Donor and recipient graft sites
- Skin loss wounds, lacerations and abrasions
- Leg ulcers

SMTH7403	5CM X 5CM	50/BX
SMTH7404	10CM X 10CM	10/BX
SMTH7477	10CM X 40CM	10/BX

## BACTIGRAS



Bactigras is indicated for adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers.

SMTH7457	10CM X 10CM	10/BX
SMTH7461	15CM X 20CM	10/BX

## IODOSORB



### PASTE

SMTH66060633	4CM X 6CM	5G	5/BX
SMTH66060634	6CM X 8CM	10G	3/BX

### OINTMENT

SMTH66060630	10G TUBE	4/BX
SMTH66060631	20G TUBE	2/BX
SMTH66060632	40G TUBE	1/BX

## INADINE



INADINE® PVP-I Non-Adherent Dressing consists of a low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base containing 10% Povidone Iodine; equivalent to 1.0% available iodine.

JOJOPO1481	5CM X 5CM	25/BX
JOJOPO1512	9.5CM X 9.5CM	25/BX

## DRESSING TRAYS



- 1 Ea - Waste Bag W/ Tie
- 1 Ea - Drape Towel 45x55cm
- 1 Ea - Safety Pin
- 1 Ea - 3 Compartment Tray
- 1 Ea - Plastic Transfer Fcp
- 2 Ea - Metal Forceps
- 8 Ea - Gauze 10x10cm
- 8 Ea - Gauze 5x5cm

BOWTRDR101278B 40/CS



- 1 Ea - Drape Towel 45x55cm
- 2 Ea - Gauze 10x10cm
- 1 Ea - Waste Bag W/ Tie
- 1 Ea - 3 Compartment Tray
- 1 Ea - Plastic Transfer Fcp
- 6 Ea - Gauze 5x5cm
- 1 Ea - Plastic Locking Fcp

BOWTRDR6588 30/CS

**\*\*Custom trays are available upon request**

## TEGADERM DRESSING



A transparent film dressing designed for protecting skin and wound sites and securing devices to skin. It consists of a patterned coated film using Comfort Adhesive Technology for gentle adhesion with increased breathability and moisture exchange.

3M1622W	4.4CM X 4.4CM	100/BX
3M1624W	6CM X 7CM	100/BX
3M1626W	10CM X 12CM	50/BX
3M1627W	4CM X 10CM	20/BX

## TEGADERM ROLL



3M16002	2" X 11YD	1 EA
3M16004	4" X 11YD	1 EA

## DUODERM EXTRA THIN



Thin polyurethane film provides a bacterial<sup>1</sup> and viral<sup>2</sup> barrier. Film provides a waterproof barrier over the dressing  
Can be removed without damaging newly formed tissue  
Designed to reduce the risk of further skin breakdown due to friction by preventing contact with clothes/bed linen

CONV187955	4"X4"	10/BX
AADL: M776		

## ADAPTIC



ADAPTIC® Non-Adhering Dressing is a primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion.

It is designed to help protect the wound while preventing the dressing from adhering to the wound<sup>1</sup> and to minimize pain and trauma upon removal.

JOJO2012	3"X3"	50/BX
----------	-------	-------

## AQUACEL



Provides excellent absorption and retention capabilities for moderate to highly exuding wounds<sup>2</sup>

Conforms to the wound surface to form an intimate contact,<sup>3</sup> Helps reduce wound pain while the dressing is in situ and upon removal.<sup>4-7</sup>

Supports wound healing by providing a moist wound healing environment

May help minimize cross-infection during dressing removal.<sup>1</sup>

CONV177902	4"X4"	10/BX
CONV177901-CW	2"X2"	10/BX

## PROMOGRAN®


PROMOGRAN® is sterile, freeze-dried composite of 55% collagen and 45% oxidised regenerated cellulose (ORC). PROMOGRAN® is designed to promote an optimal healing environment<sup>1</sup>

With PROMOGRAN® stalled wounds close:

- Faster 1,2
- Cost Effectively 1,3,4

In the presence of exudate the PROMOGRAN® protease modulating matrix transforms into a soft and conformable, biodegradable gel: this allows contact with all areas of the wound.

JOJO772028	28CM SQ	10/BX
JOJO772123	123CM SQ	10/BX



## PROMOGRAN PRISMA®


PROMOGRAN PRISMA® matrix is a freeze-dried composite of collagen, oxidised regenerated cellulose (ORC) and silver. Made from a matrix of 55% collagen, 44% ORC and 1% ORC/Silver.

PROMOGRAN PRISMA® is designed to 'kick start' the healing process while providing protection from infection.<sup>1-2</sup>

In the presence of exudate PROMOGRAN PRISMA® matrix transforms into a soft and conformable, biodegradable gel; this allows contact with all areas of the wound. With PROMOGRAN PRISMA® stalled wounds heal:

- Faster<sup>3-5</sup>
- Cost Effectively<sup>3,6</sup>
- While providing protection from infection<sup>7-8</sup>

JOJOPS2028	28CM SQ	10/BX
JOJOPS2123	123CM SQ	10/BX





## SILVERCEL



SILVERCEL® Antimicrobial Dressing is intended for use in the management of all moderate to heavily exudating partial and full-thickness chronic wounds.

JOJOCAD230	2.5CM X 30.5CM	5/BX
JOJOCAD050	5CM X 5CM	10/BX
JOJOCAD011	11CM X 11CM	10/BX
JOJOCAD020	10CM X 20CM	5/BX

## ACTISORB SILVER



ACTISORB® Silver 220 dressing is suitable as the first therapeutic step in the management of all chronic wounds. It is indicated for: fungating carcinomas, ulcerative, traumatic and surgical wounds where bacterial contamination, infection or odour occurs.

JOJOMAS065	6.5CM X 9.5CM	10/BX
JOJOMAS105	10.5CM X 10.5CM	10/BX
JOJOMAS190	19CM X 10.5CM	10/BX

## AMD PACKING STRIP



Curity™ AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA. This exciting new product provides healthcare facilities a low cost, prophylactic method of treatment without changing nursing protocol.

KEND7831AMD	1/4" X 1YD
KEND7832AMD	1/2" X 1YD
KEND7833AMD	1" X 1YD

## CICA-CARE



Cica-Care is a patented technically advanced form of silicone gel sheet for scar care. Cica-Care is designed for use in the management of both existing and new hypertrophic and keloid scars (red and raised). Also as a preventative therapy on closed wounds to prevent hypertrophic and keloid scars (red and raised).

SMTH66320703	12CM X 6CM	1 SHEET
SMTH66250706	12CM X 15CM	1 SHEET



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687

**SURGIFIX**

SIZE 1	FRA1
SIZE 2	FRA2
SIZE 3	FRA3
SIZE 4	FRA4
SIZE 5	FRA5
SIZE 6	FRA6
SIZE 7	FRA7
SIZE 8	FRA8
SIZE 9	FRA9
SIZE 10	FRA10

**TUBIGRIP**

SIZE B	SANC1436
SIZE C	SANC1443
SIZE D	SANC1437
SIZE E	SANC1434
SIZE F	SANC1438
SIZE G	SANC1439

**SURGITUBE**

1 5/8"	GL105W
5/8"	GL110W

**STOCKINETTE**

2"	WEBSTO2
3"	WEBSTO3
4"	WEBSTO4
6"	WEBSTO6

**OPSITE FLEXIFIX**

The film is waterproof and a barrier to bacteria. Its high moisture vapour permeability allows the skin to breathe. The acrylic adhesive minimises the risk of skin damage on removal, especially after long periods of wear.

SMTH66000040	5CM X 10M
SMTH66000041	10CM X 10M

**OPSITE FLEXIGRID**

OpSite Flexigrid is a transparent, adhesive film, with a unique wound mapping grid. The film is moisture vapour permeable, conformable and extensible. It is clinically proven to promote moist wound healing for superficial wounds; it is also ideal for use as a secondary dressing and to help prevent pressure sores.

**FLEXIGRID**

SMTH4630	10CM X 12CM	50/BX
----------	-------------	-------

**FLEXIGRID IV**

SMTH4628	6CM X 7CM	100/BX
SMTH4631	15CM X 20CM	10/BX

## GYPSONA

The traditional rich and creamy texture of Gypsona® allows intimate molding and creates a smooth, appealing finished cast. The unique interlocking leno weave greatly enhances stability allowing more plaster to stay on the gauze. The result is a high quality premium plaster of Paris making Gypsona® the world's leading brand.

### Features and Benefits

- Leno gauze for extra stability
- High mass weight resulting in low plaster loss
- Creamy texture
- Central plastic core
- Latex Free



### GYPSONA® - PLASTER BANDAGES - EXTRA FAST SETTING / 2 MINUTES

BSNLPL23	5 CM X 2.7 M	24 RL/BX	3 BX/CS
BSNLPL33	7.5 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL43	10 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL63	15 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL83	20 CM X 2.7 M	12 RL/BX	4 BX/CS
BSNLPL45	10 CM X 4.5 M	12 RL/BX	6 BX/CS
BSNLPL65	15 CM X 4.5 M	12 RL/BX	4 BX/CS
BSNLPL85	20 CM X 4.5 M	12 RL/BX	4 BX/CS

### GYPSONA® - PLASTER SLABS - EXTRA FAST SETTING / 2 MINUTES

BSNL5430	7.5 CM X 75 CM	50 SPLINTS	6 BX/CS
BSNL5530	12.5 CM X 75 CM	50 SPLINTS	6 BX/CS

### GYPSONA® - PLASTER BANDAGES - FAST SETTING / 4 MINUTES

BSN4000	5 CM X 2.7 M	24 RL/BX	3 BX/CS
BSN4010	7.5 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4040	10 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4045	15 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4020	10 CM X 4.5 M	12 RL/BX	6 BX/CS
BSN4030	15 CM X 4.5 M	12 RL/BX	4 BX/CS
BSN4050	20 CM X 4.5 M	12 RL/BX	4 BX/CS

**PROFORE**

The multi-layer compression bandage system has been specifically designed for the management and treatment of venous leg ulcers and associated conditions. The system pack can be used on patients with ankle circumferences of greater than 18cm (padded), but will give best results on ankle circumferences of 18-25cm.

**PROFORE LITE**

SMTH66000016  
SMTH66000415

PROFORE 4 LAYER  
PROFORE LITE 3 LAYER

**SUREPRESS COMPRESSION BANDAGE**

CONV650947

SurePress® High Compression Bandage is an elastic wrap that provides sustained graduated compression in an easy to use format. This reusable bandage has a yellow center line and two rectangular extension indicators for easier application and appropriate amount of compression for a variety of leg sizes. It is indicated for the management of venous leg ulcers and associated conditions. It can be used alone or in combination with SurePress® Absorbent Padding.

**SUREPRESS PADDING**

CONV650948

**COMPRILAN**

Comprilan is made of 100% cotton and is therefore highly skin friendly. This short stretch bandage is durable and extremely sturdy. Comprilan can be washed many times at 40°C without losing its specific features. The bandage ends do not fray and the interlocking weave provides air permeability. Support of the muscle pump for elimination of lymphatic stasis

BSN7718700 4 CM x 5 M  
BDF0102600 6 CM x 5 M  
BDF0102700 8 CM x 5 M  
BDF0102800 10 CM x 5 M  
BDF0102900 12 CM x 5 M



## Support Therapy Features

Provides relief from a range of common symptoms:

- Tired, aching legs
- Swollen feet and ankles
- Pregnancy
- Wound care
- Lymphedema
- Travelling
- Varicose and spider veins
- Infectious disease
- Sedentary lifestyle
- Surgery or trauma



## Styles Include:

- Athletic
  - Casual Cotton
  - Dress Socks (for Men & Women)
  - Relief
  - Sheer & Ultra Sheer
  - Fashionable dress
- Also available in a variety of colours and some styles offer the options of open toe or closed toe.

Available in a variety of compressions.  
Bowers also does custom measurements and fittings  
Some styles available in *Latex Free* materials

## ACCESSORIES



### Rubber Gloves

Practical and safe way to apply and remove compression stockings



### Footslip

Applicator for open toe compression stockings



### Stocking donner

To aid the process of putting on compressions stockings



### Doff N' Donner

The revolutionary donning device for all compression socks and stockings

**\*\*Please note that all stockings require measurements and a fitting, which are done in the morning before 11 am, to ensure proper sizing.**



## SIGVARIS PERFORMANCE SOCKS

A growing trend among competitive athletes and weekend warriors is the use of graduated compression socks to help improve their athletic performance and recovery. As the global leader in graduated compression therapy, SIGVARIS creates truly superior sports products featuring true graduated compression for athletes and active people. Athletes who wear SIGVARIS Sports Socks and Sleeves will benefit from:

- Increased circulation and improved blood flow
- Increased oxygen to muscle tissue
- Reduced pulled muscles and less exercised-induced muscle soreness in the legs
- Less lactic acid buildup in the legs during exercise
- Improved energy, performance and recovery
- Injury prevention and reduced shin splints & calf cramps

### ***Peak Your Performance***

- Improved blood circulation and oxygen delivery
- Less aching & cramping
- Thermal & odor control
- Achilles tendon protector reduces vibration on ligaments



### ***Available Styles & Colours***

- Men & Women
- Black or White
- Closed Toe
- 20-30mmhg

## MERINO OUTDOOR PERFORMANCE SOCKS

A growing trend among competitive athletes and weekend warriors is the use of graduated compression socks to help improve their athletic performance and recovery. As the global leader in graduated compression therapy, SIGVARIS creates truly superior sports products featuring true graduated compression for athletes and active people. Athletes who wear SIGVARIS Sports Socks and Sleeves will benefit from:

- Increased circulation and improved blood flow
- Increased oxygen to muscle tissue
- Reduced pulled muscles and less exercised-induced muscle soreness in the legs
- Less lactic acid buildup in the legs during exercise
- Improved energy, performance and recovery
- Injury prevention and reduced shin splints & calf cramps

### ***Peak Your Performance***

- 20–30mmHg graduated compression optimizes athletic performance
- Fine Australian Merino Wool Socks constructed with itch-free fibers feels soft and comfortable against the skin
- Naturally durable, anti-microbial and thermo-regulating for all seasons
- Absorbs 35% of its weight in water, keeping skin dry longer
- Cushioned padded sole keeps feet blister-free and mesh flex-zones provide ventilation



### ***Available Styles & Colours***

- Men & Women
- Charcoal or Olive
- Closed Toe
- 15-20mmhg
- 20-30mmhg



## JOBST ACTIVEWEAR

JOBST® ActiveWear is effective leg therapy in an energizing athletic sock, ideal for various activities, from walking to playing sports.

- The ideal combination of therapeutic efficacy and Dri-release® yarn for superior moisture control
- Fits an active lifestyle with a seamless toe and 360° cushioned foot for increased wearing comfort
- All-day comfort knee band that keeps socks up without binding or pinching
- Reinforced heel resists abrasion and provides greater durability

### ***Available Styles & Colours***

Men & Women

Black or White

Closed Toe

20-30mmhg & 30-40mmhg



- A high performance sock that's durable enough for the most active lifestyles
- Lightweight construction
- Knit-in heel and ample toe are ensure extra comfort and a perfect fit
- Nylon/acrylic blend helps wick away moisture
- Sporty design looks like a regular athletic sock

### ***Available Styles & Colours***

Men & Women

White

Closed Toe

15-20mmhg

## PERFORMANCE SLEEVES

A growing trend among competitive athletes and weekend warriors is the use of graduated compression socks to help improve their athletic performance and recovery. As the global leader in graduated compression therapy, SIGVARIS creates truly superior sports products featuring true graduated compression for athletes and active people. Athletes who wear SIGVARIS Sports Socks and Sleeves will benefit from:

- Increased circulation and improved blood flow
- Increased oxygen to muscle tissue
- Reduced pulled muscles and less exercised-induced muscle soreness in the legs
- Less lactic acid buildup in the legs during exercise
- Improved energy, performance and recovery
- Injury prevention and reduced shin splints & calf cramps

### ***Peak Your Performance***

- Extra-wide comfort top band keeps sleeve comfortably in place
- Less aching & cramping to optimize performance potential
- Stabilizes & diminishes vibrations of active muscles
- Ideal for triathletes—fabric doesn't deteriorate in salt or chlorine water



Available colours:



## BANDIT® TENNIS ELBOW STRAP

Helps relieve pain without constricting movement

- Latex-free.
- Useful for treating medial or lateral epicondylitis.
- Unique design compresses across the flexor, extensor, and supinator muscles.
- Allows movement without compromising circulation.
- Hand wash in cool, soapy water. Air dry.
- One adjustable size fits forearms with a minimum circumference of 8" (20 cm).



REMRM15890

## ANATECH TENNIS ELBOW STRAP



Available in sizes XS-XL

3mm AnaPrene elbow strap with Velcro® closure

- Pressure pad applies localized pressure over affected tendons
- Manufactured in AnaPrene for warmth and compression
- Velcro® cinching strap for extra adjustable pressure
- Optional extra pad available, if both tendons affected (Model #5100 - sold separately)

## SEROLA TENNIS ELBOW SUPPORT

The Serola Gel Arc Elbow Brace is superior in elbow support, designed to stop pain more effectively with:

Dual-Sided Gel Arc - concentrated or diffused pressure for maximum relief

Secondary Attachment Point - significantly reduces pull at epicondyle and allows better muscle function

Location - the arc affects more muscles and tendons as they converge near the epicondyle

Absorption - the gel is a superior barrier in absorbing muscle pull, shock, and vibration.

Foam Pad - hypoallergenic pad under buckle enhance comfort, reduces irritation

Uniquely interchangeable for right or left arm, lateral or medial epicondyle.



ONE SIZE FITS ALL



## ANAPRENE WRIST SUPPORT



- Manufactured in 3mm breathable neoprene for firm support and warmth
- Removable malleable palmar and dorsal aluminum bars
- Adjustable Velcro® closure
- Low profile design allows full finger movement
- Available in beige, black or royal blue
- Specify left or right

Available in sizes XS-XL

## ABDOMINAL BINDER

- Versatile, use as an ab binder, a back support, a rib belt or a post-operative belt
- Two heights available to accommodate most patients' needs
- Comfortable, multi-tiered plush elastic contours easily and helps prevent rolling
- CorEdge® finish for comfort
- Provides gentle compression and light support to abdominal area without restricting expansion
- Available in 9" or 12" widths (6" custom)
- Retail packaging
- Available in white only



Standard Size

\*\*Also available in 2XL & 3XL

## SOFT CERVICAL COLLAR



Soft cervical collars for support of the head and neck.

BOWCVXSS	EXTRA SMALL
BOWCVSS	SMALL
BOWCVSM	MEDIUM
BOWCVSL	LARGE
BOWCVXS	EXTRA LARGE



**BOWERS**  
HOME HEALTH

[www.bowersmedical.com](http://www.bowersmedical.com)

Richmond: 1-800-663-0047  
Edmonton: 1-800-561-3687

## ANAPRENE KNEE SLEEVE

- Manufactured in lightweight, breathable 2mm AnaPrene with a patella hole
- Made 19½" long for use under a functional knee brace
- May be specially ordered in 3mm or 5mm AnaPrene or in a different length to accommodate a specific brace
- Custom buttress options and extra features available upon request
- Available in beige, black, red or navy (other thicknesses may not have same colour choices)
- Custom sizing available upon request



Available in sizes XS-XXL

## HERNIA AID

Designed to eliminate discomfort associated with common reducible inguinal hernias Anatomically shaped pads apply comfortable direct pressure



UNI02280196

SMALL

UNI02262376

MEDIUM

UNI02264240

LARGE

## ARM SLING

Canvas arm sling with plastic snap hooks.  
Universal size.



DALCO126730P UNIVERSAL

## ARM SLING WITH IMMOBILIZING STRAP

Arm sling with immobilizing swathe strap, thumb loop, and contact closures.



DALCO73432

SMALL

DALCO73433

MEDIUM

## OBUS FORME BACK SUPPORT

Designed to promote proper posture and good body mechanics to help maintain a healthy back, reduce back pain and prevent back pain injuries due to lifting, pushing, pulling, walking, exercising and repetitive bending or twisting. The ObusForme Back Belt comes in four sizes specifically tailored to the shape of the male or female physique.

### FEATURES

- Conforms and supports the lower back, promoting muscle relaxation and proper alignment of the back
- Increases the flexibility of the stomach and back muscles by keeping them warm
- Provides intra-abdominal support to help support the spine take weight off the lower back
- Built-in, ventilated lumbar pad provides additional support to the lumbar area
- Breathable nylon mesh wrap with adjustable hook and loop front closure holds the belt in place
- Lightweight; Sleek, stylish, low profile. Wear on top of or under clothing
- Designed for industrial and/or personal usage; any profession and task where proper posture and good body mechanics are important
- Ideal for mild low back pain, strains, or sprains
- Beneficial when used in conjunction with a rehabilitation program
- All materials are hypoallergenic; hand-wash

### FEMALE BACK SUPPORT



### MALE BACK SUPPORT



Small-Medium: S-M (26"-32" / 66-81 cm)  
 Medium-Large: M-L (30"-36" / 76 cm-92cm)  
 Large-Extra Large: L-XL (37"-43" / 94-109cm)  
 Extra Extra Large: XXL (44"-60" / 112-152cm)

#### FEMALE SIZES

UNI00423632	SMALL/MEDIUM
UNI00423640	MEDIUM/LARGE
UNI00423673	LARGE/XLARGE

#### MALE SIZES

UNI00423681	SMALL/MEDIUM
UNI00423699	MEDIUM/LARGE
UNI00423715	LARGE/XLARGE

## THERA-BAND® EXERCISE BANDS

Thera-Band® latex Exercise Bands are available in 8 color-coded levels of resistance.

- Improve strength, range of motion and cooperation of muscle groups
- Color-coded progressive resistance gives at-a-glance documentation of progress from one level to the next
- Resistance levels are determined by thickness of the band, increasing in difficulty from tan (Extra Thin) to gold (Max)

**\*\*\*Also available in latex free**

REMTH-20110	TAN - EXTRA THIN	50YD
REMTH-20120	YELLOW - THIN	50YD
REMTH-20130	RED - MEDIUM	50YD
REMTH-20140	GREEN - HEAVY	50YD
REMTH-20150	BLUE - EXTRA HEAVY	50YD
REMTH-20160	BLACK - SPECIAL HEAVY	50YD
REMTH-20170	SILVER - SUPER HEAVY	50YD
REMTH-20180	GOLD - MAXIMUM	50YD



## THERA-BAND® TUBING

Thera-Band® Exercise Tubing includes natural rubber tubing available in 7 color-coded levels of resistance.

- Provides both positive and negative force on the muscles, improving strength, range of motion and cooperation of muscle groups.
- Different resistance levels are determined by thickness of the tubing.
- Color-coded levels increasing in difficulty from tan to silver.

REMTH-21110	TAN - EXTRA THIN	50YD
REMTH-21120	YELLOW - THIN	50YD
REMTH-21130	RED - MEDIUM	50YD
REMTH-21140	GREEN - HEAVY	50YD
REMTH-21150	BLUE - EXTRA HEAVY	50YD
REMTH-21160	BLACK - SPECIAL HEAVY	50YD
REMTH-21170	SILVER - SUPER HEAVY	50YD



## THERA-BAND® ACCESSORIES FOR TUBES AND BANDS

Accessories Kit

Contains Door ancor and two soft handles and a loop assist

REMTH-22135 1/KIT



## THERA-BAND® TUBING WITH HANDLES

Thera-Band® Resistance Tubing with Handles provides the same resistance, system of progression, and preferred quality as authentic Thera-Band Resistance Bands and Tubing with the added convenience of attached handles for a ready-to-use clinical or home solution.

- Available in 6 Thera-Band colors and resistances
- Soft handle grips, or standard PVC handle grips.
- Each tube is a minimum of 48" in length
- Supplied in attractive retail packaging for patient re-sale or economical bulk packaging
- Complete with usage and safety instructions



REMTH-21701	YELLOW	48" WITH PVC HANDLE
REMTH-21702	RED	48" WITH PVC HANDLE
REMTH-21703	GREEN	48" WITH PVC HANDLE
REMTH-21704	BLUE	48" WITH PVC HANDLE
REMTH-21705	BLACK	48" WITH PVC HANDLE
REMTH-21706	SILVER	48" WITH PVC HANDLE

## THERA-BAND® EXERCISE BALLS

The Thera-Band® Standard Exercise Balls are inflatable balls used to strengthen muscles to improve posture and help prevent back pain. Each exercise ball comes with two plugs and an inflation adapter, along with the Thera-Band® Exercise Ball Guide Poster detailing 24 exercises.


Five sizes are provided to allow for user height variations.

Ball Size: 45cm Height range: less than 5ft

Ball Size: 55cm Height range: 5ft - 5ft 7"

Ball Size: 65cm Height range: 5ft 8" - 6ft

Ball Size: 75cm Height range: 6ft 1" - 6ft 3"



REMTH-23010	YELLOW	45CM
REMTH-23020	RED	55CM
REMTH-23030	GREEN	65CM
REMTH-23040	BLUE	75CM





## THERA-BAND® SHOULDER PULLEY

The Thera-Band® Shoulder Pulley is specifically designed to assist patients undergoing shoulder rehab in regaining and maintaining range of motion. Unlike conventional pulleys, the Thera-Band Shoulder Pulley features a unique rope that's divided into sections by distinct black marks.

By referencing the position of these marks during exercise, patients get immediate feedback on progress. As they move from one black mark to the next during the course of rehab, they can literally see how much their range of motion is improving!

The black marks can also be used by clinicians to clearly communicate goals and limits, so that patients know how far to go and where to stop



REMTH-22160

## THERA-BAND® HAND EXERCISERS

Thera-Band® Hand Exercisers are specially formulated exercisers provided in a progressive sequence of five translucent colors and in two sizes. They're used for hand, finger and forearm strengthening, as well as hot and cold therapy. The softer yellow and red models fill the gap between putty and spring devices and they're not as messy as putty.

Progressive resistance in different levels of compression allow each exercise to be geared to the patient. Use the standard size hand exerciser for smaller hands and the extra large size for larger hands or for people who have limited hand mobility. Place a Hand Exerciser into the microwave or refrigerator for hot or cold use.



REMTH-26020	YELLOW
REMTH-26030	RED
REMTH-26040	GREEN
REMTH-26050	BLUE
REMTH-26060	BLACK



## THERA-BAND® FLEXBAR

Thera-Band® FlexBar® is used for wrist, forearm, and hand rehabilitation, upper extremity stabilization, neuromuscular re-education, massage therapy, martial arts, and sports grip strengthening.

- A flexible, durable resistance device with a ridged surface for enhanced grip during use
- It is used to improve grip strength and upper extremity stabilization by bending, twisting, or oscillation movement
- Research-proven effective for Tennis Elbow, offering a cost-effective treatment
- Thera-Band FlexBar is 12" long and is made from dry natural rubber in 4 progressive resistance levels to match user capability



REMTB-26107	YELLOW
REMTB-26100	RED
REMTB-26101	GREEN
REMTB-26102	BLUE
REMTB-26100	KIT - 3 OF EACH RED, GREEN & BLUE



## PROACTIVE PEDAL EXERCISER

For an easy, impact-free cardiovascular workout in the comfort of your home!

- Very compact retail box: 460 mm x 80 mm x 340 mm (18.1" x 3.1" x 13.4")
- Improves circulation, coordination and muscle strength
- Ideal for beginners, convalescents and seniors
- Can be placed on the floor for leg exercise or on a table for upper body exercise
- Adjustable resistance allows you to increase effort as you progress
- Compact and portable, it can be stored under a bed or chair



AMG740742

## MOIST HEAT PACKS


Moist heat packs are a simple and effective method of applying moist heat to a localized area.

- Provides 30 minutes of heat at a minimum.
- Packs can be applied with a terry or a foam filled cover and comes in various sizes and shapes to accomodate all applications.

STANDARD 10"X12"  
REMRM-09-111310

CERVICAL/NECK 9"X24"  
REMRM-09-111311

OVERSIZE 15"X24"  
REMRM-09-111312



## MOIST HEAT PACK COVERS

Terry Cloth hot pack cover allows therapeutic moist heat to penetrate without burning the skin. Terry covers are filled with foam or terry. Velcro fasteners keep moist heat packs in place. Covers are washable and can be used again.

### Terry Fill

REMRM-09-111366  
REMRM-09-111367  
REMRM-09-111368

STANDARD 20"X24"  
CERVICAL/NECK 24"X30"  
OVERSIZE 15"X24"



## THERATHERM® DIGITAL ELECTRIC MOIST HEATING PADS



Draws moisture from the air and retains it in a 100% cotton cover. Features a digital hand control allowing the operator to program the following parameters: treatment temperature (88°–166°), treatment time (1–60 min.) and on/off. Lockout mode protects against inconsistent changing of control settings throughout treatment.

Small	7" x 15"
Medium	14" x 14"
Large	14" x 27"
Cervical	20" x 23"

## ADHESIVE FELT



- 1/8 in Thickness
- White
- 5 1/2 in x 2 1/2 Yard Roll

SCHOLL58924

1RL/CS

- 1/4 in Thickness
- White
- 5 1/2 in x 2 1/2 Yard Roll

SCHOLL68926

1RL/CS

## TEN20 PASTE CONDUCTIVE PASTE



Ten20 Conductive Paste will hold the electrodes in place so there is no need to disrupt a test to reattach the electrodes. This reliable paste will conduct the signal that the technologist/technician is relying on for clear results.

### **BENEFITS:**

- Adhesive
- Washable and non-drying
- Easy clean up
- Convenient sizes

WEAVER10-20-4 4 OZ JAR

WEAVER10-20-4T 4 OZ TUBE

WEAVER10-20-8 8 OZ JAR

## NUPREP SKIN PREP GEL



To achieve quality tracings, proper skin preparation at the electrode site is crucial. Nuprep's mild abrasiveness strips away the top layer of skin and moistens the underlying skin layer. This action lowers the skin impedance with minimal skin irritation and patient discomfort.

### **MOST EFFECTIVE USES:**

- Audiology
- Cardiac catheter monitoring
- Cardiac rehabilitation monitoring
- ECG telemetry
- EEG exams
- Evoked potential procedures
- Holter monitoring
- Sleep tests
- Stress tests

Nuprep Gel is valuable whenever a reduction of skin impedance would enhance performance of a monitoring electrode. Nuprep virtually eliminates problems such as diaphoresis and unwanted artifacts.

WEAVER10-30 4 OZ TUBE

WEAVER10-61 25G TUBE



## ULTRASOUND GEL



PARKER01-08 250ML BLUE  
PARKER01-34 1 LT BLUE  
PARKER01-50 5 LT BLUE



PARKER03-08 250ML CLEAR  
PARKER03-34 1 LT CLEAR  
PACH200-C502 5 LT CLEAR

## BOTTLER WARMERS



PARKER82-01 SINGLE BOTTLE  
PARKER82-03 TRI-BOTTLE

## PROACTIVE TENS ELECTRO-STIMULATOR



- A CSA certified compact easy-to-use dual channel TENS nerve stimulator for active pain response that offers variable pulse width, frequency and amplitude controls.
- Choice of 3 versatile stimulation modes: burst, normal and modulated waveform settings to respond to most treatment protocols.
- Choice of 3 different timer settings: 30 minute, 60 minute and continuous for personal customization.
- 2 lead wires with built-in guard to protect against wire splitting
- 4 easy to apply 2" round re-useable self-adhesive comfort gel electrodes
- Custom designed sleek on-the-go pouch pack for compact storage and use on-the-go!

AMG715500

## TENS LEAD WIRES 2/PK

AMG715530



## TENS 2-IN-1 PHYSIOTHERAPY DEVICE ALEVIA™



AMG715425

The TENS 2-in-1 Physiotherapy Alevia™ by ProActive™ is a 2-in-1 (TENS and EMS) device that gives you the ability to fully customize treatment options or choose one of the 17 pre-set programs. It has a large, easy-to-read LCD and includes all accessories.

The Alevia™ comes complete with 4 self-adhesive, pre-gelled and reusable electrodes, one 9 V battery, lead wires, instructions and a carrying case.

ProActive™ offers a complete line of physiotherapy devices, clinically proven to relieve pain without drugs, including a full line-up of electrodes.

The Alevia™ Dual channel TENS Electro-Stimulator relieves pain in 2 ways:

- TENS (Transcutaneous Nerve Stimulation) stimulates or numbs nerves
- EMS (Electrical Muscle Stimulation) prevents atrophy

## TENS 3-IN-1 PHYSIOTHERAPY DEVICE THERA3™

The TENS 3-in-1 Physiotherapy Thera3™ by ProActive™ is a 3-in-1 (TENS, EMS & massage) device with 78 pre-set programs and the ability to fully customize treatment options. It has a large, easy-to-read LCD and includes all accessories.

The Thera3™ come complete with 4 self-adhesive, pre-gelled and reusable electrodes, 4 AAA batteries, lead wires, instructions and a carrying pouch.

ProActive™ offers a complete line of physiotherapy devices, clinically proven to relieve pain without drugs, including a full line-up of electrodes.

The Thera3™ Dual channel TENS Electro-Stimulator relieves pain in 3 way:

- TENS (Transcutaneous Nerve Stimulation) stimulates or numbs nerves
- EMS (Electrical Muscle Stimulation) prevents atrophy
- Massage soothes sore areas



AMG715430

## TENS ELECTRO STIMULATOR DEVICE PULSE™



AMG715420

The Pulse™ TENS by ProActive™ gives you a choice of 3 stimulation modes, a choice of 3 different timer settings and comes complete with all accessories so you can use it right away.

The Pulse™ comes complete with a durable plastic carry case, 4 easy to apply reusable electrodes, 1 9V battery, 2 lead wires and instructions. ProActive™ offers a complete line of physiotherapy devices, clinically proven to relieve pain without drugs, including a full line-up of electrodes.

The Pulse™ Dual channel TENS Electro-Stimulator relieves pain by giving you a choice of 3 stimulation modes:

Burst, normal and modulated waveform settings to respond to most treatment protocols.

This TENS also gives you the choice of 3 different timer settings: 15 minutes, 30 minutes and continuous.

## PROACTIVE™ SELF ADHESIVE ELECTRODES

Self-adhesive, pre-gelled & reusable Electrodes for use with most TENS and muscle electro stimulators

- Built-in conductive gel surface remains effective for up to 30 uses
- Flexible cloth backing allows the electrode to easily and durably conform to all body surfaces
- Pre-wired and pre-gelled for easy and convenient application
- Long lasting, hypo-allergenic gel keeps skin soft
- Rounded edges for prevention of electrode curl-up
- Anchored pigtail prevents wire pull-out and increases usage life
- Available in a variety of shapes and sizes for application to multiple treatment zones
- Latex-Free



AMG715710



AMG715715



AMG715720

## PROACTIVE™ SELF ADHESIVE ELECTRODES

For use with most TENS and muscle electro stimulators

- Flexible cloth backing allows the electrode to easily and durably conform to all body surfaces
- Pre-wired and pre-gelled for easy and convenient application
- Long lasting, hypo-allergenic gel keeps skin soft
- Rounded edges for prevention of electrode curl-up
- Anchored pigtail prevents wire pull-out and increases usage life
- Available in a variety of shapes and sizes for application to multiple treatment zones



AMG715725

*Note: Not sold in retail packaging.*

## CARBON ELECTRODES



Available with or without adhesive.

REM-NS5002P	SQUARE (1.75" X 1.75")	NON-ADHESIVE
REM-NS5008P	RECTANGLE (1.75" X 3.75")	NON-ADHESIVE
REM-S5002	SQUARE (1.75" X 1.75")	ADHESIVE
REM-S5008	RECTANGLE (1.75" X 3.75")	ADHESIVE

## KENDALL™ 200 SERIES FOAM ELECTRODES



KEND31478376

Kendall™ 200 series electrodes are a cost effective general adult monitoring electrode. The foam substrate repels fluids and conforms to body contours. The Kendall™ 200 series electrodes have a smaller footprint for convenient lead placement and a teardrop shape for easy removal from the liner and patient. Many packaging configurations are offered to satisfy all customer needs. The patented conductive adhesive hydrogel is designed to stay fresh up to 45 days out of the package.

## 3M™ RED DOT™ MONITORING ELECTRODES

A round-style monitoring electrode with 3M™ Micropore™ Surgical Tape backing.

Red Dot electrodes (Catalog Nos. 2239, 2249) are round-style electrodes designed for general or long-term clinical situations. Each electrode has a unique 3M™ solid gel conductive column and a 3M™ SureSeal moisture vapor barrier cap to ensure product freshness. Each electrode also has a 3M™ Micropore™ Surgical Tape backing that is comfortable, breathable, flexible and gentle to skin. This product is not radiolucent.

Red Dot 2249 electrodes are packaged 3 or 50 each per bag with an abrader pad on each electrode.



3M2249



# AVAILABLE AT BOWERS!



For people who suffer from chronic mucus congestion caused by COPD, asthma, emphysema, bronchitis, pneumonia, and other respiratory ailments, finding a way to dislodge the deep, thick mucus in the chest can seem impossible. The Lung Flute by Medical Acoustics is an effective, non-invasive device that is easy to use and drug-free – simply blow into it with no more force than is needed to blow out a single birthday candle.

The Lung Flute for Home Care offers these advantages:

- Easy to use
- No adjustments
- Used twice daily, improves quality of life
- Increased mucociliary clearance
- Low-pressure operation

MDCAC1003-01

MDCAC1002-01

LUNG FLUTE

REED REPLACEMENTS



## HOW TO ORDER

Visit one of our two Bowers Home Health stores:

### British Columbia

Unit 7 - 3691 Viking Way, Richmond

Tel: (604) 278-7566 or toll free 1-800-663-0047

Monday - Thursday 8 am to 5:30 pm - Friday 8 am to 4:30 pm

### Alberta

11783 - 186th Street, Edmonton

Tel: (780) 454-1666 or toll free 1-800-561-3687

Monday - Friday 8 am to 4:30 pm

Also, feel free to call or email us at [homehealth@bowersmedical.com](mailto:homehealth@bowersmedical.com)

**\*\*Please note Bowers Medical only ships within Canada. Prescriptions are not required in Canada. We accept all major credit cards. In stores we accept cash, credit & debit. Shipping charges & applicable taxes will apply on all deliveries.**

## INSTRUCTIONS FOR USE

The LUNG FLUTE® (Therapeutic) is indicated for Positive Expiratory Pressure (PEP) therapy. PEP therapy, when combined with "huff" coughing, will help remove secretions from airways, reduce the amount of air that may be trapped in airways, keep airways open and clear secretions, and improve delivery of bronchodilator medication. As with all therapies, you will need to begin slowly and build up the length and repetitions of your therapy sessions. Please read these Instructions for Use carefully before starting your new therapy. Consult with your healthcare professional to develop a personalized therapy regimen that best fits your needs.

### WARNINGS

- Do not attempt to inhale through the LUNG FLUTE®.
- If you experience shortness of breath and/or dizziness, discontinue use of the LUNG FLUTE® and consult your healthcare professional.
- A transient throat irritation, lasting less than 24 hours, was noted in some participants in clinical trials.
- Individuals who cannot follow a healthcare professional's verbal instructions should not use the LUNG FLUTE®. Specifically, young individuals who cannot follow verbal instructions, and adults who, due to physical or mental challenges, cannot follow verbal instructions from a healthcare professional.

### Getting the most benefit from the LUNG FLUTE®

Although your healthcare professional may have trained you in the use of the LUNG FLUTE®, these Instructions for Use will assist you in using this device. Using the LUNG FLUTE® on a regular basis in a manner that is comfortable for you will provide benefits that will help you breathe easier. If you have any questions about the LUNG FLUTE® or these instructions,

### The LUNG FLUTE® clears secretion build up

In many chronic cardiopulmonary conditions there is a build up of secretions in the airways. Bronchial hygiene therapy with the LUNG FLUTE® can help clear these secretions so they can be swallowed or coughed up.

The LUNG FLUTE® is a device used to help loosen, mobilize, and eliminate secretions from your airways. It consists of a mouthpiece and a reed inside a horn.



### How the LUNG FLUTE® works

When you blow gently into the LUNG FLUTE® in a series of repetitions, your breath moves the reed inside (see Figure 2 below). This causes acoustic vibrations, which thin and loosen secretions deep in your lungs and results in the secretions moving progressively up your airways until they collect at the back of your throat.



### Proper use of the LUNG FLUTE®

Proper use of the LUNG FLUTE® is important for successful therapy. Although the technique presented here works well for most users, "individualizing" the technique for your specific condition may be necessary in order to obtain the best results.

Please meet with your healthcare professional for complete personal training before beginning your LUNG FLUTE® therapy. Continue to review your technique on a regular basis, especially if you are not getting the expected results or if you notice a change in how your therapy is progressing.

### Preparing for your LUNG FLUTE® therapy

1. You may want to have a glass of water available to drink after your therapy session. Before you begin, think about what you will be doing, how you are going to breathe, and the goals of the therapy. Visualize your airways being vibrated by acoustic waves, and secretions being loosened and mobilized along with your exhaled air to progressively larger airways, where they can be coughed up.
2. Get into a relaxed position. Sit up straight so that your back is not touching the back of the chair. Tilt your head slightly downward so your throat and windpipe are wide open (see Figure 2, left). This allows the acoustic waves produced by your breath to flow from the LUNG FLUTE® into your lungs.
3. If you must perform therapy in bed, sit as upright as possible, in a position that will not restrict your smooth breathing effort.

#### A NOTE ABOUT ACOUSTIC WAVE TECHNOLOGY

The acoustic waves produced by the LUNG FLUTE® cannot be heard by human ears. The only sound you will hear is the fluttering of the reed when you exhale through the LUNG FLUTE® mouthpiece. You will not feel any vibration of the LUNG FLUTE® during operation.

### Stage One: Secretion Loosening and Mobilization

Hold the LUNG FLUTE® pointing down, as shown in Figure 2. Inhale a little deeper than normal, place your lips completely around the mouthpiece, and gently blow out through the LUNG FLUTE® as if you were trying to blow out a candle. As you blow you will hear the reed inside the LUNG FLUTE® make a fluttering noise as it moves. Concentrate on making this sound. (If the reed does not move right away, practice with the LUNG FLUTE® for several sessions. If you are still unable to move the reed, consult your healthcare professional.) Remove the mouthpiece from your mouth and quickly inhale again. Put the mouthpiece back in your mouth and blow gently through the LUNG FLUTE®.

Remove the mouthpiece and wait five seconds, taking several normal breaths.

As you adapt to using the device, your goal should be to complete up to 20 sets of 2 blows each.

Depending on your condition, you may wish to perform only a few repetitions of 2 blows, one or two times a day, until you are able to do more.

After several sessions, you should be able to determine for yourself the appropriate number of 2-breath repetitions needed to adequately clear secretions.

Remember, you only need to blow through the mouthpiece with as much force as you would to blow out a candle. Do not force a cough or use your diaphragm or stomach muscles to try to force out more air.

### Stage Two: Secretion Elimination

Once you have worked up to the required number of 2-blow repetitions, the LUNG FLUTE® has thinned and loosened your secretions. Wait five minutes after your session for secretions to collect at the back of your throat. Several minutes of coughing should bring up most of the secretions.

You may notice thinned secretions collecting at the back of your throat for several hours after your session. This is normal. A drink of water can prevent minor throat irritation.



LUNG FLUTE is a registered trademark of Medical Acoustics, LLC  
©2010 Medical Acoustics, LLC U805 Rev. A (07 /2011)

### Recommended length and frequency of sessions

Your healthcare professional will advise you on how often you should use the LUNG FLUTE® and how long each session should last. Generally, morning and late afternoon/evening sessions are recommended. The average session lasts 5 to 10 minutes, depending on the severity of the pulmonary condition. To avoid becoming overly tired from therapy, it is better to add a session rather than extending a session to the point of discomfort.

### Important Tips

- When inhaling, remove the LUNG FLUTE® from your mouth.
- Pointing the LUNG FLUTE® slightly toward the floor may make it work more efficiently.
- To avoid dizziness or shortness of breath while using the LUNG FLUTE®, take more time between each set of two blows.
- Beginning LUNG FLUTE® users or those with thick mucus may experience throat irritation caused by large quantities of secretions collecting at the back of the throat. Drinking water after a session can help reduce this irritation.
- Performing 3 “huff” coughs after each 5 sets of 2 blows may further assist you in clearing your secretions.

### Condensation build up inside the LUNG FLUTE® is normal

If condensation builds up inside the horn place it in a warm spot, such as a sunny window. The condensation should be gone before you start your next session. If condensation persists, remove the mouthpiece and let the horn air dry.

### Replacing the reed

The reed inside your LUNG FLUTE® should be replaced about every two weeks depending on the frequency of use.

**Replacements reeds are available by calling any of our Bowers Medical locations:**

**BC - (800) 663-0047**

**AB - (800) 561-3687**

**ON - (888) 216-4833**

**www.bowersmedical.com**

**email: homehealth@bowersmedical.com**

### Cleaning the LUNG FLUTE®

When you replace the reed, wash the inside of the horn and mouthpiece with mild dish soap and hot water. Allow the pieces to air dry or dry them manually before reassembling. If a film build up is noted inside the LUNG FLUTE® or if it does not dry between sessions, wash it more frequently.

### Product Performance

Frequency 16–25Hz

Minimal flow rate 128.4 L/min

Minimal pressure 1.0 cm H<sub>2</sub>O

Sound output 68 dBA1

Pressure resistance 1.0 cm H<sub>2</sub>O

1 Measured with a standard General Radio 1933

Precision Sound Level Meter. OSHA limit over a 15 min. interval is 115dBA CFR 29.1910.95 (b)(2)



**LUNG FLUTE is a registered trademark of Medical Acoustics, LLC**

**©2010 Medical Acoustics, LLC U805 Rev. A (07 /2011)**

Instructions for use obtained from Medical Acoustics/Lung Flute website

•Therapeutic Lung Flute for Home Care IFU for Patients [PDF]

<http://www.lungflute.com/PDF/LungFlute-Therapeutic-HomeCare-IFU.pdf>

## ELONGATED MEDIUM CONCENTRATION OXYGEN MASK



Transparent mask with soft anatomical shape ideal for long-term users and edentulous patients. Over the ear style fits comfortably to hold the mask snugly in place. Exceptionally comfortable for even hard-to-fit faces. Three-channel crush resistant tubing standard.

Elongated medium concentration mask with 7' (2.13m) three channel safety tube, elastic headstrap

### ADULT

SALT8110 50/CS

### PEDIATRIC

SALT1122 50/CS

## THREE-CHANNEL OXYGEN SUPPLY TUBING



Special design has three inside channels that help prevent the tubing from becoming occluded. Ideal for home patients where long lengths can cause a potential kinking problem.

Available in lengths from 4' to 100'. Special or custom lengths available...a Salter service oriented benefit.

Note: Oxygen supply tubing also available in smooth bore configuration.

### Innovative Tubing End Fitting\*

- Vinyl material with sufficient clarity to view patency
- Designed with "Ribbed Body" to permit easy gripping when making connections
- Rounded and tapered edges to prevent snagging or catching on various surfaces



Clear 14' (4.27m) oxygen supply tube, three channel safety - 3/16" I.D.\*

SALT2014



## INSPIRATION 626 COMPRESSOR NEBULIZER SYSTEM

Inspiration 626 was designed with compliance in mind. Multi-media patient education support materials and built-in conveniences make this efficient nebulizer system the ideal choice for encouraging adherence to aerosol medication treatment plans. Inspiration 626 Compressor Nebulizer is ideal for patients of all ages..



**Replacement Filters**  
UNI02213148 4/PK

UNI00334573

## COMPRESSOR NEBULIZER

Nebulizes liquid medicine using compressed air technology, for treatment of the upper and lower respiratory tract.  
Ideal for treating: Colds • Asthma • Respiratory diseases

Kit includes:

- Adult and child mask
- Mouthpiece
- Nebulizer unit
- Tubing
- Nebulizer chamber
- Replacement air filters



AMG705470

## NONIN GO2 PULSE OXIMETER

The same technology trusted by physicians, clinicians and therapists is now available in a small, portable and easy-to-use personal fingertip oximeter. Simply insert your finger into Nonin GO2 brand fingertip pulse oximeter for almost instant, highly accurate oxygen saturation and pulse rate information.

UNI02313393



## DIAGNOSTIX™ 2100 FINGERTIP PULSE OXIMETER

A portable, feature rich pulse oximeter that provides spot readings of SpO2 and heart rate for the busy healthcare professional.

The 2100 features:

- Reads SpO2 and Pulse Rate
- Accuracy: SpO2  $\pm 2\%$  / HR  $\pm 5$  BPM
- 6 changeable display modes
- 10 Display brightness settings
- 25% smaller than competing models
- Fits fingers with a thickness of .3" to .8" (6.6mm to 19.7mm)
- Dual color OLED display
- Hinged battery cover is more secure
- Complete with Carry Lanyard and Black Rubberized Bumper. (Optional orange bumper sold separately)
- 2 "AAA" Batteries
- Instruction Booklet
- 2 year limited warranty



ADC2100

## ADVANTAGE FINGERTIP PULSE OXIMETER

Reads SPO2 and Pulse Rate

- Accuracy: SPO2  $\pm 2\%$  / HR  $\pm 5\%$  BPM
- 2 changeable display modes
- Bright red LED display
- Fits fingers with a thickness of .3" to .8" (6.6mm to 19.7mm)
- Carry Lanyard
- 2 "AAA" Batteries
- Instruction Booklet
- 1 year limited warranty



ADC2200

## CLINICIAN STETHOSCOPE



\*\*Available in a variety of colours

ADC603

Combination diaphragm/bell adult chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability

- Extra large bell (1 3/8" diameter) for unsurpassed low frequency response
- Ultra-sensitive diaphragm (1 3/4" diameter) for greater amplification and crisper high frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Clinician headset boasts reinforcing yoke molded into the flexible 22" PVC tubing and stainless steel binaurals fixed at 15° angle for greater comfort
- Adsoft™ threaded PVC eartips for the ultimate in wearing comfort and acoustic seal
- Includes accessory kit with spare diaphragm, extra Adsoft™ and standard eartips
- Scope ID Tag included
- Weighs 5.77 oz. (165 gr)
- Overall length 31"
- Assembled, inspected and packaged in the USA
- Lifetime Warranty covers all parts

## ADSCOPE™ SPRAGUE STETHOSCOPE

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings - an adult and pediatric diaphragm, and three bells (1 1/4", 1", 3/4") - for complete diagnostic capability on most patients

- Latex-free, double tube configuration
- Proprietary valve mechanism eliminates acoustic leakage
- Adjustable chrome-plated brass binaurals
- Convenient 22" tubing length in choice of 29 fashion tubing colors
- Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips, and one pair of silicon eartips
- Weighs 7 oz.
- Overall length 30"
- 5 Year Warranty
- Inspected and packaged in the U.S.A.



\*\*Available in a variety of colours

ADC641

## ADIMALS™ PLATINUM PEDIATRIC STETHOSCOPE

Distinctive molded resin animal face chestpiece optimally sized (1½") for pediatric use

- Multi-frequency design switches between bell-like low frequency and diaphragm-like high frequency performance simply by altering pressure
- Non-chill bell and diaphragm retaining rim for patient comfort
- Clinical headset boasts reinforcing yoke molded into the flexible 22" PVC tubing
- Lightweight aluminum binaurals in matching satin finish
- Adsoft™ threaded PVC eartips for the ultimate in wearing comfort and acoustic seal
- 7 interchangeable animal face snap-ons included:
  - Panda
  - Frog
  - Koala
  - Deer
  - Monkey
  - Tiger
  - Bear

- Weighs: 3 oz
- Overall length 30"
- 5 year Warranty



ADC618

Available colours:



## LITTMAN LIGHTWEIGHT



3M2450	28 in. (71.1 cm)	Black
3M2451	28 in. (71.1 cm)	Burgundy
3M2452	28 in. (71.1 cm)	Caribbean Blue
3M2453	28 in. (71.1 cm)	Lilac
3M2454	28 in. (71.1 cm)	Ceil Blue
3M2455	28 in. (71.1 cm)	Seafoam Green

### 3M LITTMAN STETHOSCOPE COLOUR CHART



## LITTMAN CLASSIC II



3M2201	28 in. (71.1 cm)	Black
3M2203	28 in. (71.1 cm)	Gray
3M2205	28 in. (71.1 cm)	Navy Blue
3M2206	28 in. (71.1 cm)	Caribbean Blue
3M2208	28 in. (71.1 cm)	Hunter Green
3M2209	28 in. (71.1 cm)	Purple
3M2210	28 in. (71.1 cm)	Raspberry
3M2211	28 in. (71.1 cm)	Burgundy
3M2812	28 in. (71.1 cm)	Orange
3M2813	28 in. (71.1 cm)	Ceil Blue
3M2814	28 in. (71.1 cm)	Seafoam Green

## DUAL-HEAD STETHOSCOPE



ALM10-1310	BLACK
ALM10-1320	GRAY
ALM10-1330	NAVY BLUE
ALM10-1350	HUNTER GREEN
ALM10-1360	RED
ALM10-1370	PINK
ALM10-1380	LIGHT BLUE
ALM10-1400	TEAL
ALM10-1410	PURPLE
ALM10-1420	BURGUNDY

## TEACHING STETHOSCOPE



ALM10-5000



## WRIST BLOOD PRESSURE MONITOR



AMG106962EA

Track your blood pressure anywhere!

With fuzzy logic technology that "personalizes" your inflation level for a more comfortable reading!

Measurements obtained with this unit are equivalent to those obtained by a trained person using the cuff/stethoscope auscultation method, as per the American National Standard, Electronic or automated sphygmomanometers.

One-Touch Operation

Operates automatically at the touch of a button, for ease of use.

Stores up to 85 readings, because trending is important.

- Date and time indications
- Low Battery Indicator
- Auto Power-Off 150 seconds after last button operation

## PRISMA BLOOD PRESSURE MONITOR



AMG106950

The Prisma Blood Pressure Monitor has a unique design.

This fully automatic blood pressure monitor offers you easy to use features for optimal accuracy and comfort:

• Large backlit screen

- Auto memory for up to 60 memories
- Flex-Fit cuff (fits 8.7 in to 16.6 in arm)
- Hypertension indicator
- Date and time indicator
- Powered with 4 "AAA" batteries (not included)
- One-touch operation
- Easy-to-read LCD display with simultaneous readout of systolic, diastolic and pulse measurements
- 2 year limited warranty

## ESSENTIA BLOOD PRESSURE MONITOR



AMG106930

- Large easy to read screen
- Multi-user (4 x 30 memories)
- Heart rate indicator
- Date and time indicator

Meets all required Canadian & US FDA Standards including ANSI/AAMI SP-10 & ESH International Protocol

Includes:

Blood pressure monitor and medium size arm cuff (fits 8.7" to 14.2" circumference)  
Requires 4 AA batteries (not included)

## ADVANTAGE™ SEMI-AUTO DIGITAL BP MONITOR

Measures systolic, diastolic pressure, and pulse

- Manual inflation
- Self adjusting latex-free D-Bar cuff
- Automatic deflation valve
- Auto off

ADC6012



## ADVANTAGE™ AUTOMATIC DIGITAL BP MONITOR



ADC6021

Measures systolic, diastolic pressure, and pulse

- Smart Logic™ Technology determines ideal inflation level
- Irregular Heartbeat detector
- Automatic inflation and deflation control
- Self adjusting latex-free D-Bar cuff
- Recalls last 99 readings with date and time stamp
- Auto-off
- Optional AC adapter (sold separately)

## DISPOSABLE PENLIGHTS



Sealed concentric beam, disposable diagnostic penlight

- Activated by depressing pocket clip
- Available with or without pupil gauge (white barrel)
- Measures 5" long, 1/2" diameter
- Sold individually or in packs of 6/pk
- Available in 8 colors

PROMAR5816 6/PK  
WHITE

ADC356BK	BLACK
ADC356OR	ORANGE
ADC356R	RED
ADC356RB	ROYAL BLUE
ADC356TL	TEAL
ADC356V	PURPLE
ADC356W	WHITE
ADC356WP	WHITE W/PUPIL GAUGE



## SELF TAKING HOME BLOOD PRESSURE KIT



This easy-to-use, self-taking blood pressure kit is ideal for home use. It features a certified no-stop pin, a 300 mm Hg aneroid gauge, a self-adjust D-ring cuff and a built-in stethoscope. Comes complete with vinyl zipper case, instruction booklet, and recording chart.

ALM14-1100

## SPHYGMOMANOMETER

This B.P. kit features an aneroid with no-stop pin 300 mm Hg gauge, high quality inflation system and nylon cuff with carrying case.



ALM14-2200

## PROSPHYG™ HOME BP UNIT



Black enamel 300mmHg manometer

- Chrome-plated deflation valve
- Nylon D-Ring self-adjusting cuff for one-handed cuff application
- Diaphragm only blood pressure scope with black PVC tubing (#660BK)
- Zippered carrying case and operating instructions
- Assembled, inspected, and packaged in the U.S.A
- Meets or exceeds ANSI/AAMI/ISO 81060-1

ADC780-BP

## PROSPHYG POCKET ANEROID SPHYG



Contemporary black enamel 300mmHg manometer provides day-in, day-out dependability exceeds industry standards

- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing
- Chrome-plated brass air release valve
- Leatherette carrying case with durable nylon zipper
- Assembled, inspected, and packaged in the U.S.A. from foreign components
- Meets or exceeds ANSI/AAMI/ISO 81060-1

ADC760

## PALM ANEROID SPHYG

A versatile palm style aneroid sphygmomanometer boasting precision crafted performance and rugged durability.

The 703 series features:

Revolutionary palm style manometer (#804) with integral latex-free bulb and trigger style air release valve for one-hand convenience  
Precision crafted 300mmHg Japanese engineered movement tested to 50,000 cycles

Movement housed in a zinc cage that suspends within the polymer housing for enhanced shock resistance

Ambidextrous infinite control trigger release air valve

Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing

Quick release luer slip connector

Complete with leatherette carrying case with self repairing nylon zipper

Assembled, inspected, and packaged in the U.S.A. from foreign components

Meets or exceeds ANSI/AAMI/ISO 81060-1



ADC703

## MULTIKUF™ GENERAL PRACTICE MULTICUFF KIT

A multicuff kit designed for the busy office, clinic or ER.

The 705GPK series features:

ADC 804 palm gauge PLUS 3 cuffs with single tube bladders

Rugged gauge features precision crafted Japanese engineered movement tested to 50,000 cycles

Ambidextrous trigger release air valve

Small Adult, Adult, and Large Adult calibrated nylon Adcuffs™

Quick release luer slip connectors for rapid selection of appropriate cuff

Four section smoke gray plexiglass tray features a small 14" x 3 1/2" footprint

Tray may be wall mounted with included hardware

Assembled, inspected, and packaged in the U.S.A. from foreign components

Meets or exceeds ANSI/AAMI/ISO 81060-1



ADC705GPK-MCC



## WALL ANEROID SPHYG

Large, easy-to-read 6" luminescent dial for high visibility in virtually any light condition

- 300mmHg no-pin stop movement withstands the heaviest use
- Adcuff™ inflation system includes the Adcuff™ nylon cuff, Adflow™ bulb and valve, and eight foot length (when extended) of coiled latex-free tubing
- Contemporary Euro-style ivory casing compliments any decor
- Wall bracket that swivels through 110° rotation for easy viewing with integral cuff storage. Optional cuff basket available
- Complete with all necessary hardware and assembly instructions
- Inspected and packaged in the U.S.A. from foreign components
- Meets or exceeds ANSI/AAMI/ISO 81060-1



ADC750W

## BRONZE SERIES DS44 INTEGRATED ANEROIDS

Gear-free DuraShock™ design is much more durable than traditional aneroids with gears. Part of our Bronze Series, the DS44 Integrated features a FlexiPort cuff-mounted design that is compact and lightweight for ease of use.

Learn more about DuraShock technology.

- Integrated design
- Premium inflation system
- 5-year calibration warranty
- Latex-free



WELCHDS44-\*\*

## 767 WALL AND MOBILE ANEROIDS



Improved ergonomics and a new modern design make the 767 Aneroid an easy-to-use and attractive blood pressure aneroid. Based on Welch Allyn's near century of experience building and designing innovative blood pressure technology, it's the smart choice for accurate and reliable readings in your facility.

- Laser-engraved dial face provides unmatched accuracy
- Recessed dial face and high contrast red pointer increase visibility and reduce parallax errors
- Lifetime calibration warranty
- New 5-leg stand provides improved stability and maneuverability
- Simple color changes give the aneroid a new look
- Latex and Mercury Free

767 Wall Aneroid with Reusable One-Piece Adult Cuff

WELCH7670-01



767 Mobile Aneroid on 5-Leg Stand with Reusable One-Piece Adult Cuff

WELCH7670-03

## BLOOD PRESSURE CUFFS

BP Cuffs available in:  
1 tube or 2-tube  
Adult, Child & Pediatric sizes

WELCH5082-\*\*



## SPOT VITAL SIGNS® DEVICES



Vital signs in seconds with hospital-grade technology that's easy to use, accurate and affordable.

- Measures noninvasive blood pressure (NIBP), pulse rate, temperature and SpO2 in about 30 seconds
- IR Connectivity to Electronic Medical Records
- Lightweight and portable with a variety of mounting options
- All configurations include NIBP, pulse rate and MAP
- Optional—SureTemp® thermometry, Masimo® or Nellcor pulse oximetry

Spot Vital Signs with NIBP only (pulse rate and MAP)

\*\*Other configurations available.

WELCH4200B-E1

## SPOT VITAL SIGNS® LXI



WELCH450T0-E1

Take vital signs in about half the time. Take patient care to a whole new level.

- SureBP® technology measures blood pressure in about 15 seconds as the cuff is inflating—resulting in time savings for you and improved comfort for your patients, pulse rate and MAP included
- Choice of SureTemp® Plus or Braun ThermoScan® PRO 4000 ear thermometer
- Choice of Masimo® or Nellcor® pulse oximetry
- Ability to manually enter weight, height, respiration rate, and pain level; calculates Body Mass Index
- Connects to select Healthometer® or Tanita® digital scales
- 50-reading memory allows retrieval of previous data for review, to print or to download later
- Connectivity-ready

- Error messages give you descriptions and real solutions—not just error codes
  - Bigger, brighter LCD display provides greater visibility from farther distances
  - Now available with wireless connectivity option to transmit vitals directly to your hospital EHR! (wireless option available in U.S. and Canada only)
  - Weight Management USA and Canada: Spot Vital Signs LXI can be connected to a broad range of Health o meter Professional Scales. For more information, please go to [http://www.hom scales.com/innovation/WA\\_HOM\\_Connectivity\\_List.aspx](http://www.hom scales.com/innovation/WA_HOM_Connectivity_List.aspx) or call Pelstar/Health o meter Professional Scales at: 800-815-6615
  - Weight Management Outside USA and Canada: Spot Vital Signs LXI can be connected to a broad range of co-branded Welch Allyn/Health o meter Professional Scales.
- For more information go to: Healthometer Scales Page

## ADVIEW VITAL SIGNS

The complete Advview™ diagnostic station includes:

- An automated BP monitor that measures systolic and diastolic pressure, pulse rate and MAP in about 30 seconds using the oscillometric method.
- A temperature module using Kendall's Filac™ Fastemp™ technology that provides an oral or axillary temperature in about 10 seconds. Built-in storage well for the probe and box of 20 disposable probe covers.
- An SPO2 module that features industry leading Nellcor compatible technology for accurate readings in about 10 seconds.
- Optional EMR (electronic medical records) with BLUETOOTH® connectivity for wireless integration to your computer system.

Unlike EVERY other monitor on the market, the Advview™ is completely modular. You can purchase the Advview™ configured as you want from the factory OR upgrade the unit in the field at any time. No tools needed. With connect-and-go capability, adding optional temperature and pulse oximetry modules is effortless.

### Start With Blood Pressure

Designed with the idea that no two patients are alike, the ADC Advview™ Diagnostic Station continues the innovation first developed for the e-sphyg 2 by offering a manual measurement mode. By combining the best of our nearly half century years of blood pressure measurement experience with the capability to apply traditional methods of blood pressure measurement, the ADC Advview™ provides the most reliable and versatile means of capturing blood pressure data.

Easily capture automated blood pressure measurements at the touch of a button or choose to measure with your stethoscope by using the manual feature for your patients with arrhythmias or other complex needs. We recognize that all patients are different. Now, everything you need for accurate blood pressure is in one intuitive device. With 5 cuff sizes, ranging from child to thigh, the ADC Advview™ caters to all of your patients.\*

\* The ADC Advview™ is intended for patients ages 3 and up.

#### What's Included:


- Adult & Large Adult All Purpose Cuffs
- 8' Blood pressure hose
- Mains/Power Supply
- Mains/Power Cable
- Resource CD
- Quick Start Guide
- Reference Card

#### Modular Options:

- Automated Temperature Readings
- SpO2 Readings in an Instant

#### Mounting Options

- Wall mount
- Table Stand
- Mobile Stand



## COMPLETE DIAGNOSTIC INSTRUMENT SET



ADC5215

Replacement Parts Available

Full featured 2.5v set with 5 different instruments in a fitted case.

The complete set features:

- Bayonet locking mounts for each instrument head to allow secure connection to the battery handle
- "C" cell powered, satin finish battery handle with knurled finish
- Rheostatic on/off switch
- Complete in zippered case with fitted insert that is scuff and scratch resistant and cleans easily
- Weighs 23 1/2 oz
- Inspected and packaged in the U.S.A. from foreign components
- Batteries NOT included

Otoscope features:

- Chrome-plated Otoscope head with removable 4x magnifying lens, and insufflator fitting
- Auto-clavable, polypropylene ear specula in 2.5mm, 3.5mm, and 4.5mm sizes
- Adapter to allow use with Welch Allyn or compatible disposable specula
- Nasal speculum adapter with manual screw opening for instrumentation connects to the otoscope head
- Welch Allyn disposable polypropylene ear specula in 2.75mm and 4.25mm sizes are convenient, safe, and economical

Ophthalmoscope features:

- Ophthalmoscope head with lens wheel featuring 24 corrective lenses from -25 to +40 diopters

Bent Arm Illuminator features:

- Chrome-plated bent arm larynx and nasopharynx illuminator for lighting the oral cavity and pharyngeal area
- One 3.0mm and one 4.0mm mirror with integral tongue depressor holder that attaches to the illuminator
- One tongue depressor attachment for the illuminator
- One disposable tongue depressor adapter for the illuminator



## POCKET OTOSCOPE



Fiber optics for cool, obstruction-free light transmission

- 2.5v Xenon lamp for true tissue color
- Wide angle swivel 2.5x viewing lens permits instrumentation
- Insufflator port (fitting and insufflator sold separately) for pneumatic otoscopy
- "AA" battery handle (batteries included) with on/off power switch integrated into pocket clip
- Includes tube with 5 each Welch Allyn 2.75mm and 4.25mm specula
- Choice of fitted hard case or nylon pouch with belt loop
- Oscope in hard case measures 7 1/4" x 4 1/2" x 1 1/2", weighs .78 lbs.
- Oscope in soft case measures 1 1/4" x 7" x 3 3/4", weighs .42 lbs.
- 1 year instrument warranty, lifetime warranty on optics

ADC5111N

## POCKETSCOPE™ OTOSCOPE WITH "AA" HANDLE AND SOFT

Small, convenient, and easy to use with bright, white halogen light and a built-in throat illuminator.

- Halogen light for true tissue color and consistent, long-lasting illumination
- Simply detach the otoscope head to use the throat illuminator or all-purpose penlight
- Fiber optics provide cool light with no reflections or obstructions
- Wide-angle viewing lens allows for instrumentation under magnification
- Sealed system for pneumatic otoscopy

WELCH22821



## 3.5 V COAXIAL OPHTHALMOSCOPE



WELCH11720

The patented Welch Allyn Coaxial Vision System facilitates ophthalmoscopy by enabling easier entry into the eye, a larger field of view, and reduced glare compared to standard ophthalmoscopes.

- Halogen HPX™ lamp provides light for true tissue color and long-lasting performance
- Coaxial optics produce a shadow-free spot, easier entry into undilated pupils, and a larger field of view versus standard ophthalmoscopes
- Detect corneal abrasions with cobalt blue filter
- Polarizing filter virtually eliminates corneal reflection
- Sealed optics keep out dust and dirt
- 18 unique aperture/filter combinations for greater versatility
- Red-free filter may be used with any aperture
- 28 focusing lenses with a range of -25 to +40 diopters

## POCKET OPHTHALMOSCOPE



5 aperture selections (small spot, large spot, semicircle, fixed and red-free filter)

- 19 viewing lenses from -20 to +20 diopters for better resolution
- Illuminated magnified lens dial
- Rubber brow rest
- 2.5v Xenon lamp for true tissue color
- "AA" battery handle (batteries included) with on/off power switch integrated into pocket clip
- Choice of fitted hard case or nylon pouch with belt loop
- Ophthalmoscope in hard case measures 7 1/4" x 4 1/2" x 1 1/2", weighs .82 lbs.
- Ophthalmoscope in soft case measures 1 1/4" x 7" x 3 3/4", weighs .38 lbs.
- 1 year instrument warranty, lifetime warranty on optics

ADC5112N

## OTOSCOPE SPECULAS



4.25 mm KleenSpec® Disposable Ear Specula for Diagnostic Oscopes

WELCH52432      850/BG 10BG/CS

2.75 mm Universal KleenSpec® Pediatric Disposable Ear Specula for Diagnostic Oscopes

WELCH52434      850/BG 10BG/CS

## BRAUN THERMOSCAN® PRO 4000



This new model provides the speed and accuracy you need for all types of patients. The exclusive ExacTemp™ technology provides an active user feedback system to ensure proper positioning and improved accuracy.

- Improved accuracy and repeatability through a preheated probe tip and sensor
- Probe cover detection system ensures precise infrared heat transfer for accurate temperatures
- Automatic probe cover eject button for quick and easy disposal of used probe covers
- Easy-to-read, LCD displays temperatures in Fahrenheit or Celsius
- Three-year limited warranty

WELCH04000-200

## DISPOSABLE PROBE COVERS



Braun ThermoScan PRO 4000 Disposable Probe Covers

WELCH05075-005      200/BX 25/BX/CS

## SURETEMP® PLUS 690

Easy to use with intuitive operation, reduces the risk of cross-contamination with removable probe well and provides more reliable readings with improved accuracy. Takes approximately 4-to-6 second oral temperatures, 10-second axillary and rectal temperatures and 15-second adult axillary temperatures.

- Interchangeable, removable probe wells reduce the risk of cross-contamination
- Large LCD displays temperatures in Fahrenheit or Celsius
- Last temperature recall
- Convenient storage for 25 probe covers
- Takes approximately 6,000 readings on 3 AA batteries



SureTemp® Plus 690 Electronic Thermometer, 4 ft Cord and Oral Probe with Oral Probe Well

WELCH01690-200

SureTemp® Plus 690 Electronic Thermometer, One Room, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well

WELCH01690-300

## DISPOSABLE PROBE COVERS



WELCH05031-101

250/BX

## ADTEMP TYMPANIC THERMOMETER

1 second measurement using scan/peak technology

- Large 6mm illuminated display
- Multiple Use (wide range 89.6°F - 108°F +/- .4°F (32°C - 42.2°C +/- .2°C)
- °F/°C switchable display
- Fever alarm at 99.5°F (37.5°C)
- Last reading memory
- Audible beep when measurement is complete
- Auto off conserves battery life
- Compact design - 5 1/2" (13.97cm) long, weighs 1.6oz. (45.3g) including battery
- Powered by 1.55v (CR2032) Lithium battery (included)
- Includes 20 probe covers in convenient dispenser
- Desktop storage caddy
- Meets ASTM E1965-98 standards for accuracy



ADC425

### ADTEMP SHEATHS

ADC425-40      40/PK

## ADTEMP 6 SECOND CONDUCTIVE THERMOMETER

Non-invasive forehead temperature measurement using patented conductive technology delivers clinically accurate readings in about 6 seconds

- Reduces the risk of cross contamination and avoids the need for costly probe covers
- Ideal for the hospital, clinic, or physician office
- Auto off function conserves battery life
- Audible tone when measurement is complete
- Last reading memory
- Dual Scale (°F/°C)
- Water Resistant
- Large display
- Uses 2 "AAA" batteries (included)  
(Minimum 1000 full duty cycles over 24 month period)
- Temp range 95°F to 107.6°F +/- .2°F (35°C to 42°C +/- .1°C)
- Compact design measures just 4.7" x 1.17" x 1" (12.0 cm x 3.0 cm x 2.4 cm ) weighs 1.3 oz (37g) including batteries
- Meets the accuracy requirements for intermittent thermometers under ASTM E1112





## FEVER BUGZ



### Features

- Ideal for kids of all ages
- they love to wear them
- can be worn while bathing
- Hypo allergenic
- Fast, Gentle & Safe application
- Fuss-free for parents
- 48-hour continuous monitoring
- Visible at a glance
- Accurate within 0.5°C
- Clinically proven accuracy for readings you can trust
- Latex-free

AMG016700

## ACCUFLEX PRO, DIGITAL ORAL THERMOMETER

- Accurate within  $\pm 0.1^{\circ}\text{C}$  in the  $35^{\circ}\text{C}$  to  $39^{\circ}$  range
- Faster response time: 5.6 sec. for oral readings
- Waterproof, washable for easy cleaning
- Gentle-flex tip for patient comfort
- Auto shut-off feature
- Memory recall / display of last reading
- Includes easy to use, pop-out storage case
- Convenient color codes to differentiate uses (blue for oral and red for rectal)



AMG016638

## KLEENSPEC® VAGINAL SPECULA

Welch Allyn's single-use specula eliminate the time and expense of reprocessing and reduce the chance of cross-contamination. Each speculum is designed and manufactured to deliver durability and reliability.



WELCH59000  
SMALL 24/BX



WELCH59001  
MEDIUM 24/BX

- Newly designed light pipe provides a homogeneous spot and eliminates glare back
- New wider handle provides improved ergonomics, better balance, and is much easier to manipulate during exams
- Smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during the exam
- The new 590 Series KleenSpec® Vaginal Speculum is compatible with the new cordless illumination system (catalog #79910) as well as the current illumination system (catalog #78810)

## DISPOSABLE VAGINAL SPECULUM

Recommended for routine gynaecological exams.

- Made from hard, clear, non-toxic plastic resin of high strength.
- Allows examination of the vaginal mucosa at 360° angle without rotation.
- Provides wide range of adjustability, operation and manipulation.
- Lock-mechanism holds in place in various open positions.
- Adjusts for elevation and angle.
- Can be manipulated with one hand.
- Crystal clear transparency ensures maximum visibility, with no tissue discolouration.



AMG015300

SMALL

AMG015305

MEDIUM

AMG015310

LARGE

## KLEENSPEC® ECONOMY SERIES

- No sharp edges with this quality economy line
- Expendable illuminator can be re-sterilized up to 50 times before disposal
- No need for bulky exam lights



WELCH58600


SMALL

25/BX

WELCH58601

MEDIUM

25/BX

**KELLY FORCEPS**

AMG570848  
5.5" STRAIGHT


AMG570856  
5.5" CURVED

**DRESSING FORCEPS**

YMK-279-5 5"  
YMK-279M 4"

**LISTER BANDAGE SCISSORS**

AMG570308  
5.5" S/STEEL

**UNIVERSAL PARAMEDIC SCISSORS**

This paramedic scissor features hardened stainless steel blades and autoclavable plastic handles.

AMG120214  
ECONOMY 7.5"

AMG120210  
UNIVERSAL 8"

**STERILE SCISSORS**

BOW-119-316A  
4.5" SHARP/BLUNT

## MAGNA DISPOSABLE SCALPELS



10/bx


Magna Disposable Scalpels are individually packaged. Each scalpel blade is attached to a plastic handle, and comes with a blade cover for user protection.

ALMM92-10	SIZE 10
ALMM92-11	SIZE 11
ALMM92-15	SIZE 15
ALMM92-20	SIZE 20
ALMM92-21	SIZE 21
ALMM92-22	SIZE 22

## STERILE DISPOSABLE SCALPELS

ALMA6-510	SIZE 10
ALMA6-511	SIZE 11
ALMA6-515	SIZE 15
ALMA6-521	SIZE 21
ALMA6-523	SIZE 23

Almedic disposable Scalpels are manufactured from top quality surgical steel, sterilized by gamma radiation. Each scalpel blade is individually packed and supplied with a blade cover for user protection.



10/bx

## SCALPEL BLADES ONLY



Almedic Scalpel Blades are manufactured from the finest surgical quality stainless and carbon steel. Each blade is individually sealed in aluminum foil.

100/bx

## SCALPEL HANDLES



Reusable handles also available in stainless steel. Some sizes available in chrome, magna or utility.



## DEX4 GLUCOSE TABLETS

Dex4 Glucose Tablets contain 4 grams of fast-acting carbs per tablet that raises mild episodes of low blood sugar fast.

Great tasting orange, grape, raspberry, watermelon, and sour apple flavours to satisfy even the most discerning taste buds. What's more, they're fat-free and contain no sodium, caffeine or cholesterol!

Dex4® glucose tablets come in easy opening, crush-proof, travel tubes of 10 tablets, with a re-sealable snap-cap, so you can take them anywhere! Also available in economy-size refill bottles of 50.



AMG006318



AMG006313



AMG006315



AMG006314



AMG006316



AMG006312

## DEX4® GLUCOSE GEL

Raise Blood Sugar FAST, without chewing!


Great for moderate to severe lows, each tube contains 15g of fast-acting carbs (d-glucose) as recommended by the Canadian Diabetes Association (CDA).

- Glucose gel that raises blood sugar FAST & boosts Energy!
- Each tube contains 15 grams of fast-acting carbs (d-glucose) as recommended by CDA
- Incremental 5g dosing with gel level indicator
- Great for moderate & severe emergency lows
- Re-sealable, screw-on cap ensures sterility for multiple dosing
- Easy to swallow



AMG006337

## SAFETY LANCETS



Loris Safety Lancets offer a ready to use, easy and safe method for obtaining blood samples. The one step activation technique requires little training or technique and the automatic retraction of the needle protects patients and health care professionals from needle-stick injuries and eliminates the risk of cross contamination.

LERNA201-01	1.4MM	25G	PURPLE
LERNA202-01	1.6MM	28G	ORANGE
LERNA203-01	1.8MM	26G	YELLOW
LERNA204-01	2.2MM	21G	GREEN

## REGULAR LANCETS



LERNA205-01      200/BX

## MEDI+SURE BLOOD GLUCOSE METER

- Very small samples required - only 0.6uL - means less pain and no re-tests.
- Fast Results - Only 6 seconds.
- A Large Memory holding up to 960 recordings for historical perspective.
- No Coding Required.
- Simple and really easy to use.
- Large buttons and large visual screen for ease of use and reading.
- A lifetime replacement warranty on all meters.
- Easy to use computer uplink for storage of test results and sharing with health care providers.
- Rugged and Impact resistant while being small and portable.

**Meters are free upon the purchase of 100 Medi+Sure test strips.**



The Medi+Sure Test Strips are easy to handle and have highly accurate assays that (unlike other manufacturers' test strips) are not affected by interferences in other test strips such as high uric acid levels, or Peritoneal patients in dialysis, or higher blood concentrations of galactose, maltose, maltotriose or maltotetraose.

Our test strips are conveniently packaged in boxes of 100 containing two Vials of 50 test strips. Everyone Saves...

Quality and Performance that meets or exceeds the current Glucose meters, and strips in Canada at up to half the price of current test strips.

## COTTON TIPPED APPLICATORS NON-STERILE

3 inch wooden applicator stick with firmly wound cotton bud.

- Single Use Only
- 100% Cotton Tip

### 3" CTA'S NON-STERILE

AMG018450	100/PACK	10PK/BX
DUKAL9003	100/PACK	10PK/BX

### 6" CTA'S NON-STERILE

DUKAL9006	100/PACK	10PK/BX
-----------	----------	---------



## COTTON TIPPED APPLICATORS STERILE

The DUKAL applicator features a soft cotton tip that is absorbent and won't lint. Made from a high quality hard wood these applicators offer superior performance. The sterile 2's are available in an easy to peel package.



### 3" CTA'S STERILE

AMG018460	2/PK	100PK/BX
DUKAL9013	2/PK	100PK/BX

### 6" CTA'S STERILE

DUKAL9016	2/PK	100PK/BX
-----------	------	----------

### PLASTIC SHAFT

DUKAL9026	2/PK	100PK/BX
-----------	------	----------

## MEDPRO TONGUE DEPRESSORS STERILE



Individually wrapped tongue depressors. Sterilized and wrapped in peel open envelope.

AMG118425  
JUNIOR ST 100/BX

## MEDPRO TONGUE DEPRESSORS NON-STERILE



Ultra-smooth finish, 100% natural wood



AMG118414  
AMG118403

JUNIOR N/S  
SENIOR N/S

100/BX  
100/BX

Retail Packaging

AMG118415 SENIOR N/S 500/BX

## TONGUE DEPRESSORS



These hard wood products are clean and smooth and offer no odor or taste. They are exceptional for oral use, spreading ointments, mixing medications and even splints.

DUKAL9001 JUNIOR N/S 100/BX  
DUKAL9002 SENIOR N/S 500/BX

## FLEXIBLE STRAWS - WRAPPED



MedPro Flexible Straws, Wrapped  
19 cm x 6 mm

AMG018550	500/BX
BOWB02500W	500/BX

## FLEXIBLE STRAWS - UNWRAPPED



MedPro Flexible Straws,  
Unwrapped 19 cm x 6 mm

AMG018540	500/BX
-----------	--------

## COTTON BALLS



BOWW4010	MEDIUM	2000/BG
PACIFCL40	LARGE	1000/BG



**TABLE/EXAM PAPER****CREPE**

17"	WEST0101	125'	12/CS
18"	BOWB0118	125'	12/CS
18"	GRAHAM70002N	125'	12/CS
21"	GRAHAM70004N	125'	12/CS
24"	GRAHAM70007N	125'	12/CS

**SMOOTH**

18"	BOWB0218	225'	12/CS
18"	GRAHAM70016N	225'	12/CS
21"	GRAHAM70018N	225'	12/CS
24"	GRAHAM53216	225'	12/CS

**FACIAL TISSUE 2 PLY**

BOWB08300	100/BX	36BX/CS
BOWB08500	80/BX	135BX/CS

**WAXED PAPER CUPS**

BOWB113	4OZ	100/BG	10BG/CS
BOWB115	5OZ	100/BG	10BG/CS

**MEDICATION TEASPOONS****MEDICATION CUPS****Plastic**

BOWB2006			
1OZ	100/PK	50PK/CS	

**Paper**

BOWB100			
1OZ	250/PK	20PK/CS	



BOWB70034			
ECONO		1000/CS	

## ISOPROPYL ALCOHOL



Spill free dispenser design avoids risk of contamination and is spigot compatible for quick and easy pouring.

Ergonomic grip and compact bottle for easy handling and optimal storage.

70% and 99% concentrations available to suit all applications.

## LATEX TUBING



PRIME035R	LATEX TUBING 1/4x3/32	50'/RL
PRIME036R	LATEX TUBING 1/4x1/8	50'/RL
PRIME039R	LATEX TUBING 5/16x3/32	50'/RL
PRIME039RB	LATEX TUBING 5/16x3/32 BLACK	50'/RL
PRIME078R	LATEX TUBING 1/8x7/32	50'/RL
PRIME099BA-3	LATEX TUBING 3/16x1/16"	3ft/EA
PRIME102BA-6	LATEX TUBING 3/16x3/32"	6ft/EA
PRIME120XA	LATEX TUBING 1/4"x1/16"	100'/RL
PRIME122BA-14	LATEX TUBING 1/4"x3/32"	14ft/EA
PRIME122BA-6	LATEX TUBING 1/4"x3/32"	6ft/EA
PRIME122BA-9	LATEX TUBING 1/4"x3/32"	9ft/EA

## TOILET TISSUE



ENTSCO5131  
2 PLY 48RL/CS

## PAPER TOWELS



ENTKIM01960  
Multifold  
175/PK 25PK/CS

## TORK ELEVATION MATIC HAND TOWEL ROLL DISPENSER



Tork Elevation dispensers simplify and complement any environment, providing a clean, hygienic appearance. Our new design ensures easy and intuitive maintenance, refilling and use. The Tork Matic® Dispenser has a high capacity that is ideal for busy washrooms. Hands free dispensing provides a hygienic, one-at-a-time solution for your washroom. The user touches only the towel needed - reducing the risk of cross-contamination. Reliable and durable, Tork Matic® delivers one towel, time after time, always ready to service the next patron. Neutral, minimalist shape - a seamless capsule that fits in any environment. Semi-transparent cover and indication window makes it easy to see when a refill is needed saving maintenance time. Push bar for loading and emergency feeds. Utilizes any Tork H1 Roll Towel.

TORK5510202      WHITE  
TORK5510282      BLACK

## TORK UNIVERSAL HAND ROLL TOWEL

Designed specifically for the Tork Matic® and Tork Intuition® roll towel dispensers, this Tork Universal towel offers a truly best-in-class product for no touch, hygienic dispensing systems. Soft, strong, and highly absorbent this towel offers the perfect combination of quality, performance, and value. Embossing enhances hand feel and maximum absorbency; fewer towels are used per hand dry. Less packaging means source reduction reducing your overall labor and maintenance costs. Long lasting rolls minimize maintenance costs while maximizing service. Towels are made from 100% recycled fiber and a minimum of 50% Post Consumer fibers providing environmental benefits, and a positive image builder.

EcoLogo Certified and Green Seal Certified.

Compostable Verified by Independent Verification.



TORK290088

## TORK ADVANCED HAND ROLL TOWEL

Designed specifically for the Tork Matic® and Tork Intuition® roll towel dispensers, this Tork Advanced towel offers a truly best-in-class product for no touch, hygienic dispensing systems. Soft, strong, and highly absorbent this towel offers the perfect combination of quality, performance, and value. Embossing enhances hand feel and maximum absorbency; fewer towels are used per hand dry. Less packaging means source reduction reducing your overall labor and maintenance costs. Long lasting rolls minimize maintenance costs while maximizing service. Towels are made from 100% recycled fiber provides environmental benefits, and a positive image builder. EcoLogo Certified.



TORK290089

## TORK ADVANCED HAND TOWEL INTERFOLD



The unique 4 panel interfolded system, creates ready-to-use towels with no towel "Fall Out". Portion control, delivers one towel every time. Tork Advanced towels are easy to refill, reducing maintenance time & costs. Towels are made from 100% recycled fibers providing environmental benefits, and a positive image builder.

EcoLogo Certified.

TORK101290



The unique 3 panel interfolded system, creates ready-to-use towels with no towel "Fall Out". Portion control, delivers one towel every time. Tork Advanced Xpress™ towels are easy to refill, reducing maintenance time & costs. Towels are made from 100% recycled fibers (EPA compliant) providing environmental benefits, and a positive image builder.

EcoLogo Certified.

TORK101291



Tork Elevation Mini Xpress  
Hand Towel Interfold  
Dispenser



Tork Elevation Xpress  
Hand Towel Interfold  
Dispenser



Tork Xpress Hand Towel  
Interfold Dispenser



Tork Xpress Hand Towel  
Portable Interfold  
Dispenser

## TORK ELEVATION BATH TISSUE JUMBO ROLL MINI TWIN DISPENSER

Tork Elevation dispensers simplify and complement any environment, providing a clean, hygienic appearance. Our new design ensures easy and intuitive maintenance, refilling and use. The compact, yet high-capacity Tork Mini Twin Bath Dispensing System is the ideal choice for high-traffic washrooms as well as environments where efficiency and reduced cost are priorities. Unique tear feature eliminates risk of cutting or scratching hands. Rolls are protected so there is less waste, improved hygiene and reduced pilferage. Durable polystyrene dispensers provide low maintenance and trouble free service.



TORK5555200

WHITE

TORK5555290

BLACK

## TORK ADVANCED BATH TISSUE MINI JUMBO ROLL, 2-PLY, 7.36 INCH DIA.

Tork Advanced jumbo bath tissue's high capacity provides fewer refills and reduces labor. Large size is pilfer proof, hard to conceal. Can't use at home. Reduced stub roll means less waste. Fewer cores and no wrappers leads to source reduction. Core is made from 100% recycled fibers. Less storage space is needed. Rapid break-up capability is ideal for all plumbing systems. Made from 100% recycled fiber and 25% Post Consumer fiber providing environmental benefits, and a positive image builder.



EcoLogo Certified and Green Seal™ Certified.

TORK120246



## TORK ELEVATION HIGH CAPACITY BATH TISSUE ROLL DISPENSER

Tork Elevation dispensers simplify and complement any environment, providing a clean, hygienic appearance. Our new design ensures easy and intuitive maintenance, refilling and use. The Tork High-Capacity Bath Tissue Dispensing System is a perfect fit for washrooms with limited space. The innovative auto-transfer feature promotes hygiene and reduces refill maintenance time. Innovative and Hygienic Design - Promotes hygiene with covered roll and auto-transfer features. Prevents waste by dropping the reserve roll only after the primary roll is depleted. Keeps used roll core in the unit, keeping washroom floors clear of debris. High Capacity - Holds up to 2,000 two-ply sheets or 4,000 one-ply sheets. Contains 60 percent more tissue than conventional bath tissue dispensers. Reduced Costs - Permits quicker, easier loading and fewer refills, reducing maintenance time. Prevents pilferage and reduces waste with enclosed design. Requires rolls with core plugs.



TORK555620  
TORK555628

WHITE  
BLACK

## TORK ADVANCED HIGH-CAPACITY BATH TISSUE ROLL

Ideal for image to everyday washrooms, Tork Advanced High-Capacity Bath Tissue provides hygienic dispensing while reducing waste and maintenance expense. High capacity rolls eliminate service outages. Embossed rolls enhance bulk and softness. The banded rolls create less waste and source reduction. The easy start tail seal prevents waste. Tork bath tissue has rapid breakup, ideal for all plumbing systems. Made from 100% recycled fibers an environmental benefit, and a positive image builder.

EcoLogo Certified.



TORK110292A

## TORK TOILET SEAT COVER DISPENSER



This Tork toilet seat cover dispenser is made from white ABS Plastic, providing a bright, clean image.

TORK344080	WHITE
TORK344088A	BLACK

## TORK UNIVERSAL TOILET SEAT COVER, 1/2 FOLD



Tork Universal Seat Cover is a 1/2 fold design and is an easily dispensable, ready for use format. Hygienic, protects users from cross contamination.

TORKTC0020	25/BX
------------	-------

## TORK PREMIUM FACIAL TISSUE FLAT BOX



Tork Premium ultra soft, premium 2-ply tissue tells the user that you care about their comfort. Different color indicator sheet shows when the tissue box is nearly empty. Soft protective poly shield opening protects the tissues and provides easy, one-at-a-time dispensing. Dispenser box is manufactured from 100% recycled materials with a minimum of 35% post-consumer materials, a positive image builder.

TORKTF6920	100/BX
------------	--------

## TORK XPRESSNAP® TABLETOP NAPKIN DISPENSER



Since 2003 Tork Xpressnap® Classic has served up built-in value from foodservice counters across the country. Xpressnap Classic greatly increases personal hygiene and cuts waste by delivering one napkin at a time, every time. Customers touch only the napkins they take and take only what they use. Now available in black, granite, blue, red and yellow, Xpressnap® Tabletop dispensers also feature the AD-a-Glance message panels that can be used for in-store promotions or to generate ad revenue. The dispensers have a generous capacity of 200 napkins and measure 7.8 inches long by 5.8 inches wide and 6.2 inches tall.

TORK32XPT

## TORK ADVANCED XPRESSNAP® DISPENSER NAPKIN, INTERFOLD



Tork Advanced Xpressnap® system dramatically reduces napkin usage. Its one-at-a-time napkin dispensing design ensures that patrons use fewer napkins and touch only the napkin they use improving personnel hygiene. Plus, one-at-a-time dispensing helps decrease consumption, so you'll reduce the number of cases you need to store. Xpressnap napkins are larger than major competitor interfolded napkins offering strength, absorbency and quality to users. One napkin-at-a-time dispensing provides up to a 25% reduction in napkin usage. Smaller than standard napkin cases offers storage savings of up to 50%. Inner bag packaging ensures superior product quality while protecting unused napkins. This product contains 100% recycled fibers.

TORKDX900

500/BG

EcoLogo Certified.

**HEMECARE SUPPLIES**

AMG Medical  
Drive Medical  
Healthcraft  
Bowers Medical

<http://www.amgmedical.com/>  
<https://www.drivemedical.com/catalog/>  
<http://www.healthcraftproducts.com/>  
<http://www.bowersmedical.com/>

**MOBILITY AIDS**

Evolution Walkers

<http://www.evolutionwalker.com/>

**PATIENT SAFETY**

Secure Personal Safety  
Skilcare  
Posey

<http://www.padalarm.com/>  
<http://www.skil-care.com/>  
<http://www.posey.com/>

**BRACES & SUPPORTS**

Trulife  
Trainers Choice  
Anatech

<http://trulife.com/>  
<http://www.trainerschoice.ca/>  
<http://www.anatechinc.com/>

**COMPRESSION STOCKINGS**

Sigvaris  
Jobst

<http://www.sigvaris.ca/>  
<http://www.jobst.com/>

**PRESSURE RELIEF**

EHOB

<http://www.ehob.com/>

**EQUIPMENT**

Invacare  
MC Healthcare

<http://www.invacare.ca/>  
<http://mchealthcare.com/>

**INCONTINENCE**

Tena

<http://www.tena.ca/>

**AIDS FOR DAILY LIVING**

Parsons ADL

<http://www.parsonsadl.com/>

**WOUND CARE & FIRST AID**

Smith & Nephew  
Systagenix  
Molnlycke  
BSN  
Dukal  
Coloplast  
North Safety  
Wasip

<http://www.smith-nephew.com/canada/>  
<http://www.systagenix.ca/>  
<http://www.molnlycke.com/ca-en/>  
<http://www.bsnmedical.ca/>  
<http://www.dukal.com/>  
<http://www.coloplast.ca/en-ca/>  
<http://www.northsafety.ca/>  
<http://www.wasipltd.com/>

**SYRINGES & NEEDLES**

Terumo  
BD

<http://www.terumomedical.com/>  
<http://www.bd.com/ca/>

**INFECTION CONTROL**

Metrex  
Virox  
Ansell  
Medicom  
AMD Ritmed  
Bowman Dispensers

<http://www.metrex.com/>  
<http://www.virox.com/>  
<http://www.ansellcanada.ca/>  
<http://www.medicom.ca/en>  
<http://www.amdritmed.com/>  
<http://www.bowmandispensers.com/>

**CATHETERS & SUPPLIES**

Bard Medical  
Med-Rx  
Rochester Medical  
Rusch  
Amsino

<http://www.bardmedical.com/>  
<http://www.med-rx.ca/>  
<http://www.rocm.com/>  
<http://www.myrusch.com/frontend/bin/start.htm>  
<http://www.amsino.com/>

**DIAGNOSTICS**

American Diagnostics  
Welch Allyn

<http://www.adctoday.com/>  
<http://www.welchallyn.com/default.htm>



**REHAB & FITNESS THERAPY**

Remington Medical  
Patterson Medical  
Parker Labs

<http://www.pattersonmedical.ca/>  
<http://www.remingtonmedical.com/>  
[www.parkerlabs.com/](http://www.parkerlabs.com/)

**MISCELLANEOUS**

Covidien  
Mic-Key (Kimberly-Clarke)  
Qosina  
3M  
Medical Action  
Salter Labs  
Xodus Medical  
Key Surgical  
Texpro

<http://www.covidien.com/>  
<http://www.mic-key.com/home.aspx>  
<http://www.qosina.com/>  
<http://www.3m.com/intl/ca/>  
<http://www.medical-action.com/>  
<http://www.salterlabs.com/>  
<http://www.xodusmedical.com/>  
<http://www.keysurgical.com/>  
<http://www.texpro.net/>

## NOTES

[illegible]The logo for Bowers Home Health features the word "BOWERS" in a large, white, outlined serif font. Below it, the words "HOME HEALTH" are written in a smaller, white, solid sans-serif font. A white horizontal line runs across the background behind the text.

**Richmond: 1-800-663-0047**  
**Edmonton: 1-800-561-3687**

**www.bowersmedical.com**

## NOTES

[illegible]

**BOWERS**  
HOME HEALTH

**Richmond: 1-800-663-0047**  
**Edmonton: 1-800-561-3687**

**www.bowersmedical.com**

## NOTES

[illegible]

**BOWERS**  
HOME HEALTH

**Richmond: 1-800-663-0047**  
**Edmonton: 1-800-561-3687**

**www.bowersmedical.com**



**BOWERS**  
HOME HEALTH

(A DIVISION OF RADION LABORATORIES)

[www.bowersmedical.com](http://www.bowersmedical.com)  
[homehealth@bowersmedical.com](mailto:homehealth@bowersmedical.com)

**Head Office**

Unit 7-3691 Viking Way  
Richmond, BC V6V 2J6  
Tel: (604) 278-7566  
Toll Free: (800) 663-0047  
Fax: (604) 278-7525

**Alberta**

11783-186th Street  
Edmonton, AB T5S 2Y2  
Tel: (780) 484-1666  
Toll free: (800)561-3687  
Fax: (780) 451-6762