



Roho Dry Flotation Mattress Overlay Section



- Helps heal advanced stage (III & IV) ischemic ulcers
- Friction and shear relief surface
- Zoned: four sections each 20" X 34" / 51 cm X 86 cm
- 3.5" high air cells can be adjusted to the patient's body shape and size
- Durable neoprene rubber is flame-resistant and can be easily cleaned
- Reusable covers are available

Leveling pads and cushions are also available

ADVHEA3400

LTC T-Style Configurable Mattress



- Best choice for individuals at risk for pressure sore development at the coccyx and heels
- Four removable foam sections to accommodate the placement of Roho Neoprene (or like) inserts
- Top surface comprised of medium-firm viscoelastic foam
- Medium-firm polyurethane foam
- Transfer borders add firmness along sides
- Fluid-proof, non-allergenic, antibacterial, anti-microbial cover

EMPLTCAIS-TS-50



Made in Canada

Ultimate Comfort Premium Therapeutic Mattress



- Side safety borders add to safety and security
- Ultra-soft memory foam heel slope sculpted to redistribute pressure and reduce shear
- Fluid proof cover: non allergenic, antibacterial and phthalate-free
- Mattress cover: wraps fully down the sides and includes hidden zipper to eliminate moisture penetration
- Ca Tech Bulletin #117: meets 2013 flammability requirements
- 5 year mattress warranty

BOWESS-3580



Made in Canada

Intensive Therapeutics HeelZup Cloud Reusable



- Suspends the heels a minimum of one inch off of the surface
- Provides a high degree of comfort in both the supine and 30-degree side laying positions
- Stays neatly in place and does not require constant maintenance and adjustment
- No-Slip bottom material helps the cushion to stay neatly in place
- Manufactured of high-density, anti-microbial, hypo-allergenic VE foam, latex free materials and thermobonded polyester fiber
- Standard width is 30"

INT-CLOUD-XSR EXTRA SMALL
INT-CLOUD-SR SMALL
INT-CLOUD-MR MEDIUM
INT-CLOUD-LR LARGE

Head Office Tel: (604) 946-7712 Toll Free: (800) 663-0047

Alberta Tel: (780) 454-1666 Toll Free: (800) 561-3687

Ontario Tel: (905) 315-7455 Toll Free: (888) 216-4833

info@bowersmedical.com
www.bowersmedical.com



Single Medium Duty Steel Hamper



- Available in 37", 35" and 32" height
- Can accommodate nylon and poly liners
- Lid and foot pedals
- Solid 7/8" round tubing
- Bright chrome finish
- Durable 3" casters

PAC692

Non-Contact Infrared Thermometer



- Can be used in household environments as well as doctors in clinic settings
- Response Time: 1 second
- Measurement Range: 32-43°C
- Temperature Unit: °C/°F
- Measure Distance: 0-5cm
- Automatic Power- Off: 30 seconds
- Dimension: 149 X 38 X 42mm
- Batteries: 2 X AAA

PM-YT-1

Edan Vital Sign Monitor



- NIBP, EDAN SP02
- Portable and lightweight
- Large high resolution colour TFT display
- Pitch tones
- L-ion battery
- Real-time measurements display
- Nurse call function
- Size: 5.43" L x 3.7" W x 1.47" H

Mobile stand and accessories are also available

SOLIC-302112584

Edan H10B Finger Pulse Oximeter



- Small and light weight design
- Dual colour LED
- (2) AAA batteries
- Low power consumption automatic power off
- Low battery indicator
- SpO2, PR, bar graph, plethysmogram

SOLIC-H10B

Head Office Tel: (604) 946-7712 Toll Free: (800) 663-0047

Alberta Tel: (780) 454-1666 Toll Free: (800) 561-3687

Ontario Tel: (905) 315-7455 Toll Free: (888) 216-4833

info@bowersmedical.com
www.bowersmedical.com



Aquatec Ocean Ergo Tilt Shower Commode



- Piston free tilting mechanism reduces effort required to position center of gravity
- Rust proof, stainless steel frame
- Tool-free seat, height adjustment range of 5" to fit most toilets

INVA-OCEANERGOVIP

Aquatec Ocean Ergo VIP Tilt Shower Commode Seat



- Features key-shaped cut out
- Compatible with all Ocean Ergo shower and commode chair models
- White

INVA-AP1603142

Invacare Perfecto2 V Oxygen Concentrator



This 5-liter concentrator is small, quiet, light weight and comes standard with the Invacare SensO2 Oxygen Sensor. The Perfecto2 V unit incorporates many of the best features of the market-leading Invacare Perfecto2 Concentrator, including it's contemporary, patient preferred external design and colour.

- Top handle for ease of mobility
- Simplified preventative maintenance schedule
- 21" H x 13" W x 11.5" D
- Alarm and safety systems
- Compatible with many Invacare Platinum Oxygen concentrator parts
- Low flow alarm
- Battery-free power loss alarm
- Compressor 35 PSI pressure relief valve

IRC5PO2VC

Vacu-Aide Suction Unit



- High performance- vacuum adjustments allow for 80-55mm HG and a free flow of 27 LPM
- Adjustable flow regulator
- The Vacu-Aide series meets IOS- 10079 1;1999 suction

DRIVE7305P-D

Drive Medical Non Tilt Overbed Table



- Walnut, wood grain low-pressure laminate top
- Table top can be lowered to infinite positions between 28" - 48"
- H base provides security and stability
- Table top locks securely when height adjustment handle is released
- Can be raised with slightest upward pressure
- 2" swivel casters allow for easy maneuverability

DRIVE13067

Head Office Tel: (604) 946-7712 Toll Free: (800) 663-0047

Alberta Tel: (780) 454-1666 Toll Free: (800) 561-3687

Ontario Tel: (905) 315-7455 Toll Free: (888) 216-4833

info@bowersmedical.com
www.bowersmedical.com



BOWERS
MEDICAL SUPPLY
(A DIVISION OF RADION LABORATORIES)

ISOLATION CARTS

Storage For Personal Protective Equipment (PPE).



HIPRO-202843 (4-Drawer Cart) PREMIUM

- 20" D x 28" W x 43" H (with dolly, marine edge and casters)
- White body – HDPE
- 4 drawers : (1) 3" (2) 8" (1) 14"
- Yellow face
- 5" casters (2 swivel/ 2 brake in the front)

**NATIONWIDE
DELIVERY
AVAILABLE**



HIPRO-172236 (3-Drawer Cart) ECONOMY

- 18" D x 22" W x 36" H (with marine edge and casters)
- White body – HDPE
- 3 drawers : (2) 8" (1) 12"
- Yellow face
- 5" casters (2 swivel/ 2 brake in the front)



Utility Cart

- Ideal for an office, restaurant, healthcare facility or hotel
 - Strong sturdy handles
- Accommodates up to 200lbs of total weight
 - 40"L x 19.75"W x 37" H

ENTGCP5002

Ask us about our **NEW** Bowers Canadian N95 Equivalent Respirators (MLL-1970E),
and the Bowers Clearway Sanitizer Holders (BOWCLW-SH).

 **MADE IN CANADA**

Head Office Tel: (604) 946-7712 Toll Free: (800) 663-0047

Alberta Tel: (780) 454-1666 Toll Free: (800) 561-3687

Ontario Tel: (905) 315-7455 Toll Free: (888) 216-4833

info@bowersmedical.com
www.bowersmedical.com