



SOLUTIONS FOR TODAY...

...STRATEGIES FOR TOMORROW

FOOTCARE CATALOGUE

2019 - 2020



BOWERS
MEDICAL SUPPLY
(A DIVISION OF RADION LABORATORIES)

BOWERS LOCATIONS



BRITISH COLUMBIA

7198 Progress Way
Delta, BC V4G 1J2
Tel: 604.946.7712
Toll Free: 1.800.663.0047
Fax: 604.946.7715

ALBERTA

11783 - 186th STREET
Edmonton, AB T5S 2Y2
Tel: 780.454.1666
Toll Free: 1.800.561.3687
Fax: 780.451.6762

ONTARIO

1124 Northside Road
Burlington, ON L7M 1H4
Tel: 905.315.7455
Toll Free: 1.888.216.4833
Fax: 905.315.7422

Hours of Operation:

BC

Monday - Friday 8 am - 4:30 pm

Closed weekends & statutory holidays

Alberta



Monday - Friday 8 am - 4:30 pm

Closed weekends & statutory holidays

Ontario

Monday - Friday 8 am - 4:30 pm

Closed weekends & statutory holidays



Email us at:
info@bowersmedical.com

How to Order:

Placing an order with Bowers Medical Supply is easy and convenient. We have customer service agents across Canada waiting to assist you by phone, fax, email, online, B2B and EDI. Bowers Medical Supply has over 47 years of experience in the health care industry and we are pleased to provide the foot care community with their product and service requirements.

No account required or minimum order amount to meet.
Some orders might be eligible for no charge delivery.

Questions? Call or email us!
info@bowersmedical.com

Terms & Conditions:

- Returns:** Return requests must be authorized through your local Bowers Home Health location. All products must be returned within 30 days of receipt of authorization. Products may be subject to restocking fees. Custom and special order products are not returnable.
- Defects & Damages:** Claims for defects or shortages must be made within 24 hours of receipt of goods
- Payment:** Equipment and special orders may be subject to prepayment or deposit. All delivered orders must be paid prior to shipping. Applicable freight and taxes are extra. We accept all major credit cards. Debit & cash accepted for pick ups.



Follow us for the latest deals, events and tradeshow information.

TABLE OF CONTENTS

GENERAL SUPPLIES

Cotton Tipped Applicators.....	1
Tongue Depressors & Cotton Balls.....	2

WOUND CARE

Gauze (ST & N/S) & Conforming	3
Adhesive Bandages.....	4
Gypsona.....	5
Allevyn.....	6
Jelonet, Bactigras, Iodosorb, Inadine & Dressing Trays.....	7
Intrasite, Ichthopaste & Acticoat.....	8
Tegaderm, Duoderm, Adaptic & Aquacel.....	9
Tubular Gauze & Opsite.....	10

TAPES

	11
Paper, Transparent & Hypoallergenic.....	12
Athletic & Specialty Tapes.....	13-14
Cohesive Bandages.....	

SKIN CARE

Skin Preps & Antiseptics.....	15 - 16
Wound Cleansers.....	17
Barrier Creams & Moisturizers.....	18

GLOVES

Exam Gloves - Sterile & Non-Sterile.....	19 - 24
Glove Box Dispensers.....	25 - 27

PERSONAL PROTECTION

Masks.....	28
------------	----

HAND HYGIENE

Manorapid, Dermotan & VioNex.....	29
Hand Sanitizers & Accessories.....	30 - 35

DISINFECTANTS

Surface Disinfectants.....	36 - 38
----------------------------	---------

INSTRUMENTS

Instrument Lubricant	39
Scalpel Blades & Handles.....	40
Files & Rasps.....	41
Spatulas, Packers & Nail Splitters/Shavers.....	42
Curettes, Scissors & Forceps.....	43 - 45
Nail Nippers/Clippers, Brushes & Accessories.....	46 - 47

INSTRUMENT TRAYS

Instrument & Sterilization Trays.....	48 - 51
---------------------------------------	---------

STERILIZATION

Sterilization Pouches & Reels.....	52
Indicator Strips.....	53
CSR Wrap, Distilled Water, Indicator Incubator & Strips.....	54
Autoclaves.....	55 - 56

DIAGNOSTICS

Dopplers.....	57
---------------	----

BRACING & SPLINTS

Post-op Shoes & Cast Boots.....	58 - 59
Night Splints & Hallux Valgus Support.....	60

PRESSURE RELIEF

Heel Offloading Boots.....	61 - 62
----------------------------	---------

FOOT ACCESSORIES

Gel Lined Tubing & Toe Separators/Spreaders.....	63 - 64
Felts & Moleskin.....	65

COMPRESSION THERAPY

Compression Garments & Accessories.....	66 - 68
Profore, Surepress & Comprilan.....	69

ODOR CONTROL

Odor Eliminators.....	70
-----------------------	----

TOWELS & UNDERPADS

Podiatry & Tissue Towels.....	71
Disposable Absorbant Underpads.....	72

PAPER PRODUCTS

Exam Paper, Facial/Toilet Tissue & Paper Towels.....	73
--	----

COTTON TIPPED APPLICATORS NON-STERILE



3 inch wooden applicator stick with firmly wound cotton bud.

- Single Use Only
- 100% Cotton Tip

3" CTA'S NON-STERILE

AMG018450	100/PACK	10PK/BX
DUKAL9003	100/PACK	10PK/BX

6" CTA'S NON-STERILE

DUKAL9006	100/PACK	10PK/BX
-----------	----------	---------

COTTON TIPPED APPLICATORS STERILE

The DUKAL applicator features a soft cotton tip that is absorbent and won't lint. Made from a high quality hard wood these applicators offer superior performance. The sterile 2's are available in an easy to peel package.



3" CTA'S STERILE

AMG018460	2/PK	100PK/BX
DUKAL9013	2/PK	100PK/BX

6" CTA'S STERILE

DUKAL9016	2/PK	100PK/BX
-----------	------	----------

PLASTIC SHAFT

DUKAL9026	2/PK	100PK/BX
-----------	------	----------

MEDPRO TONGUE DEPRESSORS STERILE

Individually wrapped tongue depressors. Sterilized and wrapped in peel open envelope.



AMG118425
JUNIOR ST 100/BX

MEDPRO TONGUE DEPRESSORS NON-STERILE

Ultra-smooth finish, 100% natural wood



AMG118414	JUNIOR N/S	100/BX
AMG118403	SENIOR N/S	100/BX

Retail Packaging

AMG118415	SENIOR N/S	500/BX
-----------	------------	--------

TONGUE DEPRESSORS

These hard wood products are clean and smooth and offer no odor or taste. They are exceptional for oral use, spreading ointments, mixing medications and even splints.



DUKAL9001	JUNIOR N/S	100/BX
DUKAL9002	SENIOR N/S	500/BX

COTTON BALLS

BOWW4010	MEDIUM	2000/BG
PACIFCL40	LARGE	1000/BG





GAUZE

NON-STERILE

BOWB6112	2"X2"	200/BG
BOWB6113	3"X3"	200/BG
BOWB6114	4"X4"	200/BG

STERILE

DUKAL6122	2"X2"	2/PK	25PK/BX
DUKAL6123	3"X3"	2/PK	25PK/BX
DUKAL6124	4"X4"	2/PK	25PK/BX



TELFA NON-ADHERENT PADS

KEND1961	2"X3"	100/BX
KEND2132	3"X4"	100/BX
KEND2138	3"X8"	50/BX

ABD PADS



KEND9190A	5"X9"	36/TRAY
KEND9194A	8"X10"	18/TRAY

COMBINE ROLL



DUKAL3001 8"X20YDS

CONFORMING GAUZE



DUKAL602D	2" CONFORM	12/BG
DUKAL603D	3" CONFORM	12/BG
DUKAL604D	4" CONFORM	12/BG
DUKAL606D	6" CONFORM	6/BG

ADHESIVE BANDAGES

- Air vented, non-stick nylon pad
- Sterilized
- Individually wrapped
- 100 per box

AMG018600	3/4" X 3"	100/BX
BOWW018600	3/4" X 3"	100/BX
AMG018605	1" X 3"	100/BX



ASSORTED SHEER

ACME1260033 80/BX

FLEXIBLE FABRIC STRIPS



ACME1595033
1 X 3" 100/BX

SHEER SPOT



ACME1307033
7/8" 100/BX

SHEER PATCHES



ACME1308033
1.5" X 1.5" 100/BX

FABRIC KNUCKLE



ACME1602033
100/BX

FABRIC PATCHES



ACME1617033
2 X 3" 100/BX

GYPSONA

The traditional rich and creamy texture of Gypsona® allows intimate molding and creates a smooth, appealing finished cast. The unique interlocking leno weave greatly enhances stability allowing more plaster to stay on the gauze. The result is a high quality premium plaster of Paris making Gypsona® the world's leading brand.



Features and Benefits

- Leno gauze for extra stability
- High mass weight resulting in low plaster loss
- Creamy texture
- Central plastic core
- Latex Free

GYPSONA® - PLASTER BANDAGES - EXTRA FAST SETTING / 2 MINUTES

BSNLPL23	5 CM X 2.7 M	24 RL/BX	3 BX/CS
BSNLPL33	7.5 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL43	10 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL63	15 CM X 2.7 M	12 RL/BX	6 BX/CS
BSNLPL83	20 CM X 2.7 M	12 RL/BX	4 BX/CS
BSNLPL45	10 CM X 4.5 M	12 RL/BX	6 BX/CS
BSNLPL65	15 CM X 4.5 M	12 RL/BX	4 BX/CS
BSNLPL85	20 CM X 4.5 M	12 RL/BX	4 BX/CS

GYPSONA® - PLASTER SLABS - EXTRA FAST SETTING / 2 MINUTES

BSNL5430	7.5 CM X 75 CM	50 SPLINTS	6 BX/CS
BSNL5530	12.5 CM X 75 CM	50 SPLINTS	6 BX/CS

GYPSONA® - PLASTER BANDAGES - FAST SETTING / 4 MINUTES

BSN4000	5 CM X 2.7 M	24 RL/BX	3 BX/CS
BSN4010	7.5 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4040	10 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4045	15 CM X 2.7 M	12 RL/BX	6 BX/CS
BSN4020	10 CM X 4.5 M	12 RL/BX	6 BX/CS
BSN4030	15 CM X 4.5 M	12 RL/BX	4 BX/CS
BSN4050	20 CM X 4.5 M	12 RL/BX	4 BX/CS

ALLEVYN ADHESIVE

ALLEVYN ◊ Adhesive - Hydrocellular Dressing

- Low allergy adhesive
- Tapered construction with micro-thin edges
- Unique trilaminate structure
- Perforated wound contact layer
- Highly absorbent central hydrocellular layer
- Waterproof outer film layer

SMTH66000043	7.5CM X 7.5CM	10/BX
SMTH66000044	12.5CM X 12.5CM	10/BX



ALLEVYN NON-ADHESIVE

ALLEVYN◊ Hydrocellular Foam Dressing

- Unique trilaminate structure
- Non adherent wound contact layer
- Highly absorbent central hydrocellular layer
- Waterproof outer film layer
- Non adhesive

SMTH66007643	5CM X 5CM	10/BX
SMTH66007637	10CM X 10CM	10/BX
SMTH66000093	15CM X 15CM	10/BX



ALLEVYN GENTLE

ALLEVYN◊ Gentle Gel Adhesive Hydrocellular Foam Dressing

- Advanced triple-layered construction
- Absorbent hydrocellular pad
- Perforated soft gel adhesive wound contact layer
- Highly permeable waterproof outer film

SMTH66800247	5CM X 5CM	10/BX
SMTH66800248	10CM X 10CM	10/BX



ALLEVYN GENTLE BORDER

ALLEVYN ◊ Gentle Border Gel Adhesive Hydrocellular Foam Dressing with Border

SMTH66800269	7.5CM X 7.5CM	10/BX
SMTH66800270	10CM X 10CM	10/BX
SMTH66800272	12.5CM X 12.5CM	10/BX





JELONET

Indications

- Minor burns and scalds
- Donor and recipient graft sites
- Skin loss wounds, lacerations and abrasions
- Leg ulcers

SMTH7403	5CM X 5CM	50/BX
SMTH7404	10CM X 10CM	10/BX
SMTH7477	10CM X 40CM	10/BX



BACTIGRAS

Bactigras is indicated for adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers.

SMTH7457	10CM X 10CM	10/BX
SMTH7461	15CM X 20CM	10/BX



IODOSORB

PASTE

SMTH66060633	4CM X 6CM	5G	5/BX
SMTH66060634	6CM X 8CM	10G	3/BX

OINTMENT

SMTH66060630	10G TUBE	4/BX
SMTH66060631	20G TUBE	2/BX
SMTH66060632	40G TUBE	1/BX



INADINE

INADINE® PVP-I Non-Adherent Dressing consists of a low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base containing 10% Povidone Iodine; equivalent to 1.0% available iodine.

JOJOPO1481	5CM X 5CM	25/BX
JOJOPO1512	9.5CM X 9.5CM	25/BX

DRESSING TRAYS



- 1 Ea - Waste Bag W/ Tie
- 1 Ea - Drape Towel 45x55cm
- 1 Ea - Safety Pin
- 1 Ea - 3 Compartment Tray
- 1 Ea - Plastic Transfer Fcp
- 2 Ea - Metal Forceps
- 8 Ea - Gauze 10x10cm
- 8 Ea - Gauze 5x5cm

BOWTRDR101278B 40/CS



- 1 Ea - Drape Towel 45x55cm
- 2 Ea - Gauze 10x10cm
- 1 Ea - Waste Bag W/ Tie
- 1 Ea - 3 Compartment Tray
- 1 Ea - Plastic Transfer Fcp
- 6 Ea - Gauze 5x5cm
- 1 Ea - Plastic Locking Fcp

BOWTRDR6588 30/CS

****Custom trays are available upon request**

INTRASITE CONFORMABLE

IntraSite* Conformable is a conformable hydrogel dressing which combines the advantages of IntraSite Gel with a non-woven dressing to aid in the gentle packing of deep, open or undermined wounds. The non-woven substrate ensures the gel remains in close contact with the entire wound surface. It is suitable for use in necrotic, sloughy and granulating wounds, and helps create a moist wound environment to facilitate healing.



SMTH6600324	10CM X 10CM	10/BX
SMTH6600325	10CM X 20CM	10/BX
SMTH6600326	10CM X 40CM	10/BX

INTRASITE GEL



SMTH7308	8G
SMTH7311	15G
SMTH7313	25G

ICHTHOPASTE



SMTH4959 7.5CM X 6M

ACTICOAT

Acticoat* (with SILCRYST** Nanocrystals) Dressings are an effective barrier to bacterial penetration. The barrier function of the dressing may help reduce infection in partial and full-thickness wounds including pressure ulcers, venous ulcers, diabetic ulcers, surgical wounds, first and second degree burns and graft and dermal substitutes recipient sites. Acticoat* dressings may be used over debrided and grafted partial thickness wounds.



ACTICOAT

SMTH20101	4" X 4"	12/BX
-----------	---------	-------

ACTICOAT 7

SMTH20141	4" X 5"	5/BX
SMTH20241	6" X 6"	5/BX

ACTICOAT FLEX

SMTH66800396	5CM X 5CM	5/BX
SMTH66800399	10CM X 10CM	12/BX
SMTH66800409	10CM X 20CM	12/BX
SMTH66800419	20CM X 40CM	5/BX
SMTH66800432	40CM X 40CM	6/BX



TEGADERM DRESSING



A transparent film dressing designed for protecting skin and wound sites and securing devices to skin. It consists of a patterned coated film using Comfort Adhesive Technology for gentle adhesion with increased breathability and moisture exchange.

3M1622W	4.4CM X 4.4CM	100/BX
3M1624W	6CM X 7CM	100/BX
3M1626W	10CM X 12CM	50/BX
3M1627W	4CM X 10CM	20/BX

TEGADERM ROLL



3M16002	2" X 11YD	1 EA
3M16004	4" X 11YD	1 EA

DUODERM EXTRA THIN



The thin polyurethane film provides a bacterial¹ and viral² barrier. The film provides a waterproof barrier over the dressing. It can be removed without damaging newly formed tissue. It's Designed to reduce the risk of further skin breakdown due to friction by preventing contact with clothes/bed linen

CONV187955	4"X4"	10/BX
------------	-------	-------

AMD PACKING STRIP

Curity™ AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative microorganisms including some multi-drug resistant strains such as MRSA. This exciting new product provides healthcare facilities a low cost, prophylactic method of treatment without changing nursing protocol.



KEND7831AMD	1/4" X 1YD
KEND7832AMD	1/2" X 1YD
KEND7833AMD	1" X 1YD

SURGIFIX



SIZE 1	FRA1
SIZE 2	FRA2
SIZE 3	FRA3
SIZE 4	FRA4
SIZE 5	FRA5
SIZE 6	FRA6
SIZE 7	FRA7
SIZE 8	FRA8
SIZE 9	FRA9
SIZE 10	FRA10

TUBIGRIP



SIZE B	SANC1436
SIZE C	SANC1443
SIZE D	SANC1437
SIZE E	SANC1434
SIZE F	SANC1438
SIZE G	SANC1439

SURGITUBE



1 5/8"	GL105W
5/8"	GL110W

STOCKINETTE

2"	WEBSTO2
3"	WEBSTO3
4"	WEBSTO4
6"	WEBSTO6

OPSITE FLEXIFIX

The film is waterproof and a barrier to bacteria. Its high moisture vapour permeability allows the skin to breathe. The acrylic adhesive minimises the risk of skin damage on removal, especially after long periods of wear.

SMTH66000040	5CM X 10M
SMTH66000041	10CM X 10M



OPSITE FLEXIGRID

OpSite Flexigrid is a transparent, adhesive film, with a unique wound mapping grid. The film is moisture vapour permeable, conformable and extensible. It is clinically proven to promote moist wound healing for superficial wounds; it is also ideal for use as a secondary dressing and to help prevent pressure sores.

FLEXIGRID

SMTH4630	10CM X 12CM	50/BX
----------	-------------	-------

FLEXIGRID IV

SMTH4628	6CM X 7CM	100/BX
SMTH4631	15CM X 20CM	10/BX



MICROPORE TAPE



3M1530.5	1/2"	24/BX
3M15301	1"	12/BX
3M15302	2"	6/BX
3M15303	3"	4/BX

TRANSPORE TAPE



3M1527	1/2"	24/BX
3M15271	1"	12/BX
3M15272	2"	6/BX

PAPER TAPE



DUKALP50	1/2"	24/BX
DUKALP110	1"	12/BX
DUKALP210	2"	6/BX
DUKALP310	3"	4/BX

TRANSPARENT TAPE



DUKALT50	1/2"	24/BX
DUKALT110	1"	12/BX
DUKALT210	2"	6/BX
DUKALT310	3"	4/BX

MEFIX TAPE



SANC310250-05	2.5CM X 10M
SANC310500-05	5CM X 10M
SANC311000-04	10CM X 10M
SANC311500-05	15CM X 10M
SANC312000-01	20CM X 10M

HYPAFIX TAPE



BSN7144300	2.5CM X 10M
BSN7144301	5CM X 10M
BSN7144302	10CM X 10M
BSN7144303	15CM X 10 M

LEUKOTAPE

Leukotape® K is an elastic tape used to support the sensory functions. The stretchable cotton material works ideally with the skin's natural elasticity. The support is coated with a skin-friendly polyacrylate adhesive. It is permeable to air reducing the risk of maceration. The wave pattern of the adhesive can support the therapeutic effect. Leukotape® K is very durable, sports and showering present no problem.



LEUKOTAPE K

BSN7297811 5CM X 5M
BSN7297812 7.5CM X 5M

LEUKOTAPE P

BSN7616800 3.8CM X 13.7M

MEPITAC TAPE

Mepitac is highly conformable and easy to apply and remove. The Safetac layer assures secure fixation and gentle removal from skin and bodily hair. If inspection or repositioning is necessary, Mepitac can be lifted and re-applied. Mepitac is an ideal fixation solution for patients with fragile and sensitive skin. Mepitac is suitable for patients requiring repeated application and removal of tape over the same area e.g. neonates and dialysis patients.



SANC298300 2CM X 300CM
SANC298400-00 4CM X 1.5M

ATHLETIC TAPE



RENF20	1"	12/BX	
RENF30	1.5"	8/BX	
RENF40	2"	6/BX	
RENF50	3"	4/BX	
RENF50	4"	3/BX	
RENF101226	1.5" SECONDS		32/CS
JAYB300X-1515	1.5" SECONDS		100/CS

FOAM UNDERWRAP



JAYB50-27530
48/CS

3M™ COBAN™ SELF-ADHERENT WRAP

A self-adherent elastic wrap that functions like a tape, but sticks only to itself. Available in sterile and nonsterile styles, and in a variety of widths and colors to meet your application needs.



5 YDS/RL

3M1581	1"	30/CS
3M1582	2"	36/CS
3M1583	3"	24/CS
3M1584	4"	18/CS
3M1586	6"	12/CS

CO-PLUS®/CO-PLUS® LF ELASTIC COHESIVE BANDAGE

Co-Plus® flexible elastic cohesive bandage provides non-slip light support and compression. Because of its unique construction, Co-Plus® adheres to itself and not to skin, aiding in patient comfort. No clips or fasteners required. Available in latex and latex free versions.



Features and Benefits

- Adheres only to itself
- Non-slip support
- Breathable
- Hand-tearable (Co-Plus® LF only)
- Latex Free (Co-Plus® LF only)

BSN7210014	2.5 cm x 4.5 m* Beige	60 RL/BX
BSN7210016	5 cm x 4.5 m* Beige	36 RL/BX
BSN7210017	7.5 cm x 4.5 m* Beige	24 RL/BX
BSN7210018	10 cm x 4.5 m* Beige	18 RL/BX
BSN7210019	15 cm x 4.5 m* Beige	12 RL/BX
BSN7210033	5 cm x 4.5 m* Mixed Colors	36 RL/BX

COHESIVE BANDAGES

A lightweight compression bandage ideal for treatment of sprains and a variety of other injuries

- Great for relief bandaging of joints and other difficult-to-wrap areas
- No clips or fasteners needed – self-adherent wrap sticks to itself, not to other materials or skin
- Will not slip – layers “bond” together for the strongest holding power
- Easily torn without scissors
- Individually poly bagged
- Available with or without latex
- Available in tan, black, or in an assorted rainbow pack (red, green, purple, dark blue, pink & light blue)



* Unless otherwise stated product contains latex.

3171	Sensi Wrap, Self-Adherent - 1" x 5 yds Tan	30/Cs	
3172	Sensi Wrap, Self-Adherent - 2" x 5 yds Tan	36/Cs	
3173	Sensi Wrap, Self-Adherent - 3" x 5 yds Tan	24/Cs	
3174	Sensi Wrap, Self-Adherent - 4" x 5 yds Tan	18/Cs	
3176	Sensi Wrap, Self-Adherent - 6" x 5 yds Tan	12/Cs	
3181	Sensi Wrap, Self-Adherent - 1" x 5 yds Rainbow (5/color)	30/Cs	
3182	Sensi Wrap, Self-Adherent - 2" x 5 yds Rainbow (6/color)	36/Cs	
3183	Sensi Wrap, Self-Adherent - 3" x 5 yds Rainbow (4/color)	24/Cs	
3184	Sensi Wrap, Self-Adherent - 4" x 5 yds Rainbow (3/color)	18/Cs	
3186	Sensi Wrap, Self-Adherent - 6" x 5 yds Rainbow (2/color)	12/Cs	
3187	Sensi Wrap, Self-Adherent - Latex Free - 1" x 5 yds Tan	30/Cs	
3188	Sensi Wrap, Self-Adherent - Latex Free - 2" x 5 yds Tan	36/Cs	
3189	Sensi Wrap, Self-Adherent - Latex Free - 3" x 5 yds Tan	24/Cs	
3190	Sensi Wrap, Self-Adherent - Latex Free - 4" x 5 yds Tan	18/Cs	
3191	Sensi Wrap, Self-Adherent - Latex Free - 6" x 5 yds Tan	12/Cs	
3211	Sensi Wrap, Self-Adherent - 1" x 5 yds Black	30/Cs	
3212	Sensi Wrap, Self-Adherent - 2" x 5 yds Black	36/Cs	
3213	Sensi Wrap, Self-Adherent - 3" x 5 yds Black	24/Cs	
3214	Sensi Wrap, Self-Adherent - 4" x 5 yds Black	18/Cs	
3215	Sensi Wrap, Self-Adherent - Latex Free - 1" x 5 yds Rainbow (5/color)	30/Cs	
3216	Sensi Wrap, Self-Adherent - Latex Free - 2" x 5 yds Rainbow (6/color)	36/Cs	
3217	Sensi Wrap, Self-Adherent - Latex Free - 3" x 5 yds Rainbow (4/color)	24/Cs	
3218	Sensi Wrap, Self-Adherent - Latex Free - 4" x 5 yds Rainbow (3/color)	18/Cs	
3219	Sensi Wrap, Self-Adherent - Latex Free - 6" x 5 yds Rainbow (2/color)	12/Cs	



SKIN PREP

SMTH59420425 WIPES 50/BX
 SMTH420279 SPRAY 118ML

NO STING (ALCOHOL FREE)

SMTH59420600 WIPES 50/BX
 SMTH66800709 SPRAY 1OZ



CAVILON NO-STING BARRIER FILM

3M3344E WIPES 30/BX
 3M3346E SPRAY 28ML



REMOVE ADHESIVE REMOVER

SMTH403120 WIPES 50/BX
 SMTH403379 SPRAY 220ML



ADHESIVE REMOVER

PDIB16400 WIPES 100/BX

ALCOHOL PREP PADS



BOW3001 MED 100/BX
 BOW3002 LRG 100/BX

BENZALKONIUM WIPES



Benzalkonium
 antiseptic towelettes.

PDID35164 100/BX

IODINE PREP PADS

Iodine Prep pads 10%

PDIB40600 100/BX



POVIDONE IODINE SWABSTICKS

PDIS41125 3'S 25/BX



PVP 10% ANTISEPTIC SWAB SOLUTION

Each pad is saturated with 10% W/V Povidine Iodine solution, USP which is equivalent to 1% available Iodine
 For external use only
 Air dries quickly requiring no wiping
 Latex-Free
 For cleansing minor scrapes, cuts and burns
 Prepares skin by cleansing dead skin, dirt, or oil
 Helps reduce bacteria

LORIS-MEDPVP 200/BX



POVIDONE IODINE PREP 10% SOLUTION

BETADINE SOLUTION

UNI00452243 100 ML

UNI00184002 500 ML

LERNA109-09 115 ML

LERNA109-08 500 ML





SEA-CLENS WOUND CLEANSER

Cost-effective, saline-based wound cleanser designed to loosen debris from the wound bed without interfering with the healing process

- Saline-based cleanser provides safe and effective irrigation
- Soften eschar to facilitate debridement
- Non-toxic, non-ionic, pH-balanced formula

SWEEN1063 178ML
SWEEN1061 355ML



SALINE MODUDOSE

NACL 0.9% for irrigation.

BALLARD116
15ML 24/BX



SALINE SOLUTION

NACL 0.9% (Normal saline) for irrigation.

BAXTJF7633P 500ML
BAXTJF7634 1000ML



STERILE WATER

Sterile water for irrigation.

BAXTJF7623P 500ML
BAXTJF7624 1000ML

SPROAM® CLEANSER

- Complete head-to-toe cleansing, shampooing, bathing and incontinent skin care.
- Can be used as a spray or foam, for maximum versatility.
- Non-aerosol dense foaming action
- Dye and fragrance free
- pH-balanced

SWEEN0897 175ml
SWEEN0901 350ml



BEDSIDE-CARE® EasiCleanse™

Bedside-Care® EasiCleanse™ Bath is a no-rinse, self-sudsing, disposable washcloth.

- For the patient satisfaction of a sudsy bath without the basin
- Non-irritating, non-sensitizing
- Preservative-, dye-, alcohol- and latex-free
- pH-balanced
- CHG compatible

COLOPL7056

100 wipes/pkg



Bedside-Care® Cleanser



Bedside-Care® is a no-rinse, all-in-one body wash, shampoo and incontinence cleanser.

- Deodorizing
- Alcohol-free
- pH-balanced
- CHG compatible

SWEEN6176

120ml

Bedside-Care® Foam



Bedside-Care® Foam is a no-rinse, all-in-one, foaming body wash, shampoo and incontinence cleanser.

- Deodorizing
- Alcohol-free
- pH-balanced

COLOPL7300

125ml

SWEEN7146

240ML

Atrac-Tain® Cream



With Urea and Alpha Hydroxy Acid (AHA) Advanced moisturizing relief for severely afflicted dry, scaly, cracked skin.

Ideal for heels, elbows, knees, lower leg areas. Cream with 10% Urea and 4% AHA cream.

Cream with 10% Urea and 4% AHA cream.

Lotion with 5% Urea and 2.5% AHA.

SWEEN1843

2gm Sachet

SWEEN1802

60gm Cream

SWEEN1814

142gm Cream

Sween® 24 Cream



Once a day, long lasting 24 hour protection. Softens and soothes skin by sealing in moisture. Applies easily and absorbs quickly.

- Once a day application
- Contains 6% dimethicone
- Fragrance, petrolatum, lanolin and alcohol-free
- CHG compatible
- Suitable for neonatal use

SWEEN7090

VINYL GLOVES

DERMASHEER

- Vinyl Medical Examination Gloves
- Powder Free
- Latex Free

BOW1000S	SMALL
BOW1000M	MEDIUM
BOW1000L	LARGE
BOW1000XL	EXTRA LARGE




LATEX GLOVES

DERMASENSE

- Latex Medical Examination Gloves
- Textured
- Powder Free

BOW10110	EXTRA SMALL
BOW10111	SMALL
BOW10112	MEDIUM
BOW10113	LARGE



NITRILE GLOVES

DERMATUFF

- Nitrile Medical Examination Gloves
- Textured
- Powder Free
- Latex Free

BOW30310	EXTRA SMALL
BOW30311	SMALL
BOW30312	MEDIUM
BOW30313	LARGE
BOW30314	EXTRA LARGE



priMED PrimaTouch® Fit Nitrile Gloves

PM6-2005A	X-SMALL
PM6-2006A	SMALL
PM6-2007A	MEDIUM
PM6-2008A	LARGE
PM6-2009A	X-LARGE



- DEHP-Free
- Purple Colour
- Smooth finish
- For low-risk applications
- Powder free
- Latex Free
- Non-sterile

200/BX 10BX/CS

MICRO-TOUCH® DERMACLEAN®

ANS3180	XSMALL
ANS3181	SMALL
ANS3182	MEDIUM
ANS3183	LARGE
ANS3184	XLARGE



- Low latex protein*
- Cream color
- Textured finish
- Beaded cuff
- Length: 241 mm (9.5")
- Powder free
- Non-sterile

100/BX 10BX/CS

**XL = 80/BX

MICRO-TOUCH® AFFINITY®

ANS3770	XSMALL
ANS3771	SMALL
ANS3772	MEDIUM
ANS3773	LARGE
ANS3774	XLARGE



- Neoprene
- Green color
- Textured fingertips
- Beaded cuff
- Length: 239 mm (9.4")
- Powder free
- Non-sterile
- Easy donning due to polymer coating
- Tested for use with chemotherapy drugs

100/BX 10BX/CS

MICRO-TOUCH® NITRILE TEXTURED

- Nitrile
 - Blue color
 - Textured finish
 - Beaded cuff
 - Length: 240 mm (9.4")
 - Excellent chemical and puncture resistance
 - Tested for use with chemotherapy drugs.
- ASTM D6978



ANS4000	XSMALL
ANS4002	SMALL
ANS4004	MEDIUM
ANS4006	LARGE
ANS4008	XLARGE
<i>100/BX</i>	<i>10BX/CS</i>

MICRO-TOUCH® NITRILE

- Nitrile
 - Blue color
 - Thinner to enhance sensitivity
 - Textured fingertips
 - Beaded cuff
 - Excellent chemical and puncture resistance
 - Tested for use with chemotherapy drugs.
- ASTM D6978
- Length: 240 mm (9.4")



ANS6034300	XSMALL
ANS6034301	SMALL
ANS6034302	MEDIUM
ANS6034303	LARGE
ANS6034304	XLARGE
<i>200/BX</i>	<i>10BX/CS</i>

MICRO-TOUCH® PLUS

- Low latex protein*
- Cream color
- Textured finish
- Beaded cuff
- Length: 241 mm (9.5")
- Easy donning due to polymer coating



ANS6015300	XSMALL
ANS6015301	SMALL
ANS6015302	MEDIUM
ANS6015303	LARGE
ANS6015304	XLARGE
<i>150/BX</i>	<i>10BX/CS</i>

MICRO-TOUCH® NITRAFREE™

ANS6034510 XSMALL
 ANS6034511 SMALL
 ANS6034512 MEDIUM
 ANS6034513 LARGE
 ANS6034514 XLARGE



- Nitrile
- Pink color
- Chemical accelerator-free*
- Textured fingertips
- Beaded cuff
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs. ASTM D6978
- Length: 240 mm (9.4")

100/BX 10BX/CS



SYNSATION®

ANS6043300 XSMALL
 ANS6043301 SMALL
 ANS6043302 MEDIUM
 ANS6043303 LARGE
 ANS6043304 XLARGE




- Polyvinyl chloride
- DEHP-Free
- Clear
- Smooth finish
- Beaded cuff
- For low-risk applications
- Length: 230 mm (9.1")

200/BX 10BX/CS



WALL BRACKET



- Universal, fits most glove box sizes
- Wall-mountable
- Keeps glove organized and accessible

PM-BRACKET

DISPOS-A-GLOVE® STERILE

- Copolymer
- Clear
- Smooth finish
- Length: 269 mm (10.6")

JOJOMDG601	SMALL		
JOJOMDG701	MEDIUM		
JOJOMDG801	LARGE	100/BX	14BX/CS



VINYL GLOVES

- Vinyl
- Non-sterile
- Pre-Powdered

AMD9993A	SMALL		
AMD9993B	MEDIUM		
AMD9993C	LARGE		
AMD9993D	XLARGE	100/BX	10BX/CS



ACCUTOUCH NITRILE GLOVES

Our Accutouch Chemo Nitrile Exam Gloves are both latex and powder free.

Nitrile gloves conform comfortably to the shape of the hands and return to original shape after removal.

They feature excellent tactile sensitivity and high tensile strength for durable barrier protection.

Textured for a secure grip, with tight cuffs to minimize accidental exposure to blood-borne pathogens and other contaminated fluids.

Approved for use with chemotherapy drugs, formaldehyde and gluteraldehyde.

MEDLMDS192084	SMALL		
MEDLMDS192085	MEDIUM		
MEDLMDS192086	LARGE		
MEDLMDS192087	XLARGE		

100/BX 10BX/CS



CONFORM® PROCEDURE



- Cream color
- Bisque finish
- Unbeaded cuff
- Length: 312 mm (12.3")

ANS786960	6.0
ANS787065	6.5
ANS787170	7.0
ANS787275	7.5
ANS787380	8.0
ANS787485	8.5

1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

DERMA PRENE® ISOTOUCH®



- Polyisoprene
- Green color
- Smooth finish
- Beaded cuff
- Ergonomic for better comfort
- Tested for use with chemotherapy drugs. ASTM D6978
- Length: 306 mm (12")

ANS20685255	5.5
ANS20685260	6.0
ANS20685265	6.5
ANS20685270	7.0
ANS20685275	7.5
ANS20685280	8.0
ANS20685285	8.5
ANS20685290	9.0

1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

PERRY® STYLE 42®



- Cream color
- Smooth finish
- Beaded cuff
- Ergonomic for better comfort
- Length: 300 mm (11.8")
- Contains natural rubber latex

PER5711100	5.5
PER5711101	6.0
PER5711102	6.5
PER5711103	7.0
PER5711104	7.5
PER5711105	8.0
PER5711106	8.5
PER5711107	9.0

1 Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

GLOVE BOX DISPENSER - SINGLE

Glove Box Dispenser - Single with Built-in Flexible Spring
Holds one box of gloves

Two-way keyholes for vertical or horizontal wall mounting
White powder coated steel

BOWMAN-GB001



GLOVE BOX DISPENSER - DOUBLE

Glove Box Dispenser - Double
Holds two boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting
White powder coated steel

BOWMAN-GB002



GLOVE BOX DISPENSER - TRIPLE

Glove Box Dispenser - Triple
Holds three boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting
White powder coated steel

BOWMAN-GB003



GLOVE BOX DISPENSER - SINGLE WITH FLEXIBLE SPRING



Glove Box Dispenser - Single with Flexible Spring
Holds one box of gloves with built-in flexible spring

Two-way keyholes for vertical or horizontal wall mounting
Clear PETG plastic

BOWMAN-GP013

GLOVE BOX DISPENSER - DOUBLE



Glove Box Dispenser - Double
Holds two boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting
Clear PETG plastic

BOWMAN-GP014

GLOVE BOX DISPENSER - TRIPLE



Glove Box Dispenser - Triple
Holds three boxes of gloves

Two-way keyholes for vertical or horizontal wall mounting
Clear PETG plastic

BOWMAN-GP015

GLOVE BOX DISPENSER - QUAD

Glove Box Dispenser - Quad
Holds four boxes of gloves horizontally

Keyholes for wall mounting
Clear PETG plastic

BOWMAN-GP340



GLOVE BOX DISPENSER - SINGLE WITH FLEXIBLE SPRING

Glove Box Dispenser - Single with Flexible Spring
Holds one box of gloves with built-in flexible spring

Two-way keyholes for vertical or horizontal wall mounting
Stainless steel

BOWMAN-GS004



GLOVE BOX CADDY - SINGLE

Glove Box Caddy - Single
Holds one box of gloves

Attaches to the bottom of protection organizer LD-037
Quartz ABS plastic

BOWMAN-GD011



EARLOOP MASKS



BOWB2115 50/BX

SURGICAL TIE MASK



PRIMAGARD® Surgical Mask with Soft Horizontal Tie, Blue Swoosh

PM-PG4-2560

50/BX 300/CS

PROCEDURE EARLOOP W/VISOR



PRIMAGARD®160 Procedure Ear Loop with attached Visor, Anti-fog Foam, Anti-glare, Indigo

PM-PG4-1573

25/BX 100/CS

SURGICAL TIE W/VISOR



PRIMAGARD®120 Surgical Tie with attached Visor, Anti-fog Foam, Anti-glare, Denim Blue

PM-PG4-2592

25/BX 100/CS

3M™ VFlex™ Healthcare Particulate Respirator and Surgical Mask 1804, N95



- NIOSH approved N95
- FDA cleared for use as surgical mask
- Fluid resistant and disposable
- 3M proprietary filter media and pleats help make breathing easier

3M1804

50/BX 8BX/CS

MANORAPID

Manorapid is a highly effective poly alcohol formula that when rubbed on your hands for 20-30 seconds kills bacteria and fungi.

That is a lot of germs! 99.999.33% is not your average sanitizer.

Manorapid has a special lanolin-derivative emollient that ensures your skins stays soft and smooth after each application even under the harshest winter conditions found in Canada.

Manorapid is a proven formula used in scores of countries and for over 40 years has helped prevent infections in hospitals and nursing homes.



REMAS-HAND-120	120ML WITH FLIP CAP
REMMANO1000P	1 L WITH PUMP
REMMANO1000B	BRACKET FOR 1 L

DERMOTAN

Dermotan is an imported soap from Germany manufactured with caring substances found in high quality cosmetic products making it very gentle on the skin.

The smallest amount of Dermotan will clean the dirtiest of hands, making this a cost effective product.

Dermotan can also be used in bath water cleansing the most fragile skin; from new born babies to the elderly even with difficult skin problems.

Dermotan is used as a pre scrub prior to using Manorapid in the OR for a perfect combination to keep your hands from drying out even in the winter months when the air is dry and cold.



REMDERM1000	1 L
-------------	-----

VIONEX HAND SOAP

VioNex Antimicrobial Liquid Soap is specially formulated to be gentle and non-irritating for health care professionals who handwash 30-50 times per day. It is recommended for use between routine patient contact, procedures, after exposure to potential pathogens and after removal of personal protective equipment.

METRMX1518	18 OZ WITH PUMP
METRMX1500	1 G REFILL





PURELL® HAND SANITIZER The Brand People Know And Trust.

PURELL Advanced Hand Rub is the first and only sanitizer allowed for Commercial Use, Hospital and Healthcare Professional Use and in Food Premises in Canada.* That's because it is the first sanitizer to be granted an Exemption by Health Canada for Hospital and Healthcare Professional use, Commercial use and in Food Premises based on the 2009 Health Canada GUIDANCE DOCUMENT Human-Use Antiseptic Drugs.*

It provides proven persistent effect and kills 99.99% of germs* including:

- Norovirus (MNV)
- Staphylococcus aureus
- Escherichia coli
- Staphylococcus epidermis
- Serratia marcescens
- Candida albicans
- Aspergillus niger
- Rotavirus
- Rhinovirus
- Adenovirus

This green certified hand rub gel formulation contains ingredients made from natural renewable resources. Fragrance free. Meets EcoLogo hand sanitizer standard for environmental leadership and proven performance. Scientifically advanced, patent-pending 70% formulation outperforms other leading hand sanitizers milliliter for milliliter(1). Includes a nourishing blend of four different skin conditioners. Clinically proven to help maintain skin health(2). North America's #1 instant hand sanitizer.

PURELL® 70 INSTANT HAND SANITIZER

The power of PURELL® in a 70% ethyl alcohol, fragrance free formula. Effective in destroying certain bacteria. Acceptable for frequent use. Contains vitamin E and glycerin.



GOJO9770-12 354ml

PURELL® 70 INSTANT HAND SANITIZER MOISTURIZING FOAM

PURELL 70 Moisturizing Foam is specially formulated to be mild to the skin, non-drying and increase skin hydration.



GOJO5795-04 535ML



PURELL® Personal Gear
Retractable Clip

GOJO9608-24



GOJO5695-24
45ML

PURELL® VERSAHOLD™ POINT OF CARE BRACKET

Drive hand hygiene compliance by making PURELL Instant Hand Sanitizer available at the point of patient care. The PURELL VERSAHOLD bracket provides flexible options for attaching to a wide variety of fixtures.

PURELL VERSAHOLD Point of Care Bracket

- Addresses a key barrier to hand hygiene compliance: access to hand hygiene products
- Helps meet WHO and PIDAC compliance standards throughout the patient zone – in patient rooms, on mobile transport, in work areas, wherever needed
- Quick and easy on and off mounting capability, no tools required



PURELL® SANITIZING HAND WIPES



Enhance guest and patient experience and satisfaction with economical towelettes that clean and sanitize. Individually wrapped, “on-the-go” amenity item for a wide variety of markets, including foodservice, healthcare, travel and hospitality. Made with 62% ethyl alcohol – the same active ingredient as original formula PURELL Instant Hand Sanitizer.

GOJO9022-10 100/BX

PURELL® ALCOHOL HAND SANITIZING WIPES



For personal hand hygiene to help prevent the spread of certain bacteria •Use anytime for quick clean-up when soap and water are not available.

GOJO9031-06 175/TUB

PURELL® 70 INSTANT HAND SANITIZER MOISTURIZING

PURELL 70 Moisturizing Foam is specially formulated to be mild to the skin, non-drying and increase skin hydration.

TFX™ 1200 mL Refill

- Fresh dispensing valve with each refill
- Refill snaps in place with audible click



GOJO5395-02

PURELL® ADVANCED MOISTURIZING FOAM HAND RUB

1200 mL Refill

- GOJO SANITARY SEALED™ refill helps lock out germs
- Durable, recyclable PET bottle uses 30% less material than standard rigid HDPE
- Patent pending GOJO CONTROLLED COLLAPSE technology keeps bottle shape longer while emptying
- Crystal clarity for visual appeal and easier monitoring of fill level
- Fresh dispensing valve with each refill
- Removable pump promotes easier recycling



GOJO1913-02

GOJO1903-02 ADVANCED HAND RUB 1200ML

PURELL® LTX-12™ DISPENSER - WHITE

- Reliable touch-free dispensing
- Smart, trouble-free electronics eliminate battery changes in most installations
- Large sight window and skylight for at-a-glance product monitoring
- 1200mL high capacity output for fewer changeouts
- Easily converts to locked dispenser
- Lifetime guarantee even includes batteries
- Fully ADA compliant
- Customization available

GOJO1920-04



GOJO® CLEAR & MILD FOAM HANDWASH



Spa-inspired, green certified foam hand soap offers luxurious foam lather in a fragrance- and dye-free formulation. Enriched with moisturizers and skin conditioners. USDA certified biobased formulation. EcoLogo certified. Delights users while supporting a healthy environment.

GOJO1311-03

700ML

GOJO® LTX-7™ DISPENSER - WHITE



Compact 700mL dispenser is ideal for tight spaces. Lifetime Performance Guarantee; including batteries. GOJO SMART ELECTRONICS technology minimizes energy use per actuation, eliminating need to change batteries over the life of the system in most installations. Large sight window, skylight and clear refill bottles make it easy to check fill status. CONTROLLED COLLAPSE refills hold their shape as they empty for a more attractive appearance. Easily converts to locking dispenser. Removable pump makes recycling easy. SANITARY SEALED™ refills made of durable, recyclable PET material, using 30% less plastic than rigid HDPE bottles. Fully ADA compliant.

GOJO1380-04

GOJO® GREEN CERTIFIED FOAM HAND CLEANER



A mild, biodegradable foam soap formula certified by Green Seal™ and EcoLogo® to meet their joint standard (GS-41/CCD-104) for hand cleaners and hand soaps. The certifications mean that the product has been thoroughly evaluated by these unbiased, third-party organizations and meet standards for being environmentally responsible without sacrificing product performance

GOJO1380-04

ISAGEL HAND SANITIZING GEL

Ethyl Alcohol Antiseptic No-Rinse Hand Gel
70% Ethyl Alcohol

Effective, fast-acting handrub kills greater than 99.9% of most transient microorganisms (including MRSA and VRE) in 15 seconds or less. Enriched with emollients, Isagel is gentle on the skin. Compatible with gloves and CHGbased products.



SWEEN1644	118 ML
SWEEN1645	621 ML
SWEEN7041	800 ML

DISPENSER BRACKET

Made of very sturdy plastic.
Made specifically to hold 621 ml sized pump bottles. Theft resistant.



SWEEN0089

DISPENSER FOR BAG



Safe and compact, easy to dispense. Designed for 800 ml and 1000 ml sized bags.

SWEEN7251

DRIP TRAY



Accessory to be use with 7251 to prevent spillage.

SWEEN7262

ACCEL CS20



Accel Instrument Products provide the perfect balance between Safety and Efficacy by providing a safer, more responsible solution to disinfection without the use of harmful chemicals such as Aldehydes.

Accel CS 20 is designed for use in Healthcare and Personal Services when sterilization in a rapid contact time is required. Accel CS 20 achieves chemical sterilization and high level disinfection in 20 minutes and has a 14 day re-use claim.

VIROXCS20-1	1 LT
VIROXCS20-4	1 GALLON

ACCEL CONCENTRATE



Accel Concentrate is a Ready-to-use One-step disinfectant cleaner for use as a daily cleaner and as a disinfectant of environmental surfaces and patient care devices. Accel Concentrate uses the Accelerated Hydrogen Peroxide (AHP) technology and is effective against a broad spectrum of pathogens in only 5 minutes!

VIROXACCDISC4 1 GALLON

ACCEL PREVENTION



Accel PREVENTion Hydrotherapy is designed for use as a cleaner and disinfectant for foot baths and hydrotherapy tubs where the requirement for broad-spectrum germicidal performance in a short contact time is required. Accel Prevention is the FIRST and ONLY concentrated 5 minute intermediate level disinfectant available in Canada.

VIROXACCPREXH4 1 GALLON

ACCEL TB



Accel TB is a Ready-to-use One-step disinfectant cleaner for use as a daily cleaner and as a disinfectant of environmental surfaces and patient care devices. Accel TB uses the Accelerated Hydrogen Peroxide (AHP) technology and is effective against a broad spectrum of pathogens in only 5 minutes!

LIQUID

VIROXACCDISR1-TB
VIROXACCDISR4-TB

1 LT
4 LT

WIPES

VIROXWIP1-TB	6X7"	160/TUB
VIROXWIP1-TB	10X10"	60/TUB
VIROXWIP1TB-SP	7.5X10"	10/PK SOFT PACK

CAVICIDE

Surface Disinfectant/Decontaminant Cleaner

CaviCide is a convenient, ready-to-use, intermediate-level surface disinfectant which is effective against TB, HCV, bacteria (including MRSA and VRE) and fungi. It is safe for all areas of the facility including NICU, operating rooms, isolation rooms, patient care areas and laboratories. When used as directed, it will also effectively clean and decontaminate critical and semi-critical instrumentation.

CaviCide is safe for use on non-porous surfaces. It is perfect for cleaning and decontamination of environmental and medical device surfaces.



CAVIWIPES

Low Alcohol (17.20% isopropanol) Surface Disinfectants

Presaturated with CaviCide®, CaviWipes is a cleaner and disinfectant, in one. CaviWipes are durable, nonwoven, nonabrasive towelettes that offer quick, easy-to-use, time-saving convenience. They're recommended for use on non-porous surfaces and fixtures.

Ideal for operating rooms, surgical centers, neonatal units, medical and dental operatories and other critical-care areas. CaviWipes should not be disposed of in the toilet.

CAVICIDE

METRMX1024	24 OZ SPRAY
METRMX1000	1 GALLON
METRMX1025	2.5 GALLON

CAVIWIPES

METRMX1100	160/TUB
METRMX1150	66/TUB (XL)
METRMX1224	45/PK FLAT PACK

METRICIDE

High-Level Disinfectant and Sterilant

MetriCide is a 2.6% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant on immersible, lensed and delicate instruments for up to 14 days.



METRICIDE 28

High-Level Disinfectant and Sterilant

MetriCide 28 is a long life, 2.5% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant for up to 28 days.

14 DAY

METRMX1400	3.78 LT
------------	---------

28 DAY

METRMX2800	3.78 LT
------------	---------

METRIZYME

Extra Strength Dual-Enzymatic Detergent



MetriZyme is a detergent cleaner comprising of protease dual-enzymes. The active ingredients offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid soils and protein soils. The protein-dissolving enzymes in MetriZyme® allow for fast, simple removal of body fluids that often cannot be accomplished manually with soaps and other detergents.

METRMX4000

1 GALLON

METRISPONGE

Sponge Saturated with Dual-Enzymatic Detergent



According to 2008 SGNA guidelines, manual cleaning of endoscopes immediately after patient use is the first and most important step in removing the microbial burden prior to automated or manual disinfection.

METRMX4025

25/BX

EMPOWER

Dual-Enzymatic Detergent



EmPower is a fresh scented, detergent cleaner containing two protease enzymes that offers the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments. Non-corrosive to metals.

METRMX4100

1 GALLON

EMPOWER FOAM

Foaming Enzymatic Spray



EmPower Foam propels thick foam of protease enzymes. Effectively cleans away blood, tissue, mucous and other protein-rich body fluids. Foaming spray is ideal for use in surgery, central processing, operating room and other departments where instrument cleaning is not immediately available.

METRMX4224

24 OZ

MetriLube™

MetriLube is an oil-based instrument lubricant that coats metal instruments to prevent spotting, staining, and rusting of surgical instruments and helps keep hinges working freely.

MetriLube is offered as a concentrate as well as a 24-oz ready-to-use spray bottle, which is suitable for alternate care facilities.

- Extends life of instruments
- Anticorrosive
- Non-greasy and silicone-free
- Economical

METRMX3425

24oz SPRAY



MILTEX Spray Lube

MILTEX Spray Lube, 8 fluid ounces (0.24 liter) bottles (case of 12). Pre-mixed, non-silicone and anti-corrosive formula in a ready-to-use 8 ounce spray pump bottle. Helps prevent spotting, staining and rusting. Lubricates instrument box locks, screw locks, scissor blades, and any other metal on metal moving part.

MILTEX3-700

8oz



MAGNA DISPOSABLE SCALPELS

Magna Disposable Scalpels are individually packaged. Each scalpel blade is attached to a plastic handle, and comes with a blade cover for user protection.



ALMM92-10	SIZE 10	
ALMM92-11	SIZE 11	
ALMM92-15	SIZE 15	
ALMM92-20	SIZE 20	
ALMM92-21	SIZE 21	
ALMM92-22	SIZE 22	10/bx

STERILE DISPOSABLE SCALPELS

Almedic disposable Scalpels are manufactured from top quality surgical steel, sterilized by gamma radiation. Each scalpel blade is individually packad and supplied with a blade cover for user protection.



ALMA6-510	SIZE 10	
ALMA6-511	SIZE 11	
ALMA6-515	SIZE 15	
ALMA6-521	SIZE 21	
ALMA6-523	SIZE 23	10/bx

SCALPEL BLADES ONLY

Almedic Scalpel Blades are manufactured from the finest surgical quality stainless and carbon steel. Each blade is individually sealed in aluminum foil.



100/bx

SCALPEL HANDLES

Reusable handles also available in stainless steel. Some sizes available in chrome, magna or utility.



NAIL OR BONE RASP



Single ended file
inside 15cm 6"

ALMA4-300



Single ended file
outside 15cm 6"

ALMA4-302



Double ended
1 file outside
1 file inside
17cm 6 3/4"

ALMA4-308

BLACKS FILE



11.5 CM 4.5"

ALMA4-306

PROBE & EXCAVATOR



Double ended
Blunt hook end
Blunt probe end
15cm 6"

ALMA4-315

FILE



Double ended
Angled end
12.5cm 5"

ALMA4-316

FILE



Double sided
Coarse side
Fine side
17cm 6 3/4"

ALMA4-320

FILE



Double sided & double ended
Coarse side & fine side
17cm 6 3/4"

ALMA4-322


DIAMOND DEB FILE



20cm 8"

ALMA4-324

SPATULA



Double ended
18cm 7"

ALMA4-310

SPATULA & PACKER



Oval spatula 7mm wide
Packer 25mm wide
14cm 5 1/2"

ALMA4-311

PACKER



Double ended
Straight end & angled end
40mm wide
16cm 6 1/4"

ALMA4-312

PACKER & EXCAVATOR



Double ended
Blunt hook end
2 x 25mm packer end
14cm 5 1/2"

ALMA4-313

NAIL SPLITTER



Angled
2mm wide
12.5cm 5"

ALMA4-342

INGROWN NAIL SHAVER



Fenestrated with
inside cutting edge
12.5cm 5"

ALMA4-346

VERUCCA CURETTE



Single ended
4.0mm cup
12.5cm 5"

ALMA4-400

CURETTE EXCAVATOR



Single ended
with hole
15cm 6"

ALMA4-410	1.0mm cup
ALMA4-412	1.5mm cup
ALMA4-414	2.0mm cup
ALMA4-416	2.5mm cup

CURETTE EXCAVATOR



Double ended
1.5mm & 2.0mm cup
with hole
14cm 5 1/2"

ALMA4-430

CURETTE EXCAVATOR



Double ended
1.5mm & 2.5mm cup
with hole 14cm 5 1/2"

ALMA4-432

LISTER BANDAGE SCISSORS



AMG570308
5.5" S/STEEL

UNIVERSAL SCISSORS




This paramedic scissor features hardened stainless steel blades and autoclavable plastic handles.

AMG120214
ECONOMY 7.5"

AMG120210
UNIVERSAL 8"

IRIS SCISSORS STRAIGHT



Fine & Stainless Steel

AMG570200 3.5"
AMG570204 4.5"

IRIS SCISSORS CURVED



Fine & Stainless Steel


AMG570208 3.5"
AMG570212 4.5"

STERILE SCISSORS



BOW-119-316A
4.5" Sharp/Blunt

KELLY FORCECP



AMG570848
5.5" Straight

AMG570856
5.5" Curved

HALSTEAD MOSQUITO FORCEPS



Halstead Mosquito
Forceps, 5 in., Curved,
stainless steel

AMG570828



Halstead Mosquito
Forceps, 5 in., Straight,
stainless steel

AMG570824

DRESSING FORCEPS



YMK-279-5 5"
 YMK-279M 4"

SPLINTER FORCEP FINE



Splinter Forceps, Stainless Steel,
 Fine, Straight

AMG570584 4"
 AMG570592 4.5"
 ALMA10-604 4.5"
 ALMA10-606 5"

ADSON DRESSING FORCEPS



Adson Dressing
 Forceps, 4 1/2 inch
 (11.4 cm), Serrated,
 Stainless Steel

AMG570630



Adson Dressing Forceps,
 4 1/2 inch (11.4 cm), 1 x
 2 teeth, Stainless Steel

AMG570634


CUTICLE NIPPER



10cm 4"

ALMA4-211 4mm Blade
 ALMA4-212 7mm Blade
 ALMA4-213 9mm Blade

TISSUE NIPPER



Double Spring

ALMA4-218 10cm 4"
 ALMA4-219 11.5cm 4 1/2"

NAIL NIPPER



Barrel Spring

ALMA4-224
ALMA4-225

13cm 5 1/4"
11.5CM 4 1/2"

NAIL SPLITTER



Single Spring

ALMA4-229
ALMA4-230

11.5cm 4 1/2"
12.5cm 5"

NAIL NIPPER



Front cutting
14cm 5 1/2"

ALMA4-232

INGROWN NAIL SPLITTING SCISSORS



One Serrated Blade
15cm 6"

ALMA4-234

TOENAIL CLIPPERS



With File
7.5cm 3"

ALMA4-238

FINGERNAIL CLIPPERS



With File
5cm 2"

ALMA4-242

NAIL SCISSORS



Curved
9cm 3 1/2"

ALMA4-204

SPLINTER FORCEP FINE



Curved
9cm 3 1/2"

ALMA4-208

EMERY BOARD



Emery boards, 2 sided,
coarse/fine

AMG018860 144/BX

MANICURE STICK



- 100% White Birch Wood
- Double-Sided
- Pointed & Angled Ends

AMG018830 144/BX

NYLON INSTRUMENT BRUSH



Plastic Handle Small Utility Brushes
are used for cleaning and scrubbing
small surface areas.

JUST1434

NYLON SCRUB BRUSH



AMG118140

- Open back hand grip
- 16 x 7 rows of 1/2" (1.3 cm) bristles

Scrub brushes are made from polypropylene and can withstand temperatures of up to 80°C or 176°F.

INSTRUMENT TRAYS & FLAT/SOLID COVERS

The wide assortment of heavy duty, seamless trays have covered corners for easy sterilization. The covers have recessed strap handles for ease in handling and storing.

ITEM CODE	DIMENSIONS: L x W x H IN. (CM)
POLAR43	5 ³ / ₈ x 4 x 2 ⁵ / ₈ (13.7 x 10.2 x 6.7)
POLAR43-2	Cover—Fits 43 Tray
POLAR63	6 ¹ / ₈ x 6 ¹ / ₈ x 2 ⁷ / ₈ (15.6 x 15.6 x 7.3)
POLAR63-2	Cover—Fits 63 Tray
POLAR952	8 ⁷ / ₈ x 5 x 2 (22.5 x 12.7 x 5.1)
POLAR953	8 ⁷ / ₈ x 5 x 3 (22.5 x 12.7 x 7.6)
POLAR952-2	Cover—Fits 952 & 953 Trays
POLAR1002	10 x 6 ¹ / ₂ x 2 (25.4 x 16.5 x 5.1)
POLAR1002-2	Cover—Fits 1002 Tray
POLAR1202	12 ¹ / ₄ x 7 ³ / ₄ x 2 ¹ / ₄ (31.1 x 19.7 x 5.7)
POLAR1202-2	Cover—Fits 1202 Tray
POLARE1082	10 x 8 ¹ / ₄ x 2 ¹ / ₂ (25.4 x 21.0 x 6.4)
POLAR1082-2	Cover—Fits 1082 Tray
POLARE1650	16 ¹ / ₂ x 10 x 2 ¹ / ₂ (41.9 x 25.4 x 6.4)
POLARE1654	16 ¹ / ₂ x 10 x 4 (41.9 x 25.4 x 10.2)
POLAR1650-2	Cover—Fits Above Two Trays
POLARE902	6 ⁷ / ₈ x 4 ¹ / ₄ x 2 (17.5 x 10.8 x 5.1)
POLARE904	6 ⁷ / ₈ x 4 ¹ / ₄ x 4 (17.5 x 10.8 x 10.2)
POLAR902-2	Cover—Fits Above Two Trays
POLARE2302	13 ⁷ / ₈ x 12 ⁷ / ₈ x 2 ¹ / ₂ (35.2 x 32.7 x 6.4)
POLARE2304	13 ⁷ / ₈ x 12 ⁷ / ₈ x 4 (35.2 x 32.7 x 10.2)
POLARE2306	13 ⁷ / ₈ x 12 ⁷ / ₈ x 6 (35.2 x 32.7 x 15.2)
POLAR2302-2	Cover—Fits Above Three Trays
Standard Weight (22 Gauge)	
POLARE06062	6 ⁷ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂ (17.5 x 16.2 x 6.4)
POLARE06064	6 ⁷ / ₈ x 6 ³ / ₈ x 4 (17.5 x 16.2 x 10.2)
POLARE06066	6 ⁷ / ₈ x 6 ³ / ₈ x 6 (17.5 x 16.2 x 15.2)
Heavy Weight (20 Gauge)	
POLARS06062	6 ⁷ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂ (17.5 x 16.2 x 6.4)
POLARS06064	6 ⁷ / ₈ x 6 ³ / ₈ x 4 (17.5 x 16.2 x 10.2)
POLARS06066	6 ⁷ / ₈ x 6 ³ / ₈ x 6 (17.5 x 16.2 x 15.2)
POLAR602-2	Cover—Fits Above Six Trays
Standard Weight (22 Gauge)	
POLARE10062	10 ³ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂ (26.4 x 16.2 x 6.4)
POLARE10064	10 ³ / ₈ x 6 ³ / ₈ x 4 (26.4 x 16.2 x 10.2)
POLARE10066	10 ³ / ₈ x 6 ³ / ₈ x 6 (26.4 x 16.2 x 15.2)
Heavy Weight (20 Gauge)	
POLARS10062	10 ³ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂ (26.4 x 16.2 x 6.4)
POLARS10064	10 ³ / ₈ x 6 ³ / ₈ x 4 (26.4 x 16.2 x 10.2)
POLARS10066	10 ³ / ₈ x 6 ³ / ₈ x 6 (26.4 x 16.2 x 15.2)
POLAR402-2	Cover—Fits Above Six Trays

Order pans and covers separately.



INSTRUMENT TRAYS & FLAT/SOLID COVERS

ITEM CODE DIMENSIONS: L x W x H
IN. (CM)

Standard Weight (22 Gauge)

POLARE12062 12³/₄ x 6⁷/₈ x 2¹/₂ (32.4 x 17.5 x 6.4)
POLARE12064 12³/₄ x 6⁷/₈ x 4 (32.4 x 17.5 x 10.2)
POLARE12066 12³/₄ x 6⁷/₈ x 6 (32.4 x 17.5 x 15.2)

Heavy Weight (20 Gauge)

POLARS12062 12³/₄ x 6⁷/₈ x 2¹/₂ (32.4 x 17.5 x 6.4)
POLARS12064 12³/₄ x 6⁷/₈ x 4 (32.4 x 17.5 x 10.2)
POLARS12066 12³/₄ x 6⁷/₈ x 6 (32.4 x 17.5 x 15.2)
POLAR302-2 Cover—Fits Above Six Trays

Standard Weight (22 Gauge)

POLARE12102 12³/₄ x 10³/₈ x 2¹/₂ (32.4 x 26.4 x 6.4)
POLARE12104 12³/₄ x 10³/₈ x 4 (32.4 x 26.4 x 10.2)
POLARE12106 12³/₄ x 10³/₈ x 6 (32.4 x 26.4 x 15.2)

Heavy Weight (20 Gauge)

POLARS12102 12³/₄ x 10³/₈ x 2¹/₂ (32.4 x 26.4 x 6.4)
POLARS12104 12³/₄ x 10³/₈ x 4 (32.4 x 26.4 x 10.2)
POLARS12106 12³/₄ x 10³/₈ x 6 (32.4 x 26.4 x 15.2)
POLAR1311-2 Cover—Fits Above Six Trays

Standard Weight (22 Gauge)

POLARE20122 20³/₄ x 12³/₄ x 2¹/₂ (52.7 x 32.4 x 6.4)
POLARE20124 20³/₄ x 12³/₄ x 4 (52.7 x 32.4 x 10.2)
POLARE20126 20³/₄ x 12³/₄ x 6 (52.7 x 32.4 x 15.2)

Heavy Weight (20 Gauge)

POLARS20122 20³/₄ x 12³/₄ x 2¹/₂ (52.7 x 32.4 x 6.4)
POLARS20124 20³/₄ x 12³/₄ x 4 (52.7 x 32.4 x 10.2)
POLARS20126 20³/₄ x 12³/₄ x 6 (52.7 x 32.4 x 15.2)
POLAR2000-2 Cover—Fits Above Six Trays

Order pans and covers separately.



STERILIZATION TRAYS

HOT/COLD

Tray units can be used for either hot or cold sterilization. The covers feature recessed strap handle for ease and safety in transporting, handling, and storing.



ITEM CODE CAPACITY DIMENSIONS: L x W x H
QT. (LITRE) IN. (CM)

POLARE20062 3¹/₂ (3.31) 20³/₄ x 6³/₈ x 2¹/₂ (52.7 x 16.2 x 6.4)
POLARE20064 5¹/₄ (4.97) 20³/₄ x 6³/₈ x 4 (52.7 x 16.2 x 10.2)
POLARE20066 8³/₄ (8.28) 20³/₄ x 6³/₈ x 6 (52.7 x 16.2 x 15.2)
POLARS20062 3¹/₂ (3.31) 20³/₄ x 6³/₈ x 2¹/₂ (52.7 x 16.2 x 6.4)
POLARS20064 5¹/₄ (4.97) 20³/₄ x 6³/₈ x 4 (52.7 x 16.2 x 10.2)
POLARS20066 8³/₄ (8.28) 20³/₄ x 6³/₈ x 6 (52.7 x 16.2 x 15.2)
POLAR2006-2 Flat, Solid Cover—Fits all above pans

Order pans and covers separately.

PERFORATED INSTRUMENT TRAYS

Constructed of heavy duty stainless steel, PMP series perforated instrument trays are finished to provide a tough, clean and corrosion-resistant surface for long term use in a variety of applications.



ITEM CODE DIMENSIONS: L X W H
IN. (CM)

Perforated Mayo Trays

POLAR10F-PMT	10 x 6½ x ¾ (25.4 x 16.5 x 1.9)
POLAR13F-PMT	13 ¹¹ / ₁₆ x 9 ¹³ / ₁₆ x ¾ (34.8 x 24.9 x 1.9)
POLAR15F-PMT	15½ x 10 ⁵ / ₈ x ¾ (38.4 x 27.0 x 1.9)
POLAR17F-PMT	17 x 11 ⁵ / ₈ x ¾ (43.2 x 29.5 x 1.9)
POLAR19F-PMT	19½ x 12 ⁵ / ₈ x ¾ (48.6 x 32.1 x 1.9)
POLAR1520-PMT	21¼ x 16¼ x 7/8 (54.0 x 41.3 x 2.2)

For excellent steam sterilization without moisture build-up.

High-Sided Perforated Trays

POLAR952-HPT	8 ⁷ / ₈ x 5 x 2 (22.5 x 12.7 x 5.1)
POLAR952-FPP*	8 ⁷ / ₈ x 5 x 2 (22.5 x 12.7 x 5.1)
POLAR1002-HPT	10 x 6½ x 2 (25.4 x 16.5 x 5.1)
POLAR1002-FPP*	10 x 6½ x 2¼ (25.4 x 16.5 x 5.7)
POLAR1202-HPT	12¼ x 7¾ x 2¼ (31.1 x 19.7 x 5.7)
POLARS12102-HPT	12¾ x 10 ³ / ₈ x 2½ (32.4 x 26.4 x 6.4)
POLARE1650-HPT	16½ x 10 x 2½ (41.9 x 25.4 x 6.4)
POLARS20122-HPT	20¾ x 12¾ x 2½ (52.7 x 32.4 x 6.4)

High sides protect fragile instruments during set-up and procedures.

*With end handles



OBLONG INSTRUMENT/DRYING TRAYS

The user of these trays will benefit from coved corners, and large, open bead, which allow for easy handling, stacking and cleaning.



ITEM CODE DESCRIPTION DIMENSIONS: L X W H
IN. (CM)

POLAR10F	With Rolled Bead	10 x 6½ x ¾ (25.4 x 16.5 x 1.9)
POLAR13F	With Rolled Bead	13 ¹¹ / ₁₆ x 9 ¹³ / ₁₆ x ¾ (34.8 x 24.9 x 1.9)
POLAR15F	With Rolled Bead	15½ x 10 ⁵ / ₈ x ¾ (38.4 x 27.0 x 1.9)
POLAR17F	With Rolled Bead	17 x 11 ⁵ / ₈ x ¾ (43.2 x 29.5 x 1.9)
POLAR19F	With Rolled Bead	19½ x 12 ⁵ / ₈ x ¾ (48.6 x 32.1 x 1.9)
POLAR1520	With Rolled Bead	21¼ x 16¼ x 7/8 (54.0 x 41.3 x 2.2)
POLAR1525	With Flat Flange	25 x 15 x 5/8 (63.5 x 38.1 x 1.6)
POLAR25F	With Flat Flange	25 x 20 x 5/8 (63.5 x 50.8 x 1.6)
POLARE12101		12¾ x 10 ³ / ₈ x 1¼ (32.4 x 26.4 x 3.2)
POLARE20121		20¾ x 12¾ x 1¼ (52.7 x 32.4 x 3.2)



PERFORATED INSTRUMENT TRAYS

1/4" (0.64 cm) diameter holes in tray bottom and sides enhance sterilizing liquid flow, reducing wet packs. Tray edges are rounded to aid in handling.



ITEM CODE	DIMENSIONS: L X W H IN. (CM)
POLARP20121	20 ³ / ₄ x 12 ³ / ₄ x 1 ¹ / ₄ (52.7 x 32.4 x 3.2)
POLARP20122	20 ³ / ₄ x 12 ³ / ₄ x 2 ¹ / ₂ (52.7 x 32.4 x 6.4)
POLARP20124	20 ³ / ₄ x 12 ³ / ₄ x 4 (52.7 x 32.4 x 10.2)
POLARP20126	20 ³ / ₄ x 12 ³ / ₄ x 6 (52.7 x 32.4 x 15.2)
POLARP12101	12 ³ / ₄ x 10 ³ / ₈ x 1 ¹ / ₄ (32.4 x 26.4 x 3.2)
POLARP12102	12 ³ / ₄ x 10 ³ / ₈ x 2 ¹ / ₂ (32.4 x 26.4 x 6.4)
POLARP12104	12 ³ / ₄ x 10 ³ / ₈ x 4 (32.4 x 26.4 x 10.2)
POLARP12106	12 ³ / ₄ x 10 ³ / ₈ x 6 (32.4 x 26.4 x 15.2)

INSTRUMENT TRAYS WITHOUT COVER & VERY SLIGHT FLANGE



ITEM CODE	DIMENSIONS: L X W H IN. (CM)
POLAR4V	8 ³ / ₄ x 4 ⁵ / ₈ x 3 ¹ / ₈ (22.2 x 11.7 x 7.9)
POLAR6V	10 ⁵ / ₁₆ x 5 ¹ / ₂ x 4 (26.2 x 14.0 x 10.2)


SAFE-SEAL DUET STERILIZATION POUCHES

Assured sterilization with a Class 4 internal multi-variable indicator and a Class 1 external process indicator for steam sterilization. Meets CDC guidelines for effective sterilization and responds to all variables: time, temperature and saturated steam.

Secure seal with extra wide chevron seal protects against instrument breaching while high quality porous medical-grade paper meets medical standard.

Perforated fold facilitates easy and accurate folding of adhesive strip to form airtight seal. Thumb notch enables easy opening of pouch.

Lead free, high-quality inks compatible with a large number of commonly used disinfectants and assure easy confirmation of sterilization process.



MEDICOM88005	3 1/2" x 5 1/4"	200/PK
MEDICOM88015	2 1/4" X 4"	200/PK
MEDICOM88010	3 1/2" x 9"	200/PK
MEDICOM88025	5 1/4" x 10"	200/PK

STERILIZATION REELS

Various sized reels allows for customized sterilization of a wide variety of items.

Individual and flexible formats save time and lower costs. Ideal for specialty items of various sizes.

Special grade, heavy-weight paper provides superior strength. Blue tinted transparent material for easy viewing and possible instrument breaching.



MEDICOM9520	3" x 100'
MEDICOM9530	4" x 100'
MEDICOM9540	6" x 100'

COMPLY™ LEAD FREE STEAM INDICATOR TAPE

As an exposure monitor, 3M™ Comply Indicator tapes securely seal packs and enables CSSD and OR personnel to know at a glance that the packs have been exposed to the sterilization process. The chemical indicator ink, printed in a diagonal striped pattern, changes color or turns darker when exposed to the sterilization process.

Benefits

Lead Free Ink means no hazardous waste issues

Trusted, high quality adhesive

stretchable backing that minimizes tape "pop-off" during sterilization

can be written on or labeled with pre-printed labels



3M1322-18MM 18MM x 55M

3M1322-24MM 24MM x 55M

COMPLY™ INDICATOR STRIP-STEAM

The 3M™ Comply™ Indicator Strip enables pack control by validating that the sterilant has penetrated to the point of placement in the pack and confirming that sufficient exposure conditions have been met. Use in steam sterilization cycles.



3M1250 240/BX

COMPLY™ (STERIGAGE™) CHEMICAL INTEGRATOR

As an internal pack monitor, the 3M™ Comply™ (SteriGage™) Class 5 Steam Chemical Integrators offer immediate, accurate and easy to read method of monitoring the sterilization process conditions inside each pack. It can be used in all steam sterilizers. The easy-to-read moving front style has been shown to have a higher reading accuracy than color interpretation style indicators.

Benefits

Moving front ink technology provides easy to read 'Accept or Reject' result

Can be used in all 250 - 275°F steam sterilization cycles

Response correlates with a Biological Indicator at 3 time/temperature relationships under ideal steam sterilization conditions

Confirmed by BSI, a leading global independent Product Testing Service Company, to meet the performance requirements of ANSI/AAMI/ISO 11140-1:2005 Class 5.



3M1243B 100/BG

CSR WRAP

CSR Autoclave Wrap

ARJ052-1212 12" x 12"



DISTILLED WATER

Ripple Effects distilled water.

CONT6375 3.85L



3M™ ATTEST™ BIOLOGICAL INDICATOR INCUBATOR 116

Biological Indicator Incubator, 14-vial (round); 56°C steam, 120 nominal voltage required; 90-132V acceptable.

Attest incubators provide optimal incubation conditions for the accurate readout of the 3M™ Attest™ Biological Indicators. Since the growth of each test organism is temperature dependent, it is critical that the appropriate Attest Biological Indicator be used with the correct Attest incubator.

3M116 1 EA



3M™ ATTEST™ INDICATOR 1262P STEAM

This 3M™ Attest™ Indicator is for 132 degrees C vacuum assisted or 121 degrees C gravity steam sterilizers. It provides assurance of sterility by challenging the sterilizer with quantifiable, highly resistant spores. Incubation is for 48 hours.

3M1262P 25 EA



UTILITY BRUSH

Brush Soft Brass Utility *autoclavable

JUST1434



PRESTIGE CLASSIC 2100 ELECTRIC AUTOCLAVE



The Prestige Classic 2100 is one of the world's largest selling and cost effective automatic autoclaves. Sold in 80 countries, it has gained an enviable reputation for reliability, durability and quality. As a commitment to this quality, each unit comes with a Lifetime Guarantee on its heating element and is issued with a test certificate showing the temperature time cycle.

Additional Features:

- Mechanical interlock system prevents lid from being removed while pressure remains in the chamber
- Gasket offset device prevents pressure build up if lid is not correctly fitted
- Net weight of 10 lbs and chamber size of approximately 8" x 9" makes the 2100 ideal for most sterilization applications
- Simple 1-step operation starts the fully automatic cycle
- Light sequence indicators show "power to the unit", "cycle in progress", "sterilization cycle complete"
- Electronic detectors turn off power if there is insufficient water charge
- Application Specific Integrated Circuit (ASIC) is used to improve performance and reliability

PACH03-21006

EZ9 - AUTOMATIC AUTOCLAVE



With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer EZ9 fully automatic autoclave will fill, sterilize, exhaust and dry at the touch of a button.

EZ9 Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying.
- Hepa filter 0.2 µm air filter (replaceable).
- Pre-loaded cycles for: unwrapped instruments, wrapped/packs, glassware and additional drying cycle.
- All program parameters can be customized and stored.
- Automatic shut off at the end of both the sterilization and dry cycles.
- Audible and visual cycle interruption alert.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir.
- Constant monitoring of temperature and pressure for sterilization.
- Power outage recovery system.
- Door design prevents steam from coming in contact with control panel.
- Optional internal printer documents: date, time, temperature and pressure.
- 510k clearance for dental handpieces.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- International certifications.
- 2 year parts and labor warranty.

TUTTEZ9P WITH PRINTER

Edan SonoTrax Vascular Doppler

Edan SonoTrax Vascular Doppler 8MHz or 4MHz probe for arterial and venous blood flow. The hand-held Doppler has been widely adopted as standard equipment for the assessment of peripheral vascular disease. The versatile, inexpensive and portable equipment allows the fast and simple but effective assessment of both arterial and venous disease. Besides, it helps physicians to diagnose some diseases, such as lower limb and foot ulceration in both diabetics and non-diabetics alike. SonoTrax Vascular offers a traditional pencil-style 4MHz or 8MHz probe. It can locate both superficial and deep vessels faster and provides stronger signals.

- Ergonomic design, compact and light
- Easy-to-use
- Accurate detection with clear sound
- High sensitivity interchangeable probe
- Support 4MHZ,5mhz ,8MHZ waterproof probes
- Long time continuous use (about 8h)
- Automatic power-off

SOLIC-MS314322 Vascular doppler with 4 MHz probe

SOLICMS314321 Vascular doppler with 8 MHz probe

SOLIC-120114346 SonoTrax probe only 4 MHz

SOLIC-120114347 SonoTrax probe only 8 MHz

Optional Accessories

Part Number	Description
12.01.14320	2MHz Waterproof Interchangeable Probe
12.01.14321	3MHz Waterproof Interchangeable Probe
12.01.14346	4MHz Waterproof Interchangeable Probe
02.01.104822	5MHz Waterproof Interchangeable Probe (Only for SonoTrax Vascular)
12.01.14347	8MHz Waterproof Interchangeable Probe
01.21.064122	9V Rechargeable Ni-MH Battery (200mAh)
01.21.064113	9V Ni-MH Battery Charger (USA Standard, ETL Mark)
11.21.64048	9V Alkaline Battery
03.28.328007	MT-503 Rolling stand for SonoTrax Series
02.04.240040	Mounting Plate Assembly for SonoTrax series
02.04.101976	Basket for Rolling Stand (in the bottom)
01.56.14268-01	Blue Carry Case
01.56.104581-12	Normal Carry case
11.13.36032	Audio Record Line-in Cable
01.14.104751	Earphone



DARCO ULTRA OFF-LOADING SHOE

The Ultra™ Off-loading Shoe is the most versatile post operative, trauma and off-loading boot available anywhere. The soft ankle strap eliminates sharp edges and allows full shoe on foot contact through gait. The PegAssist™ off-loading insole comes standard, and the boot features a roomy closed-toe box design to accommodate large dressings and keeps toes warm and dry.

- Rocker sole reduces plantar pressure under the forefoot and heel by over 25%
- PegAssist™ insole comes standard and provides an additional 60% reduction in focal pressure when off-loading wounds
- Enhanced toe box is extra roomy to accommodate very large dressings while still keeping the toes safe from bumps and the elements
- Removable / reversible forefoot closure provides security without the pressure of buckles or straps
- Available in black only



DARCO ALL PURPOSE BOOT

The APB™ All Purpose Boot's all season closed toe protection allows dressings to be kept clean and toes dry, providing more protection for post operative procedures and fiberglass casts

- Square toe design acts as a bumper when K-wires are present and provides additional room and comfort for the patient
- Removable EVA insole is fully customizable to off-load pressure points if needed
- Wider sole and deep opening gives more room to accommodate large dressings, bandages and casts
- Wide ankle strap eliminates heel slippage as well as seats and secures the foot or cast in the shoe
- Rocker bottom encourages smooth roll through gait pattern and reduces plantar pressure on the heel and forefoot by over 25%
- MetaShank™ construction provides for a reinforced sole under the Metatarsal heads to reduce motion
- Strapless forefoot closure provides uniform security without the pressure of buckles or straps



DARCO ULTRA OFF-LOADING SHOE

The Ultra™ Off-loading Shoe is the most versatile post operative, trauma and off-loading boot available anywhere. The soft ankle strap eliminates sharp edges and allows full shoe on foot contact through gait. The PegAssist™ off-loading insole comes standard, and the boot features a roomy closed-toe box design to accommodate large dressings and keeps toes warm and dry.

- Rocker sole reduces plantar pressure under the forefoot and heel by over 25%
- PegAssist™ insole comes standard and provides an additional 60% reduction in focal pressure when off-loading wounds
- Enhanced toe box is extra roomy to accommodate very large dressings while still keeping the toes safe from bumps and the elements
- Removable / reversible forefoot closure provides security without the pressure of buckles or straps
- Available in black only



DARCO ALL PURPOSE BOOT

The APB™ All Purpose Boot's all season closed toe protection allows dressings to be kept clean and toes dry, providing more protection for post operative procedures and fiberglass casts

- Square toe design acts as a bumper when K-wires are present and provides additional room and comfort for the patient
- Removable EVA insole is fully customizable to off-load pressure points if needed
- Wider sole and deep opening gives more room to accommodate large dressings, bandages and casts
- Wide ankle strap eliminates heel slippage as well as seats and secures the foot or cast in the shoe
- Rocker bottom encourages smooth roll through gait pattern and reduces plantar pressure on the heel and forefoot by over 25%
- MetaShank™ construction provides for a reinforced sole under the Metatarsal heads to reduce motion
- Strapless forefoot closure provides uniform security without the pressure of buckles or straps



SWEDE-O NIGHT SPLINT



Allows the patient to apply a consistent 90° angle stretch to the plantar fascia by holding the foot in a gentle dorsiflexion. Soft gray padded cover provides greater patient comfort and protection against irritation to the skin and breakdown

NIGHT SPLINT WITH TREAD



- The featherweight design offers superior comfort while being worn when sleeping
- Slip resistant tread offers maximum patient safety and enable them easier mobility in the house
- Adjustable strap improves patient's ability to apply steady stretch to aid in pain relief from achilles tendonitis and plantar fasciitis
- Plantar flexion contracture greater than 10° degrees
- Grey shell can be wiped clean with a damp cloth using a very mild soap
- Straps and blue liner can be hand washed in lukewarm water with a mild soap, thoroughly rinsed and laid flat to dry
- Easy release clips for freedom of use
- Fits either foot
- Three sizing options

BORT HALLUX VALGUS SUPPORT



- The BORT Valco® Hallux Valgus splint, the comfortable solution for hallux valgus
- Manufactured in an virtually unbreakable ABS plastic material, is a conservative therapy in the hallux valgus treatment. It helps to corrects the unnatural X-position of the big toe by forcing a more straightened position. Also providing optimal post-operative capsular relief following corrective surgery
- The anatomically adapted shape of the splint aids the physiological correction of the basic joint of the big toe.
- It affords much comfort thanks to the generously dimensioned foam padding, which safely protects the critical area of the toe and in the ball of the foot
- To avoid excessive heat, the splint is provided with a ventilation hole in the toe area.
- Two adjustable Velcro fasteners allow precise adjustment of the moments of leverage with the splint remaining in the same optimum position

WAFFLE® FOOT CUSHION - LINED

The Foot WAFFLE® Custom surrounds the lower limb in air and is designed to float the heel while sustaining the integrity of the foot, ankle, calf and Achilles area. The air-filled straps and foot pad offer pressure ulcer prevention, treatment and patient comfort. It can be used conveniently at home and on O.R. tables, stretchers, ICU beds, wheelchairs and general ward beds.

*36 Month Warranty



EHOB503FW	SMALL	CALF 10" TO 11.5"
EHOB504FW	MEDIUM	CALF 12" TO 14"
EHOB505FW	LARGE	CALF 14.5" TO 18"

WAFFLE® HEEL ELEVATOR

The Waffle Heel Elevator is used for prevention and healing through Stage IV pressure ulcers. It suspends the heel completely off the surface and reduces plantar flexion and foot drop. The cloth lining and air venting holes keep the patient cool and dry.



EHOB345CX ONE SIZE FITS MOST PATIENTS

WAFFLE FOOTHOLD®

A new positioning device for the lower extremity that promotes proper alignment of the ankle, knee and hip, to prevent and treat heel ulcers. The patented Hydro-Cushion floats and lifts the calf while unloading the heel when recumbent. Ideal for use in hospital, long term care and home care.

*30 Day Warranty

EHOBFH10-SMXX	SMALL
EHOBFH10-MDXX	MEDIUM
EHOBFH10-LGXX	LARGE
EHOBFH10-XLXX	EXTRA LARGE



HEEL FLOAT



- Floats the heel over an air cavity for zero pressure
- Provides air circulation to aid healing of existing ulcers
- Foam cushioned sides protect lateral and medial malleoli
- Water-based gel insert lowers pressure beneath ankle
- Cozy Cloth cover reduces skin shear and wicks away moisture
- Anti-foot-drop construction controls dorsal flexion
- Launderable, reusable, and guaranteed for six months

SKIL503034

SMALL 3" WIDE

SKIL503035

MEDIUM 4" WIDE

SKIL503036

LARGE 5" WIDE

ULTRA-SOFT HEEL CUSHION



- Filled with fluffy spun fiber
- Covered with low-friction brushed polyester fabric to reduce skin shear
- Cover wicks away moisture to eliminate skin-damaging perspiration
- Comfortable, extra-wide cloth strap won't irritate skin
- Velcro® closure keeps heel cushion in place
- Reusable and launderable
- Meets California 117B flammability standards

SKIL503030

UNIVERSAL SIZE 1 PR

ALL GEL BUNION PROTECTOR

- Ideal to protect the hallux bunion from friction and pressure
- Soft flexible medical grade gel cushions and comforts tender sensitive skin and surrounding tissue
- Comfortable toe loop, with ample stretch to prevent breakage, slides over 1st digit to hold the product in place
- Thin profile design offers comfort in any footwear



ALL GEL METATARSAL PAD

- 100% medical grade mineral oil based visco-elastic gel effectively absorbs shock & shear forces
- For protection and comfort of the entire forefoot
- Provides even pressure distribution of metatarsal head
- Flexible and stretchy toe loop fits easily over toe to stay in place
- Will not flatten out like foam or felt products
- Dermatologist tested - hypoallergenic
- Washable and reusable



GEL LINED TUBING

- Fully coated with medical grade mineral oil based visco-elastic gel
- All-round digital protection
- Durable and effective
- Uses less room than regular foam tubing and the gel will not bottom out like most foam materials
- Can help to reduce scar tissue
- 6" long tubes - May be cut to required size without fraying
- Widely used in sports for pressure and friction relief
- Can be used to protect over/under-lapping toes
- Available in two materials:
 - Stretchable ribbed fabric material - offers a mild compression depending on the size chosen
 - Mesh – a more elasticized woven material. Suited for athletes
- Oil smooths and enhances skin
- Dermatologist tested
- Washable and reusable



TOE CRESTS



- 100% Medical Grade Silicone absorbs shock and vibration and provides even pressure distribution
- Helps to reduce pressure and pain caused by hammer, claw or mallet toes
- Will not bottom out like foam products

TOE SEPARATORS



- Anatomically molded to fit comfortably between toes reducing shearing and inter-digital pressure
- May be trimmed if required
- Cost effective alternative to silicone compounds
- Also available in 100/bag (special order)

TOE SPREADERS



- Thicker than the Separators
- The polymer's soft conforming properties makes these spreaders extremely comfortable
- May be used as an inter-digital wedge for bunion deformity, overlapping toes or toe drift
- Can be trimmed with scissors
- Also available in 100/bag (special order)

CENTRE CUSHION HEEL CUP



- Dual Density Orthotic - Maximum support and shock absorption using two densities of medical grade silicone
- Firm density cup with softer centre 'shock spot' offering relief to high impact area
- Non-compressible medical grade silicone, reduces shock loads to feet, knees and back by providing even pressure distribution throughout the orthotic
- Hypoallergenic, odorless and fully washable

ADHESIVE FELT

- 1/8 in Thickness
- White
- 5 1/2 in x 2 1/2 Yard Roll

SCHOLL58924 1RL/CS



- 1/4 in Thickness
- White
- 5 1/2 in x 2 1/2 Yard Roll

SCHOLL68926 1RL/CS

KUROTUX HEAVY MOLESKIN

- Beige
- 12 in x 5 Yard Roll
- Water Repellent
- Adhesive Backed

SCHOLL58364 1RL/CS



ORTHOPEDIC FELT

- 1/4 in Thickness
- 21 in x 36 in Roll

SCHOLL58472 1RL/CS

- 1/2 in Thickness
- 21 in x 36 in Roll

SCHOLL58364 1RL/CS





Support Therapy Features

Provides relief from a range of common symptoms:

- Tired, aching legs
- Swollen feet and ankles
- Pregnancy
- Wound care
- Lymphedema
- Travelling
- Varicose and spider veins
- Infectious disease
- Sedentary lifestyle
- Surgery or trauma



Styles Include:

- Athletic
 - Casual Cotton
 - Dress Socks (for Men & Women)
 - Relief
 - Sheer & Ultra Sheer
 - Fashionable dress
- Also available in a variety of colours and some styles offer the options of open toe or closed toe.

Available in a variety of compressions.
Bowers also does custom measurements and fittings
Some styles available in *Latex Free* materials

ACCESSORIES



Rubber Gloves

Practical and safe way to apply and remove compression stockings



Footslip

Applicator for open toe compression stockings



Stocking donner

To aid the process of putting on compressions stockings



Doff N' Donner

The revolutionary donning device for all compression socks and stockings

****Please note that all stockings require measurements and a fitting, which are done in the morning before 11 am, to ensure proper sizing.**

SIGVARIS PERFORMANCE SOCKS

A growing trend among competitive athletes and weekend warriors is the use of graduated compression socks to help improve their athletic performance and recovery. As the global leader in graduated compression therapy, SIGVARIS creates truly superior sports products featuring true graduated compression for athletes and active people. Athletes who wear SIGVARIS Sports Socks and Sleeves will benefit from:

- Increased circulation and improved blood flow
- Increased oxygen to muscle tissue
- Reduced pulled muscles and less exercised-induced muscle soreness in the legs
- Less lactic acid buildup in the legs during exercise
- Improved energy, performance and recovery
- Injury prevention and reduced shin splints & calf cramps

Peak Your Performance

- Improved blood circulation and oxygen delivery
- Less aching & cramping
- Thermal & odor control
- Achilles tendon protector reduces vibration on ligaments



Available Styles & Colours

Men & Women
Black or White
Closed Toe
20-30mmhg

JOBST ACTIVEWEAR

MODERATE - HIGH COMPRESSION

JOBST® ActiveWear is effective leg therapy in an energizing athletic sock, ideal for various activities, from walking to playing sports.

- The ideal combination of therapeutic efficacy and Dri-release® yarn for superior moisture control
- Fits an active lifestyle with a seamless toe and 360° cushioned foot for increased wearing comfort
- All-day comfort knee band that keeps socks up without binding or pinching
- Reinforced heel resists abrasion and provides greater durability

Available Styles & Colours

Men & Women
Black or White
Closed Toe
20-30mmhg & 30-40mmhg



LOW COMPRESSION

- A high performance sock that's durable enough for the most active lifestyles
- Lightweight construction
- Knit-in heel and ample toe are ensure extra comfort and a perfect fit
- Nylon/acrylic blend helps wick away moisture
- Sporty design looks like a regular athletic sock

Available Styles & Colours

Men & Women
White
Closed Toe
15-20mmhg

COMPRESSION BANDAGES

PROFORE & PROFORE LITE

The multi-layer compression bandage system has been specifically designed for the management and treatment of venous leg ulcers and associated conditions. The system pack can be used on patients with ankle circumferences of greater than 18cm (padded), but will give best results on ankle circumferences of 18-25cm.

SMTH66000016 PROFORE 4 LAYER
SMTH66000415 PROFORE LITE 3 LAYER



SUREPRESS COMPRESSION BANDAGE

SurePress® High Compression Bandage is an elastic wrap that provides sustained graduated compression in an easy to use format. This reusable bandage has a yellow center line and two rectangular extension indicators for easier application and appropriate amount of compression for a variety of leg sizes. It is indicated for the management of venous leg ulcers and associated conditions. It can be used alone or in combination with SurePress® Absorbent Padding.



CONV650947

SUREPRESS PADDING

CONV650948



COMPRILAN

Comprilan is made of 100% cotton and is therefore highly skin friendly. This short stretch bandage is durable and extremely sturdy. Comprilan can be washed many times at 40°C without losing its specific features. The bandage ends do not fray and the interlocking weave provides air permeability. Support of the muscle pump for elimination of lymphatic stasis

BSN7718700 4 cm x 5 m
BDF0102600 6 cm x 5 m
BDF0102700 8 cm x 5 m
BDF0102800 10 cm x 5 m
BDF0102900 12 cm x 5 m



METRIMIST



- A freshly scented deoderizer spray for eliminating odors.
- Natural Aromatic Deodorizer
- Environmentally friendly alternative to aerosol air fresheners
- Eliminates offensive odors by cleaning the air, not masking the smell

Features:

- Biodegradable and nontoxic
- Nonflammable
- Cost effective

METRMX1158 8OZ 12/CS

ODORGON



Odorgon has a long history of success with no allergy complaints, or reactions. Still, there are situations where any inherent smell of a product is not acceptable, so Quest fills this need. Additional uses include: hunting, fishing & security needs. It's an interesting phenomenon when people use odor eliminators, they want to smell something! So using products with no scent may require training to prevent excessive use.

CHEM83006 114ML 12/CS

PODIATRY PROFESSIONAL TOWELS

Graham Medical's Podiatry Professional Towel comes in exciting blue or mauve outlined footprints to brighten exam environments!

SIZE: 13.5" x 18"
500/CS

GRAHAM70190N	3 PLY	FOOTPRINT BLUE
GRAHAM70191N	3 PLY	FOOTPRINT MAUVE
GRAHAM70192N	3 PLY TISSUE/POLY	FOOTPRINT BLUE
GRAHAM70194N	3 PLY TISSUE/POLY	FOOTPRINT MAUVE



2-PLY TISSUE EMBOSSED TOWEL

Graham Medical's Standard 2-Ply Tissue Embossed Towels are available in two weights and absorbencies to provide maximum protection against fluidings and staining medications.

SIZE: 13.5" x 18"
500/CS

GRAHAM70175N	2 PLY TISSUE	WHITE
--------------	--------------	-------




3-PLY TISSUE EMBOSSED TOWEL

Graham Medical's Standard 3-Ply Tissue Embossed Towels come in various colors and can be used in a variety of applications.

SIZE: 13.5" x 18"
500/CS

GRAHAM70171N	3 PLY TISSUE	BLUE
--------------	--------------	------



DRY WIPES

Very soft and very strong wipes designed to quickly and gently cleanse the skin. These multi-purpose dry disposable washcloths are ideal for perineal care, but are also suitable for the whole body, and are more hygienic than textile cloth alternatives.

PROM74500 50/PK



DISPOSABLE UNDERPADS

TENA Underpads for protection are designed to absorb leakages, reduce odours and maintain dryness. They consist of a polyethylene layer combined with sealed edges to ensure leakage protection. The highly absorbent material is distributed evenly throughout the core.

TENA Underpads provide far greater levels of absorption and retention than draw-sheets or wadding-filled products. Available in a comprehensive range of sizes and absorbencies.



REGULAR ABSORBENCY

CASC350	12 pkgs of 25 Pcs/Case	43 x 61 cm (17 x 24")
CASC351	8 Pkgs of 25 Pcs/Case	58 x 61 cm (23 x 24")
CASC352	6 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")

EXTRA ABSORBENCY

CASC353	12 pkgs of 25 Pcs/Case	43 x 61 cm (17 x 24")
CASC354	8 Pkgs of 25 Pcs/Case	58 x 61 cm (23 x 24")
CASC355	6 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")

ULTRA ABSORBENCY

CASC357	4 Pkgs of 25 Pcs/Case	58 x 91 cm (23 x 36")
---------	-----------------------	-----------------------

BARIATRIC

CASC361	10 Pkgs of 10 Pcs/Case	91 x 91 cm (36 x 36")
---------	------------------------	-----------------------

"SURE STAY" BACKING keeps absorbent pad in place without damaging surfaces with tape, pins or adhesives.

CASC381	6 pkgs of 25 Pcs/Case	57 x 60cm (22.5 x 23.5")
---------	-----------------------	--------------------------



TABLE/EXAM PAPER

CREPE

17"	WEST0101	125'	12/CS
18"	BOWB0118	125'	12/CS
18"	GRAHAM70002N	125'	12/CS
21"	GRAHAM70004N	125'	12/CS
24"	GRAHAM70007N	125'	12/CS

SMOOTH

18"	BOWB0218	225'	12/CS
18"	GRAHAM70016N	225'	12/CS
21"	GRAHAM70018N	225'	12/CS
24"	GRAHAM53216	225'	12/CS



FACIAL TISSUE 2 PLY

BOWB08300	100/BX	36BX/CS
BOWB08500	80/BX	135BX/CS



WAXED PAPER CUPS

BOWB113	4OZ	100/BG	10BG/CS
BOWB115	5OZ	100/BG	10BG/CS



TOILET TISSUE

ENTSCO5131	
2 PLY	48RL/CS



PAPER TOWELS

ENTKIM01960	
Multifold	
175/PK	25PK/CS





BOWERS
MEDICAL SUPPLY
(A DIVISION OF RADION LABORATORIES)

Head Office

7198 Progress Way
Delta, BC V4G 1J2
Tel: (604) 946-7712
Toll Free: (800) 663-0047
Fax: (604) 946-7715

Alberta

11783 - 186th Street
Edmonton, AB T5S 2Y2
Tel: (780) 454-1666
Toll free: (800)561-3687
Fax: (780) 451-6762

Ontario

1124 Northside Road
Burlington, ON L7M 1H4
Tel: (905) 315-7455
Toll free: (888) 216 -4833
Fax: (905) 315 -7422

www.bowersmedical.com
info@bowersmedical.com