

workFLOTM

Mobility management of hygiene and safety solutions - Hygienic Roll Stand.



Easy-to-maneuver and lightweight design provides mobile access at various entry points throughout your facility.

Interchangeable dispensing solution for masks, gloves and sanitizer.

- Main back board fastens to mounting plate of the workFLO[®] device mount roll stand.
- Holds up to 3 modules. Tabs on modules slide easily into slots on back board.
- Locking plate on rear of the back board keeps module secure while allowing for easy exchange of modules.
- Easily modify messaging as needed using the 9" x 12" clear magnetic document holder. Accommodates both US standard 8.5" x 11" and European standard A4 paper.



Learn How You Can Adapt Your Existing workFLO[®] Device Mount Roll Stand


Contact Your Local Sales Rep

Craig Alward
calward@touchpointmed.com
Tel: 604 318 5783





Main Back Board Fastened to Device Mount Roll Stand
119-002210-0



Mask Cassette

119-002219-00

- Holds masks that are bulk packaged (not tissue box style).
- Dispense by grabbing the top most mask in the cassette and pull through window.
- Holds individual masks up to 7.5"L x 4"W.
- Locks into back board, ONLY in the middle position (1 per back board max).



Hand Sanitizer Cover

119-002216-00

- "One size fits most" hand sanitizer bottles up to 32 fl oz.
- Fits cylindrical bottles up to 4" diameter.
- Will not work with gallon jug style dispensers.
- Locks into back board in ANY position (3 per back board max).



Glove Box Cover

119-002217-00

- "One size fits most" qty 100 and qty 200 glove boxes.
- Fits glove boxes up to 9.5"L X 5.5"W X 4"D.
- Locks into back board in ANY position (3 per back board max).




Mask Box Cover

119-002218-00

- One size fits most" tissue box style mask boxes.
- Fits mask boxes up to 8"L X 4.5"W X 3.5"D.
- Locks into back board in ANY position (3 per back board max).

Equipment changes over time and our line of accessories can be added or removed to fit your application needs.



DEVICE MOUNT ROLL STAND WF-RDM-FH

Height Adjustment Range: 1065mm (42in) with rear handle

Weight Capacity: 15kgs (33lbs)

Casters: 62.5mm (2.5in) – 5 locking

Mounting Options: Slide Plate

Main Material: Aluminum alloy/Stainless steel

Cart Weight: 5.4kgs (12lbs)

Cart Footprint: Ø 604.5mm (Ø23.8in) maximum



Antimicrobial Surfaces on Entire Stand

Info@touchpointmed.com touchpointmed.com

BRINGING CARE FORWARD.™