

# M U L T I - P U R P O S E E Q U I P M E N T W A S H E R

You too can reduce Bacteria by 99.9%

As demonstrated in an independent laboratory study

- 5 Minutes per Cycle • All Stainless Steel • Portable
- Washes Two Wheelchairs at a Time • Direct from Factory
- Customers in U.S.A. - Canada - New Zealand

#### DIMENSIONS

- NEW MODEL 64X

L 58" • W 34.25" • H 60.50"

- Door Opening

W 29.75" • H 47"

#### CHEMICALS

- Detergent – Approximately 300 washes per gallon
- Rinse Agent – Approximately 300 washes per gallon

#### WARRANTY

- Parts – wall to wall 5 years, labor 90 days
- Product use life – over 35 years

#### CONTROLS

- One button operation
- Electronic solid-state timer
- Door Switch – Automatic off if door is open

#### ELECTRICAL

- 115 Volt, 15 Amp Circuit

#### WATER INLET

- Standard hose connection to hot water faucet

#### WATER DRAIN

- Standard hose connection pumped to sewer drain



Manufactured by:



**Medco**<sup>®</sup>  
EQUIPMENT

**To Order in Canada:**

1 . 800 . 663 . 0047

[www.bowersmedical.com](http://www.bowersmedical.com)

# M U L T I - P U R P O S E EQUIPMENT WASHER

Over 2,000 Washers Now in Use

## Highlights

- Dramatic bacteria reduction  
Independent lab confirms over 99.9%
- All stainless-steel construction
- Portable. Heavy-duty castors
- No hassle or risk for 30 day free trial  
and free delivery
- 5-minute cycle cleans and sanitizes
- Will wash two chairs at one time
- Easy installation – 115V-15 amp plug-in  
Standard hose connections
- 300 washes per gallon of chemicals  
(600 if two at a time)
- UL listed and C/UL listed
- Government approved FSS #V797P4012a

## Warranty

- "Wall-to-Wall" parts, 5 years, Labor, 90 days
- Built to last over 40 years!

Distributed by:



**BOWERS**  
MEDICAL SUPPLY  
(A DIVISION OF RADION LABORATORIES)



Accommodates  
two wheelchairs at once.



Cleans and sanitizes commode  
chairs and shower chairs.



Washes up to four walkers  
at a time.



Large, unruly laundry/  
trash bins are a snap.

Try it today -  
You will **LOVE IT!**